

Micro-Kleen³TM

Accelerated Kill-Time FormulaTM

Fungicidal · Bactericidal · Virucidal · Tuberculocidal

A Fast 3-Minute, Ready-To-Use, Disinfectant, Decontaminant & Cleaner
for Hospital, Medical, Surgical and Institutional use on
Inanimate Environmental & Equipment Surfaces

DIN 02317524 (Solution) • DIN 02319101 (Wipes)



DESCRIPTION: Micro-Kleen³TM is a ready to use cleaner disinfectant that provides effective disinfection in hospitals, surgical facilities, neonatal & infant care units, child care facilities, plasma pheresis and blood centers, GI/pulmonary endoscopy labs, dialysis units and the many other areas where strict control of cross contamination is critical.

- Micro-Kleen³'s surfactant rich formula will quickly loosen and clean both wet & dried bodily fluids and other organic and non-organic soil from surfaces. It safely and effectively disinfects environmental and equipment surfaces such as: operating room lights & tables, lasik devices, slit lamps, tonometer devices, IV pumps, chest pumps, EMS equipment, examination room fixtures, dialysis machine surfaces, infant care equipment & incubators, oxygen hoods, respiratory therapy & anesthesia equipment surfaces, laboratory surfaces as well as other contaminated items.
- Micro-Kleen³TM will pre-clean blood-soiled surgical instruments prior to handling during the terminal cleaning and reprocessing steps by SPD and nursing personnel. To use for this purpose; spray instruments until thoroughly drenched, allow to remain wet for 3-minutes and rinse under cool running water. Follow with an appropriate terminal cleaning, wrapping or packaging and sterilization process.



BACTERIOLOGY, VIROLOGY & MYCOLOGY

Micro-Kleen³TM solution provides broad spectrum *Accelerated Quick-Kill*TM of pathogenic microorganisms including:

Gram Positive & Gram Negative Bacteria:

Acinetobacter baumannii
Mycobacterium tuberculosis (TB)
Staphylococcus aureus
Pseudomonas aeruginosa
Salmonella choleraesuis (enterica)
Methicillin Resistant Staphylococcus Aureus (MRSA)
Vancomycin Resistant Enterococcus (VRE)

Viruses:

Both Non-enveloped & Enveloped viruses including:
Rotavirus
Hepatitis B Virus (HBV)
Hepatitis C Virus (HCV)
Influenza Virus
Infectious Rhinovirus
Human Immunodeficiency virus (HIV)
Herpes Simplex Virus
Poliovirus
Avian Influenza Virus (Bird Flu)

Pathogenic Fungus:

Trichophyton mentagrophytes

For technical data and product information contact:



Micro-Scientific Industries
16228 Carrs Landing Road
Winfield, BC V4V 1A9
1-888-253-2536 www.micro-scientific.com



Micro-Kleen³TM is safe for use on inanimate surfaces including plastics (polyvinyl chloride, polypropylene, polycarbonate and polystyrene), stainless steel, brass, aluminum, chrome, acrylic, vinyl, glass, Plexiglas[®], and baked-on painted surfaces. Leaves no hazardous chemical residue when used as directed. Contains no dangerous phenols, caustic hydrogen peroxide, toxic chlorine or staining dyes. Requires no special handling, ground transportation or special storage precautions.

Unlike many other spray bottles, the specially designed Micro-Kleen³TM sprayer head produces large liquid droplets without producing fine aerosolized (airborne mist) particles that may be dispersed into the air and subsequently aspirated.

Micro-Kleen³TM exceeds blood-borne pathogen requirements for disinfection in healthcare facilities and is formulated following strict Good Manufacturing Practices.



APPLICATIONS:

Micro-Kleen³TM is a very effective:

- Surface Disinfectant and Cleaner
- Instrument Decontaminant Solution
- Ultrasonic Decontamination / Cleaning Solution



PRODUCT	ORDER #	PACKED
---------	---------	--------

Solution:

709 ml. trigger spray bottle	MK12-024	12/case
3.8L. pour bottle	MK04-128	4/case
9.46L. with spigot	MK02-320	2/case

Wipes:

100 Large Kimberly-Clark [®] Wipes Saturated with Micro-Kleen ³ TM Solution Packed in a Dispenser Canister.	MKW06-100	6/case
---	-----------	--------

Accessories:

3.85L. size- Large bulk sprayer for application to larger surfaces.	MKBS01-001	Each
Metal Wipes Canister Wall Bracket	MKWB01-001	Each

Micro-Kleen³TM
Family



Distributed by:



Valuemed Professional Products Ltd.
(780) 454-1899 * Toll Free (800) 919-8444
www.valuemed.ca