

# Multi-Function Chemical Integrator



### UNEXPOSED

Integraph integrators are designed to change color from purple to green when all parameters for steam sterilization (time, temperature, and the presence of saturated steam, are met.

### CORRECT COLOR CHANGE NORMAL CYCLE

When both of the first two bars have turned completely to green, sterilization conditions have been met.

### CORRECT COLOR CHANGE EXTENDED CYCLE

When all three bars have turned completely to green, extended cycle conditions have been met.

### INSUFFICIENT

If the first two bars of the integraph have not completely turned to green, exposure to sterilization conditions are insufficient.




### EXCESSIVE - NORMAL CYCLE

If the final bar on the integraph turns completely to green during a normal cycle, excessive exposure to sterilization processing is indicated.

### SUPERHEATED - DRY

The indicator spots will turn dark brown if the steam is superheated.



### WET STEAM

If the purple/green color appears blurry, and the plastic laminate bubbles, wet steam conditions may be present.