Preparing your unit for MicroSure Shock

- 1) Ensure no attachments are on the lines (handpieces, air/water syringes tips, ultrasonic and/or piezo scalers).
- 2) Empty any remaining water from bottle.
- 3) Add 500 mL of fresh tap water to your bottle.
- 4) Add 50 mL of MicroSure Shock concentrate to the bottle. (1:10 dilution ratio)
- 5) Gently swirl in bottle to mix concentrate with water.
- 6) Connect bottle to dental unit and run each line for 2 minutes. Ensure that the water volume on the unit is turned to "High". Leave in lines for **10 TO NO MORE THAN 20 MINUTES**.
- 7) Flush the entire bottle of MicroSure Shock through the lines.
- 8) Reconnect a bottle with fresh water.
- 9) Flush an entire bottle of tap water through the lines to ensure all solution has been purged. Turn down the volume of water to its previous level of flow.
- 10) Test with MicroSureH2Os to ensure effective shock, and ensure bacterial count is safe for human consumption.
- 11) Register for an account and *Log results* at: <http://www.duwltesting.com/microsure>

Note: If you witness particles dislodging while purging after your shock treatment, you may want to consider a second shock before re-testing with MicroSureH2Os.

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