# GlucoCheck Multi-Patient Blood Glucose Test Strip

## IMPORTANT

Please read this information and GlucoCheck Multi-Patient Blood Glucose Monitoring System Owner's Booklet before using the GlucoCheck Multi-Patient Blood Glucose Test Strips.

## INTENDED USE

The GlucoCheck Multi-Patient Blood Glucose Test Strips are for use with the GlucoCheck Multi-Patient Blood Glucose Meter to quantitatively measure glucose in fresh capillary whole blood samples drawn from the fingertips. The GlucoCheck Multi-Patient Blood Glucose Monitoring System is intended for use outside the body (*in vitro* diagnostic use) and is intended for multiple-patient use in professional healthcare settings as an aid to monitor the effectiveness of diabetes control. The system is only used with auto-disabling, single use lancing device. It is not intended for use on neonates and is not for the diagnosis or screening of diabetes.

## TEST PRINCIPLE

Glucose in blood samples reacts with the chemical in the test strip to produce a small electrical current. The GlucoCheck Multi-Patient Blood Glucose Meter detects the electrical current which reflects the amount of glucose in the blood sample.

## STORAGE AND HANDLING

- Store vial in a cool and dry place between 1-30 °C. Do not freeze.
- Keep the vial of test strips away from direct sunlight or heat.
- Store unused test strips in their original vial to avoid damage or contamination.
- Push the lid down on the vial immediately after taking out a test strip to fully close the vial and maintain air tightness.
- Handle test strips only with clean and dry hands.
- Use the test strip immediately after taking it out of the vial.
- Do not bend, cut, or alter the test strips in any way.
- Do not force a test strip into the meter. Gently push it into the meter's test strip port.
- Apply only fresh capillary whole blood for testing.
- Use all of the test strips within the expiration date printed on the test strip box and vial label.
- Dispose of test strips past the expiration date immediately. Using test strips past their expiration date can produce incorrect test results.
- Test strips in new, unopened vials and test strips in vials that have been opened can be used up until the expiration date printed on the test strip box and vial label if the test strips are used according to its storage and handling methods.

## WARNINGS AND PRECAUTIONS

- Inaccurate results may occur in patients undergoing oxygen therapy.
- Keep test strips and the test strip vial away from children. The test strips and vial cap may be choking hazards. Drying agents in the vial cap may be harmful if inhaled or swallowed and may cause skin or eye irritation.
- Test strips are for single use only. Do not reuse.
- If the test strip does not absorb the blood sample properly, please contact your authorized sales representative.

## BLOOD SAMPLE COLLECTION PROCEDURE

Wash hands and sample site with soap and warm water. Rinse and dry thoroughly before collecting the blood sample with a lancing device.

### Fingertip Site Blood Sampling

Unscrew the lancing device tip. Place the loaded lancing device against the side of the fingertip and press the release button. Massage the fingertip to obtain a round drop (at least 0.5 µL, actual size: ●) of blood. Apply test strip tip to the blood sample.

### Alternate Site Blood Sampling (forearm, palm, thigh, calf)

Select a clean, soft and fleshy sample site area free of visible veins and hair and away from bones. Gently massage the sample site to help blood circulation to minimize result differences between fingertip and alternative site sampling. Firmly press and hold the lancing device against sample site. Wait until the skin surface under the lancing device changes color. Then press the release button while continuing to apply pressure.

Keep holding the lancing device against your skin until sufficient (at least 0.5 µL) amount of blood is drawn. Carefully lift the lancing device away from your skin.

## Acceptable Situations for AST

When your blood glucose levels are stable

- Fasting period
- Before a meal
- Before going to bed

## Situations Requiring Fingertip Test

When your blood levels are unstable

- During the two (2) hours after a meal or exercise
- When sick or when glucose levels seem quite lower than test value
- When hypoglycemia is not well recognized
- When insulin has the biggest effect
- During two (2) hours after an insulin injection

## CAUTION

Alternative site and fingertip results may differ significantly due to rapid changes in the glucose level after meals or exercise, hypoglycemic symptoms, or effects of drugs such as insulin. Use a fingertip sample site if you suffer from hypoglycemia or have experienced hypoglycemic shock or symptoms. For instructions on how to obtain samples from alternative sites, please refer to the AST section of the Owner's Booklet of your meter.

## TEST PROCEDURE

- 1) Wash hands and sample site with soap and warm water. Rinse and dry thoroughly.
- 2) Insert the test strip into the port with contact bars facing upwards. Push the strip in gently until the meter beeps.
- 3) The  symbol will appear.
- 4) Use lancing device to obtain blood sample. Sample must be at least 0.5 µL (actual size: ●) to fill the test strip confirmation window. When the symbol appears on display, apply blood sample to edge of the narrow end of the test strip until meter beeps. If confirmation window has not filled, Er4 message or inaccurate test results may appear.
- 5) Meter will count down from five-to-one (5-to-1) on the display. Test result, time, and date will appear and automatically be stored in the meter's memory. Remove used test strip from port. Meter will turn off after three (3) seconds.

## TEST RESULTS

The GlucoCheck Multi-Patient Blood Glucose Meter will display results between 1.1-33.3 mmol/L.

### Normal Blood Glucose Results

A recommended target for fasting blood glucose or glucose before meals is below 7.0 mmol/L.

A recommended target for blood glucose 2 hours after eating is below 7.8 mmol/L.<sup>1</sup>

### Low Blood Glucose Results

If the test result is below 1.1 mmol/L, **Lo** will appear on the display indicating hypoglycemia (low blood glucose). You should follow the appropriate treatment recommendations of your healthcare professional.

### High Blood Glucose Results

If the test result is above 33.3 mmol/L, **HI** will appear on the display to indicate hyperglycemia (high blood glucose). If so, follow hyperglycemia treatment recommendations of your healthcare professional.

### Unexpected Results

Low or high blood glucose readings can indicate a potentially serious medical condition. If your results are unusually high or low, or do not match the way you feel, repeat the test with a new test strip. If your reading is inconsistent with your symptoms or your result is less than 3.3 mmol/L or higher than 13.3 mmol/L, contact your healthcare professional.

## Please note that:

- An abnormally high or low red blood cell count (hematocrit level over 65 % or below 15 %) may produce inaccurate results.
- Severe dehydration (excessive water loss) may cause inaccurate results. If you believe you are suffering from severe dehydration, consult your healthcare professional immediately.
- Altitude of up to 3,000 m (10,000 ft) above sea level has no effect on the performance of the test strip.
- Interferences: Paracetamol, ascorbic acid (vitamin C), uric acid and other reducing substances (when occurring in normal blood or normal therapeutic concentrations) do not significantly affect results. However, abnormally high concentrations in blood may cause inaccurate high results.
- Blood samples that contain a high concentration of dissolved oxygen may lower the test result.
- Discard used test strips properly in an appropriate container.

## METER AND TEST STRIP PERFORMANCE CHECK

The GlucoCheck Blood Glucose Control Solution (Control A and/or B) contains a known amount of glucose that reacts with the GlucoCheck Multi-Patient Blood Glucose Test Strip in combination with the GlucoCheck Multi-Patient Blood Glucose Meter to ensure they are working properly together and the correct testing procedure is being followed. You may run a check when you:

- Want to practice the test procedure using the control solution instead of blood.
- Use the meter for the first time.
- Open a new vial of test strips.
- Have symptoms that are inconsistent with your blood glucose test results.
- Believe your test results are not accurate.
- Suspect your meter and test strips are not performing properly.
- Drop or damage the meter.

If your control solution test results do not fall within the range printed on the test strip vial, repeat the test. Out of range results may be due to one or more of the following factors:

- Error in performing the test.
- Expired or contaminated control solution.
- Expired or damaged test strip.
- Failure to shake control solution bottle.
- Failure to discard first drop of control solution and wipe bottle tip clean.

If results continue to fall outside the range printed on the vial, the GlucoCheck Multi-Patient Blood Glucose Test Strip and Meter may not be working properly. If so, do not use your system and contact your authorized GlucoCheck Multi-Patient sales representative.

## CHEMICAL COMPOSITION

Each GlucoCheck Multi-Patient Blood Glucose Test Strip contains the following reagents:

- Glucose oxidase (*Aspergillus* sp.): 2.7 units
- Hexaamineruthenium (III) chloride: 45.7 µg
- Other ingredients: 1.6 µg

## PERFORMANCE CHARACTERISTICS

The performance of the GlucoCheck Multi-Patient BGM System has been evaluated in laboratory and in clinical tests.

## ACCURACY

The GlucoCheck Multi-Patient BGM System is calibrated to yield results equivalent to plasma glucose concentrations. The accuracy of the GlucoCheck Multi-Patient BGM System was assessed by comparing blood glucose results obtained by patients with those obtained using a YSI Model 2300 Glucose Analyzer, a laboratory instrument. The following results were obtained by diabetic patients at clinic centers.

Slope	0.946
Y-intercept	0.37 mmol/L
Correlation coefficient (r)	0.994
Number of sample	600
Range tested	1.6 - 27.1 mmol/L

Accuracy results for glucose concentration < 5.55 mmol/L

Within ± 0.28 mmol/L	Within ± 0.56 mmol/L	Within ± 0.83 mmol/L
100/186 (53.8 %)	169/186 (90.9 %)	180/186 (96.8 %)

Accuracy results for glucose concentration ≥ 5.55 mmol/L

Within ± 5 %	Within ± 10 %	Within ± 15 %
266/414 (64.3 %)	395/414 (95.4 %)	409/414 (98.8 %)

## PRECISION

Precision studies were performed in a laboratory using the GlucoCheck Multi-Patient BGM System.

<i>Within Run Precision</i>		
*Blood avg.	2.1 mmol/L	SD=0.1 mmol/L
	3.2 mmol/L	SD=0.1 mmol/L
	6.7 mmol/L	CV=3.6 %
	9.7 mmol/L	CV=2.8 %
	16.8 mmol/L	CV=3.2 %


<i>Total Precision</i>		
*Control avg.	2.2 mmol/L	SD=0.1 mmol/L
	6.7 mmol/L	CV=3.5 %
	17.7 mmol/L	CV=2.6 %

This study shows that there could be variation of up to 3.6 %.

## Reference

1. Canadian Diabetes Association, 2018. Clinical Practice Guidelines for the Prevention and Management of Diabetes in Canada.

## DESCRIPTION OF SYMBOL

	<i>In vitro</i> diagnostic medical device
	Do not re-use
	Consult instructions for use
	Temperature limit
	Manufacturer
	Batch code
	Catalogue number
	Use-by date
	Quantity

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🏭 **Innovatek Medical Inc.**