# GlucoCheck Blood Glucose Control Solution

## IMPORTANT

Please read the following information and the GlucoCheck Multi-Patient Blood Glucose Monitoring System Owner's Booklet before using the GlucoCheck Blood Glucose Control Solution.

## INTENDED USE

The GlucoCheck Blood Glucose Control Solutions are red in color and are to be used to check that the GlucoCheck Multi-Patient Blood Glucose Meter and the GlucoCheck Multi-Patient Blood Glucose Test Strips are working together properly and that you are performing the test correctly.

## SUMMARY AND EXPLANATION

GlucoCheck Blood Glucose Control Solutions are solutions containing a specific amount of glucose that can be tested with GlucoCheck Multi-Patient Blood Glucose Test Strips. The GlucoCheck Blood Glucose Control Solutions are for outside the body (*in vitro* diagnostic use) and for professional use only. These solutions may be used to check the performance of your meter and test strip or your testing technique. The test results you obtain should fall within the Control Solution Range printed on the vial of test strips that you are using.

## CAUTION

The control solution ranges printed on the test strip vial are for GlucoCheck Blood Glucose Control Solutions only. These ranges are used to check the meter and test strip performance.

**These are not recommended ranges for your blood glucose level.**

## Use GlucoCheck Glucose Control Solution when:

- You want to practice the test procedure, instead of using blood.
- You use your GlucoCheck Multi-Patient Blood Glucose Meter for the first time.
- You begin using a new vial of test strips.
- You suspect that the GlucoCheck Multi-Patient Blood Glucose Meter or the GlucoCheck Multi-Patient Blood Glucose Test Strips are not working properly.
- You think your test results are inaccurate or do not reflect how you feel.
- You drop your meter.

## STORAGE AND HANDLING

Store the GlucoCheck Blood Glucose Control Solutions at temperatures between 8-30°C/46-86°F tightly capped in the bottle. It must be kept at room temperature (20-25°C/68-77°F) before testing to obtain correct results. Do not refrigerate or freeze the control solution. Use before the expiration date shown on the bottle and within 3 months (90 days) of opening.

## WARNING AND PRECAUTIONS

- The GlucoCheck Blood Glucose Control Solutions are for *in vitro* diagnostic use only.
- Use only GlucoCheck Blood Glucose Control Solutions with your GlucoCheck Multi-Patient Blood Glucose Meters.
- Do not swallow or inject the control solution.
- Check the expiration date on the bottle. Do not use if expired.
- If the test results keep falling outside the specified range continually, do not use your system to test blood.
- Control solutions, meter and test strips must be brought up to room temperature (20-25°C/68-77°F) before testing.
- Shake the bottle well before use, discard the first drop of solution, and wipe the tip with a clean tissue before applying the solution to the GlucoCheck Multi-Patient Blood Glucose Test Strip.
- This solution contains a dye that may stain clothes.

## TEST PROCEDURE

### Step 1

Insert the GlucoCheck Multi-Patient Blood Glucose Test Strip into the strip port of your meter.

### Step 2

Shake the control solution vial well before using.

### Step 3

Discard the first drop of control solution you squeeze and wipe the dispenser tip clean.

### Step 4

Apply the control solution at the tip of the strip until the meter beeps.

### Step 5

Test results appear after the meter counts down from 5 to 1.

Compare your control result with the appropriate control solution range printed on the test strip vial. If the result falls within the specified range, you are ready to use the meter and the test strips to test your blood.

### Step 6

Remove the used test strip from the port. The meter will turn off automatically.

## INTERPRETATION OF OUT OF RANGE RESULTS

At room temperature, GlucoCheck Blood Glucose Control Solution test results should fall within the Control Solution Range printed on the GlucoCheck Multi-Patient Blood Glucose Test Strip vial about 95% of the time.

If your control solution test results do not fall within the expected range printed on the test strip vial, repeat the test. Out-of-range results may be due to one or more of the following factors:

- An error in the way you are doing the test.
- Expired or contaminated control solution.
- Damaged or outdated test strips.
- A problem with the meter.
- Control solution that is too warm or too cold.
- Failure to shake the control solution vial vigorously and to wipe the vial tip clean.

## IMPORTANT

If you continue to get control solution test results that fall outside the range printed on the test strip vial, the system may not be working properly. Do not use the system to test your blood. If you are unable to resolve the problem, contact your authorized sales representative.

## CHEMICAL COMPOSITION

Each 4 mL vial of GlucoCheck Blood Glucose Control Solution contains:

- D-Glucose
- Preservative
- FD&C red dye

## DESCRIPTION OF SYMBOLS

	<i>In vitro</i> diagnostic medical device
	Consult instructions for use
	Temperature limit
	Manufacturer
	Batch code
	Catalogue number
	Use-by date

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**Innovatek Medical Inc.**