MEDICAL

• 🔁 🖥

Guide Book Effective January 1, 2022

midmark

General Information

Ordering Notes	
Freight and Shipping	
Midmark [®] and Ritter [®] Medical Equipment Warranty	5
Midmark Chairs + Tables	<i>,</i> –
Midmark 641 Procedure Chair	
Midmark 630 HUMANFORM® Procedure Chair	
Midmark 647 Podiatry Chair	
Midmark 646 Podiatry Chair	
Midmark 626 Barrier-Free® Examination Chair	
Midmark 640 Pediatric Examination Table	
Upholstery Color Palette	
Ritter Chairs + Tables	
Ritter 230 Procedure Chair	
Ritter 244 Barrier-Free® Examination Table	
Ritter 225 Barrier-Free® Examination Chair	
Ritter 224 Barrier-Free® Examination Chair	
Ritter 204 Manual Examination Table	
Ritter 203 Treatment Table	
Ritter 95 Treatment Table	
Contine	
	20
Midmark 425 Physician Seating (Hand Operated)	
Midmark 427 Physician Seating (Foot Operated)	
Ritter 279 Physician Seating (Hand Operated)	
Ritter 278 Physician Seating (Foot Operated)	
Ritter 277 Physician Seating with Backrest	
Ritter 276 Physician Seating	
Ritter 275 Physician Seating with Backrest	
Ritter 274 Physician Seating	
Ritter 273 Physician Seating with Backrest	
Ritter 272 Physician Seating	
Ritter 271 Physician Seating with Backrest	
Ritter 270 Physician Seating	
Ritter 280 Patient Side Chair	
Ritter 282 Patient Bariatric Side Chair	
Ritter 281 Blood Draw Chair	
Lighting	
Ritter 255 LED Procedure Light	
Ritter 253 LED Exam Light	
Ritter 250 LED Exam Light	
Diagnostic Producto	
Diagnostic Products	20
Midmark Diagnostic Medical Products Warranty	
Midmark Diagnostic Medical Products Ordering Information	
Midmark HEINE® Diagnostic Instruments Limited Warranty	
Midmark Diagnostic Medical Products Support Information	
Midmark Digital Vital Signs Device, IQvitals® Zone™ and Supplies	
Midmark Digital Spirometer, Software and Supplies	
Midmark IQecg [®] , Software and Supplies	
Midmark HEINE Diagnostic Instrument Wall Kits and Accessories	
Midmark HEINE Diagnostic Instrument Desktop Kits and Accessories	
Midmark Synthesis Instrument Panel 11-12	
Midmark HEINE Diagnostic Instrument Replacement Parts and Supplies	
Midmark HEINE Gamma XXL LF Sphygmomanometer Kit and Accessories	
Midmark HEINE Gamma G5 Sphygmomanometer Kit and Accessories	
Midmark IQcart® Diagnostic Cart Accessories and Supplies	

HEINE is a registered trademark of Heine Optotechnik GmbH & Co KG.

Midmark[®] Workstations

Midmark Workstation Warranty	
Midmark 6211 Secure Laptop Workstation	
Midmark 6212 Flat Panel PC Workstation	60
Midmark 6213 Flat Panel Secure PC Workstation	61
Midmark 6214 Procedure Workstation	62
Midmark 6215-001 Compact Workstation	63
Midmark 6215-002 Compact Rectangular Workstation	64
Midmark 6215-003 Compact Kidney Workstation	65
Midmark 6217 Secure Laptop Workstation + Mouse Tray	
Midmark 6218 Flat Panel PC Workstation + Mouse Tray	67
Midmark 6219 Flat Panel Secure PC Workstation + Mouse Tray	
Midmark 6281 Workstation	
Midmark 6282 Workstation	70
Midmark 6283 Workstation	

Midmark Synthesis® Mobile Solutions

Midmark Synthesis Cabinetry Warranty	
Midmark Synthesis Cabinetry Color Selections	
Midmark Synthesis Mobile Treatment Cabinet M2	
Midmark Synthesis Mobile Treatment Cabinet M4	
Midmark Synthesis Mobile Treatment Cabinet M5	
Midmark Synthesis Mobile Procedure Cart P1	
Midmark Synthesis Mobile Procedure Cart P2	
Midmark Synthesis Mobile Procedure Cart P3	
Midmark Synthesis Mobile Waste Cart W2	
Midmark Synthesis Mobile Linen Cart L1	
Ritter 261 Waste Receptacle	

Instrument Processing

Ritter M11 [®] Steam Sterilizer	
Ritter M9® Steam Sterilizer	
Ritter M9D Manual Sterilizer	
Sterilizer Accessories	
VistaCool™ Direct-to-Drain System	
QuickClean® Ultrasonic Cleaner QC6/QC6R	
QuickClean® Ultrasonic Cleaner QC3/QC3R	
QuickClean® Ultrasonic Cleaner QC1	
Ultrasonic Parts and Accessories	

Midmark Customer Experience

For order status, requests and product information, contact Medical Customer Experience: 1.844.856.1233

Email: allcustservice-medical@midmark.com

To place an order, email: Orders-Medical@midmark.com

For Synthesis® Cabinetry order placement or information, email: cabinetry-medical@midmark.com

Midmark Technical Service

For service part identification, technical/operational questions, and warranty and product issues/complaints, contact Technical Service: **1.844.856.1230** Sterilizer + Ultrasonic Products: sterilizers.ts@midmark.com

Diagnostic Products: techsupport@midmark.com

Procedure + Exam Chairs, Seating, Lighting, Workstations, Diagnostic Instrument Wall/Desk Kits, Otoscopes,

Ophthalmoscopes + Sphygmomanometers: examroomeq.ts@midmark.com

Special Notes

1. ORDERING—Please provide the following information on all orders:

- Purchase order number
- Bill-to and ship-to addresses
- Product models and/or part number(s)
- Description including upholstery color and/or
- panel/countertop color
- Quantity needed
- Pricing
- Required shipping date
- Special shipping requirements (see details below)

DROP-SHIP REQUIREMENTS

- Please specify if the customer has a loading dock or has special shipping requirements.
- Receiving contact name, email and phone number is required.
- 2. SPECIAL DESIGN—Orders for equipment of special design, construction, upholstery or special finish are not subject to cancellation after your acceptance of our acknowledgment.
- 3. SPECIAL SHIPPING REQUIREMENTS—These services are available to enhance delivery to your customers. Please call for pricing.
 - APPOINTMENT REQUIRED: Carrier will call delivery contact and mutually agree upon delivery appointment. This is not a guaranteed time/date.
 - LIFTGATE: Carrier uses liftgate/ramp to get product to ground level when facility does not have a dock. Facility is responsible for taking product once at ground level.
 - INSIDE DELIVERY: Carrier unloads product and delivers inside the first set of doors only, provided entrance to facility can accommodate freight as packaged (including skids) and there are no steps.
 - MIDMARK DELIVERY SERVICE: Carrier sets product in place per the delivery contact. Carrier is responsible for assembling product and removal of all packaging, trash and skids from the delivery site (Synthesis[®] Cabinetry does not include installation).
- 4. FREIGHT DAMAGE—Merchandise arriving with external damage should be so described on the bill of lading. In case of concealed damage, carrier must be notified within five days after receipt and claims for damage instituted.
- 5. ORDER CHANGES—All changes must be submitted in writing. Any changes made within five days of the order shipping will result in the order being moved to the next available ship date.

- 6. ORDER CANCELLATION—Any order, if cancelled within 10 days before the scheduled shipping date, will be charged a 20% cancellation charge.
- 7. RETURN OF MERCHANDISE— A Return Material Authorization (RMA) must be obtained for the return of any item. If return is authorized, transportation charges must be prepaid. Risk of loss with respect to returned goods shall remain with purchaser until receipt by Midmark. Items returned for credit are subject to a charge of 25% or more to cover inspection, damage in transit, reconditioning and repacking. Items may be returned within 6 months (except service parts within 90 days) from the original ship date if properly packaged and in a good, unused condition. RMA # should be written on the outside of the packaging. Diagnostic products may be returned within 30 days if in a good, unused condition and in original sealed packaging.

Note: Certain items are non-returnable: Specials, Synthesis Cabinetry, Midmark[®] Workstations, Discontinued Items or Parts, Custom or Non-Standard Colors, Service Upholstery and Service Parts if value is less than \$50, and Diagnostic Consumables; e.g., Electrodes, Prep Kits, Pneumotach Mouthpieces, Temperature Probe Covers, Nose Clips, Thermal Paper, all Disposable Items and Diagnostic Items valued less than \$50 per item.

8. SPECIAL ORDER REQUIREMENTS—Any order requiring any special deviation from standard products will carry the following disclaimer:

We hereby accept and acknowledge the above captioned order. However, because of the special nature of this order in that this equipment is being supplied upon special design specifications, special materials and/or field modifications, as per your request, the standard warranties of Midmark and Ritter on its products are not applicable to this order, and this acceptance and acknowledgment is conditioned upon your indemnification and promise to hold Midmark and Ritter, its officers, employees and agents, harmless in respect to any and all warranty claims, damages, both direct and consequential, expenses, costs and any liability arising from the use of or relating to this special equipment.

9. IMPROVE PRODUCT QUALITY—Midmark Corporation reserves the right to change design and construction without notice.

Midmark[®] + Ritter[®] Medical Equipment (Tables, Chairs, Physician Seating, Lights, Workstations and Sterilizers) Limited Warranty

SCOPE OF WARRANTY: Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace components of the domestic and international medical products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

Midmark warrants to the original retail purchaser that during the applicable warranty period it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications.

APPLICABLE WARRANTY PERIOD: The applicable warranty period, measured from the date of invoice to the original retail purchaser of the product and shall be one (1) year for all warranted products and components. QuickClean® Ultrasonic Cleaners are warranted for a period of three (3) years. Workstations have a 5-year warranty including all parts, except power system components and coiled power cords that have a 1-year warranty.

OBTAINING WARRANTY SERVICE: Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via email at midmark.com; by phone at 1.800.MIDMARK or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380.

It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

EXCLUSIONS: This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/ Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Products considered to be of a consumable or sterile nature;
- (4) Accessories or parts not manufactured by Midmark;
- (5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (6) Costs and expenses of routine maintenance and cleaning; and
- (7) Representations and warranties made by any person or entity other than Midmark.
- (8) Matching of color, grain or texture except to commercially acceptable standards;
- (9) Changes in color caused by natural or artificial light;
- (10) Custom manufactured products;
- (11) Alterations or modifications to the product by any person or entity other than Midmark; and

(12) Products that would otherwise by covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market.

SOFTWARE; WITH RESPECT TO SOFTWARE THAT IS A PRODUCT OR COMPONENT THEREOF, MIDMARK DOES NOT WARRANT THAT THE SOFTWARE: (1) IS ERROR FREE; (2) CAN BE USED WITHOUT PROBLEMS OR INTERRUPTIONS; OR (3) IS FREE FROM VULNERABILITY TO INTRUSION OR ATTACK BY VIRUSES OR OTHER METHODS.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER: MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN.

WARRANTY DISCLAIMER: THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS: No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

NO AUTHORIZATION: No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

Warranty information included in product manuals.

Medical Guide Book Effective January 1, 2022



Midmark® 641 Procedure Chair

- Power height, back, foot and tilt
- 450 lb patient weight capacity
- 19" to 40" seat height adjustment (non-rotation models)
- 22.5" to 43.5" seat height adjustment (rotation models)
- Perimeter contouring for maximum access to patient
- Flip-out foot extension
- Hand and foot control standard
- Paper roll holder standard
- Choice of Premium/Premium UltraFree seamless upholstery
- Chair weight: 450 lb
- Shipping weight: 500 lb

Model	Description
Standard Models 641-002 002-0966-00-XXX ⁶	Procedure Chair Non-Programmable w/o receptacle, w/o top Upholstery top - Premium
641-002 002-10147-XXX ⁶	Procedure Chair Non-Programmable w/o receptacle, w/o top 641 Premium UltraFree 24" Upholstery top
641-003 002-0966-00-XXX ⁶	Procedure Chair Programmable w/receptacle, w/o top Upholstery top - Premium
641-003 002-10147-XXX ⁶	Procedure Chair Programmable w/receptacle, w/o top 641 Premium UltraFree 24" Upholstery top
Rotation Models 641-004 002-0966-00-XXX ⁶	Procedure Chair Non-Programmable w/rotation, w/o top Upholstery top - Premium
641-004 002-10147-XXX ⁶	Procedure Chair Non-Programmable w/rotation, w/o top 641 Premium UltraFree 24" Upholstery top
641-005 002-0966-00-XXX ⁶	Procedure Chair Programmable w/rotation, w/o top Upholstery top - Premium
641-005 002-10147-00-XXX ⁶	Procedure Chair Programmable w/rotation, w/o top 641 Premium UltraFree 24" Upholstery top

Accessories

Description

Arms/Armboards

9A179002-XXX^{5,6} 9A396001-XXX^{3,5,6} 9A396001-XXX^{3,6} 9A81001-XXX^{5,6} 9A82001-XXX^{5,6}

Casters 9A3300011

9A365001¹

Controls

9A408001² 9A408002² 9A408003² 9A408004²

Headrests

9A374001-XXX^{5,6} 9A374001-XXX⁶ 9A378001-XXX^{5,6} 9A382001-XXX^{5,6} 9A382001-XXX^{5,6} 9A386001-XXX^{5,6} 9A386001-XXX^{5,6} 9A390001-XXX^{5,6} 9A394001-XXX^{5,6} 9A394001-XXX^{5,6}

Instrument Trays/Work Shelves/IV Poles

9A427002⁷ 9A267003⁴ 9A77001⁴

Miscellaneous

Whiseenancous	
9A152010	3' power cord
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A406001	Security straps, plant installed
9A406002	Security straps, field installed
9A78001-XXX ^{4,5,6}	Vision block screen
9A416001	Foot section plastic cover
002-0966-00-999	Upholstery top - Premium Special Color

IV pole

¹ - Not available on rotation models

- ² Accessory supplied with hand and foot control (batteries included)
- ³ Interference with armboard accessories will likely occur
- ⁴ Requires 9A20500X base rail
- ⁵ Additional charge for Special, Premium UltraFree for accessory
- ⁶ See color palette on page 19
- ⁷ Not compatible on same side as 9A20500X base rail

Fixed armboard Chair arms Chair arms, UltraFree upholstery Articulating armboard Hand surgery armboard

Caster base, swivel wheels Casters

Wireless control, non-programmable, linear-style, plant installed Wireless control, programmable, linear-style, plant installed Wireless control, non-programmable, linear-style, field installed Wireless control, programmable, linear-style, field installed

Contoured headrest, Premium upholstery Contoured headrest, Ultrafree upholstery Flat headrest, Premium upholstery U-shaped headrest, Premium upholstery U-shaped headrest, UltraFree upholstery Oval-shaped headrest, Premium upholstery Rectangular-shaped headrest, UltraFree upholstery Rectangular-shaped headrest, UltraFree upholstery
Magnetic headrest, Premium upholstery Magnetic headrest, UltraFree upholstery

Swing arm instrument tray, extension top

Instrument tray, double arm

Midmark Chairs + Tables

7

Medical Guide Book Effective January 1, 2022



Midmark 630

HUMANFORM® Procedure Chair

- Power height, back, foot and tilt
- Motion Profiling® adjusts the chair as the patient moves to help maintain head and back support, and reduce the need for the patient to reposition as the chair back raises or lowers
- Active Sensing Technology[®] (AST) helps avoid impact or disruption while moving the chair by pausing chair movement when disrupting impact is detected
- Barrier-free design
- 650 lb patient weight capacity
- 19" to 40" seat height adjustment (non-rotation models)
- 21.5" to 43.5" seat height adjustment (rotation models)
- Hand and foot controls standard
- Available in 28" or 32" upholstery widths
- Chair weight: 450 lb
- Shipping weight: 525 lb

Model Non-Programmable	Description
630-010	Procedure Chair Non-Programmable w/receptacle, w/o top
002-10108-XXX ⁶	630 Premium 28" Upholstery top
630-010	Procedure Chair Non-Programmable w/receptacle, w/o top
002-10114-XXX ⁶	630 Premium UltraFree 28" Upholstery top
630-011	Procedure Chair Non-Programmable wireless, w/receptacle, w/o top
002-10108-XXX ⁶	630 Premium 28" Upholstery top
630-011 002-10114-XXX ⁶	Procedure Chair Non-Programmable wireless, w/receptacle, w/o top 630 Premium UltraFree 28" Upholstery top
Non-Programmable v	with Rotation
630-012	Procedure Chair Non-Programmable w/rotation, w/o top
002-10108-XXX ⁶	630 Premium 28" Upholstery top
630-012	Procedure Chair Non-Programmable w/rotation, w/o top
002-10114-XXX ⁶	630 Premium UltraFree 28" Upholstery top
630-013	Procedure Chair Non-Programmable wireless, w/rotation, w/o top
002-10108-XXX ⁶	630 Premium 28" Upholstery top

Model	Description
Programmable 630-020 002-10108-XXX ⁶	Procedure Chair Programmable w/receptacle, w/o top 630 Premium 28" Upholstery top
630-020	Procedure Chair Programmable, w/receptacle, w/o top
002-10114-XXX⁴	630 Premium UltraFree 28" Upholstery top
630-021	Procedure Chair Programmable wireless, w/receptacle, w/o top
002-10108-XXX ⁶	630 Premium 28" Upholstery top
630-021	Procedure Chair Programmable wireless, w/receptacle, w/o top
002-10114-XXX⁰	630 Premium UltraFree 28" Upholstery top
Programmable with	Potation
630-022	Procedure Chair Programmable, w/rotation, w/o top
002-10108-XXX ⁶	630 Premium 28" Upholstery top
630-022	Procedure Chair Programmable, w/rotation, w/o top
002-10114-XXX ⁶	630 Premium UltraFree 28" Upholstery top
630-023	Procedure Chair Programmable wireless, w/rotation, w/o top
002-10108-XXX ⁶	630 Premium 28" Upholstery top
630-023	Procedure Chair Programmable wireless, w/rotation, w/o top
002-10114-XXX ⁶	630 Premium UltraFree 28" Upholstery top
Upholstery Options	
002-10108-999	Upholstery top - 28" wide Premium Special Color
002-10114-999	Upholstery top - 28" wide UltraFree Special Color
002-10112-XXX ⁶	Upholstery top - 32" wide Premium
002-10112-999	Upholstery top - 32" wide Premium Special Color
002-10116-XXX ⁶	Upholstery top - 32" wide UltraFree
002-10116-999	Upholstery top - 32" wide UltraFree Special Color

9

Medical Guide Book Effective January 1, 2022

Accessories

Description

Fixed armboard

Articulating armboard

Hand surgery armboard

Chair arms, 28" wide Premium upholstery Chair arms, 28" wide I ternium upholstery Chair arms, 28" wide UltraFree upholstery Chair arms, 32" wide Premium upholstery Chair arms, 32" wide UltraFree upholstery

Wired hand control, non-programmable, field installed

Wired foot control, non-programmable, field installed

Wired hand control, programmable, field installed

Wired foot control, programmable, field installed

Rectangular-shaped headrest, Premium upholstery Rectangular-shaped headrest, UltraFree upholstery

Oval-shaped headrest, Premium upholstery Oval-shaped headrest, UltraFree upholstery

U-shaped headrest, Premium upholstery

U-shaped headrest, UltraFree upholstery

Standard headrest, Premium upholstery

Standard headrest, UltraFree upholstery

Extended headrest, Premium upholstery

Arm Systems/Armboards

9A550001-XXX^{5,6} 9A550001-XXX6 9A550005-XXX5,6 9A550005-XXX6 9A179002-XXX^{2,5,6} 9A81001-XXX^{2,5,6} 9A82001-XXX^{2,5,6}

Controls (additional)

9A552001 9A553001 9A552002 9A553002

Headrests

9A496001-XXX5,6 9A496001-XXX⁶ 9A500001-XXX5,6 9A500001-XXX⁶ 9A512001-XXX5,6 9A512001-XXX6 9A683001-XXX5,6 9A683001-XXX⁶ 9A684001-XXX5,6 9A684001-XXX6

Extended headrest, UltraFree upholstery Instrument Trays/Work Shelves/IV Poles

9A2670032 9A770012

Instrument tray, double arm IV pole

Knee Crutches 9

9A410001 ³	Standard knee crutches (crutches only)
9A410002 ³	Standard knee crutch system
9A411001 ³	Articulating knee crutches (crutches only)
9A411002 ³	Articulating knee crutch system
9A412001 ³	Knee crutch mounting brackets (brackets only)

Accessory Rails

9A521001	Seat accessory rails, surgical size, 28" wide upholstery, plant installed
9A522001	Seat accessory rails, surgical size, 32" wide upholstery, plant installed
9A567001	Seat accessory rails, surgical size, 28" wide upholstery, field installed
9A568001	Seat accessory rails, surgical size, 32" wide upholstery, field installed
9A557001	Accessory rails, 28" wide upholstery, field installed
9A558001	Accessory rails, surgical size, 28" wide upholstery, field installed
9A557002	Accessory rails, 32" wide upholstery, field installed
9A558002	Accessory rails, surgical size, 32" wide upholstery, field installed

Miscellaneous

9A5540011	Casters - set of 4
9A152010	3' power cord
9A157003 ²	Security side panels
9A5590014	Security straps, field installed
9A555001	Gynecology debris tray accessory
9A556001	Gynecology debris tray (tray only)
9A78001-XXX ^{5,6}	Anesthesia/vision block screen
AD 1779 ¹	Tamper-resistant duplex receptacle

¹ - Not available on rotation models

- ³ Not compatible with 9A521001, 9A522001, 9A567001 or 9A568001
- ⁴ Requires (1) 9A557001/9A557002/9A558001/9A558002 accessory rails and (1) 9A521001/9A522001 seat rails
- ⁵ Additional charge for Special, UltraFree color for accessory

⁶ - See color palette on page 19



² - Requires (1) 9A557001/9A557002/9A558001/9A558002 accessory rails



Midmark® 647

Podiatry Chair

- Power height, tilt and back
- 450 lb patient weight capacity
- 19" to 40" seat height adjustment (non-rotation models)
- 22.5" to 43.5" seat height adjustment (rotation models)
- 9" foot section extension
- Sealed foot control standard
- Integrated touch controls standard
- Chair arms standard
- Choice of Premium Seamless or UltraFree upholstery
- Protective foot cover standard with upholstery
- Chair weight: 433 lb (non-rotation models)
- Shipping weight: 515 lb (non-rotation models)
- Chair weight: 539 lb (rotation models)
- Shipping weight: 621 lb (rotation models)

Model	Description
647-001	Podiatry Chair non-programmable w/receptacles, w/o top
002-1119-XXX ⁴	Upholstery top - Premium
647-001	Podiatry Chair non-programmable w/receptacles, w/o top
002-10129-XXX ⁴	646/647 Premium UltraFree 28" Upholstery top
647-002	Podiatry Chair programmable w/receptacles, w/o top
002-1119-XXX ⁴	Upholstery top - Premium
647-002	Podiatry Chair programmable w/receptacles, w/o top
002-10129-XXX ⁴	646/647 Premium UltraFree 28" Upholstery top
647-003	Podiatry Chair non-programmable w/rotation, w/o receptacles and top
002-1119-XXX ⁴	Upholstery top - Premium
647-003 002-10129-XXX ⁴	Podiatry Chair non-programmable w/rotation, w/o receptacles and top 646/647 Premium UltraFree 28" Upholstery top
647-004	Podiatry Chair programmable w/rotation, w/o receptacles and top
002-1119-XXX ⁴	Upholstery top - Premium
647-004	Podiatry Chair programmable w/rotation, w/o receptacles and top
002-10129-XXX ⁴	646/647 Premium UltraFree 28" Upholstery top

11

Ritter Chairs + Tables

Seating

Lighting

Diagnostics

Accessories

Description

Instrument Trays/Work \$ 9A427002 ⁵ 9A267003 ² 9A77001 ²	Shelves/IV Poles Swing arm instrument tray, extension top Instrument tray, double arm IV Pole
Casters 9A330001 ¹ 9A365001 ¹	Caster base, swivel wheels Casters
Miscellaneous 9A418001 9A418002 9A418003 9A418004 9A422001-XXX ^{3,4} 9A420001 9A423001 9A423001 9A205003 9A205004 9A152010 002-1119-999	Wireless foot control, non-programmable, linear-style, plant installed Wireless foot control, non-programmable, linear-style, field installed Wireless foot control, programmable, linear-style, plant installed Wireless foot control, programmable, linear-style, field installed Vision block screen Plastic protective foot cover Paper roll holder Base rails, plant installed Base rails, field installed 3' power cord Upholstery top - Premium Special Color

Not available on rotation models
 Requires 9A20500X base rail
 Additional charge for Special or UltraFree accessory
 See color palette on page 19
 Not available on the second seco

⁵ - Not compatible on same side as 9A20500X base rail

Not available for sale in Canada



Midmark® 646

Podiatry Chair

- Power tilt and back
- 450 lb patient weight capacity
- 21" fixed seat height
- 9" foot section extension
- Sealed foot control standard
- Integrated touch controls standard
- Chair arms standard
- Choice of Premium Seamless or UltraFree upholstery
- Protective foot cover standard with upholstery
- Chair weight: 323 lb
- Shipping weight: 405 lb

Description

646-001 002-1119-XXX³

Podiatry Chair non-programmable w/receptacles, w/o top Upholstery top - Premium

646-001 002-10129-XXX³ Podiatry Chair non-programmable w/ receptacles, w/o top 646/647 Premium UltraFree 28" Upholstery top

Accessories

Description

Instrument Trays/Work Shelves/IV Poles

9A427002⁴ 9A267003¹ 9A77001¹ **Casters** 9A330001 9A365001 Swing arm instrument tray, extension top Instrument tray, double arm IV Pole Caster base, swivel wheels

Miscellaneous

9A422001-XXX ^{2,3}	Vision block screen
9A420001	Plastic protective foot cover
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A152010	3' power cord
002-1119-999	Upholstery top - Premium Special Color

Casters

¹ - Requires 9A20500X base rail

² - Additional charge for Special or UltraFree Upholstery color for accessory

³ - See color palette on page 19

⁴ - Not compatible on same side as 9A20500X base rail

Midmark Chairs + Tables



Midmark 626

Barrier-Free® Examination Chair

- 650 lb weight capacity
- Low height of 15.5" (uncompressed) to 37" high
- Clean Assist Roller System
- Power base and back
- Hand and foot controls
- Digital Scale (-003 and -004 models)
- Heated upholstery options (-002, -004 and -006 models)
- Quick Exam[®] and Return to Home
- Active Sensing Technology®
- Mist and Stone chair base color
- Lockout of hand and foot controls
- Chair weight: 435 lb
- Shipping weight: 485 lb

Model	Description
626-001	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2000-XXX ¹	626 Premium 28" Upholstery top
626-001	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2002-XXX ¹	626 Premium 32" Upholstery top
626-001	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2004-XXX ¹	626 Premium UltraFree 28" Upholstery top
626-002	Barrier-Free Examination Chair w/ Clean Assist Roller System,
	receptacles, drawer heater, pelvic tilt, wired controls, heated upholstery options, w/o top
002-2006-XXX ¹	626 Premium 28" Heated Upholstery top
626-002	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, heated upholstery options, w/o top
002-10201-XXX1	626 Premium 32" Heated Upholstery top

Model	Description
626-002	Barrier-Free® Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, heated upholstery options, w/o top
002-2011-XXX ¹	626 Premium UltraFree Heated 28" Upholstery top
626-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2000-XXX1	626 Premium 28" Upholstery top
626-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2002-XXX1	626 Premium 32" Upholstery top
626-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2004-XXX1	626 Premium UltraFree 28" Upholstery top
626-004	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2006-XXX1	626 Premium 28" Heated Upholstery top
626-004	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-10201-XXX1	626 Premium 32" Heated Upholstery top
626-004	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2011-XXX1	626 Premium UltraFree Heated 28" Upholstery top

BACK TO CONTENTS

15 NTENTS

Medical Guide Book Effective January 1, 2022

Model	Description
626-005	Barrier-Free® Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2000-XXX1	626 Premium 28" Upholstery top
626-005	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2002-XXX ¹	626 Premium 32" Upholstery top
626-005	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2004-XXX1	626 Premium UltraFree 28" Upholstery top
626-006	Barrier-Free Examination Chair w/ Clean Assist Roller System,
020 000	Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2006-XXX ¹	626 Premium 28" Heated Upholstery top
626-006	Barrier-Free Examination Chair w/ Clean Assist Roller System,
	Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-10201-XXX ¹	626 Premium 32" Heated Upholstery top
626-006	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer
002-2011-XXX1	heater, pelvic tilt, wired controls, w/o top 626 Premium UltraFree 28" Heated Upholstery top
	ozo rieman onarie zo rieated opholstery top

Model Arm Systems	Description
9A600002 ² 9A600001 ² 9A582001	Patient Support Rails+, field installed Patient Support Rails, field installed Patient Support Rails, receiver extension - pair; field install for 32" wide upholstery top

Knee Crutches

9A643001	Standard knee crutch kit w/ acc. receiver, field install
9A645001	Articulating knee crutch kit w/ acc. receiver, field install
9A641001	Mounting brackets, hardware w/ acc. receiver, field install

Vital Signs Accessories

9A478001	Digital Scale serial cable for Digital Vital Signs Device, 6' length, coiled
9A478002	Digital Scale serial cable for Digital Vital Signs Device, 15' length, straight
9A478003	Digital Scale serial cable for Digital Vital Signs Device, 30' length, straight
9A478004	Digital Scale serial cable for Digital Vital Signs Device, 50' length, straight
9A479001	Welch Allyn [®] Spot Vitals Lxi, 6' length coiled
9A479002 002-10064-00	Welch Allyn Spot Vitals Lxi, 15' length straight 626 chair connectivity kit

Miscellaneous

9A677001	Wireless hand and foot controls
9A610001	Stainless steel treatment pan
002-2015-XXX1	626 Premium UltraFree 32" upholstery top
002-10203-XXX1	626 Premium UltraFree heated 32" upholstery top
002-2008-XXX1	224/225/626 Flat 28" upholstery top
002-10361-XXX1	224/225/626 Flat UltraFree 28" upholstery top
002-10494-XXX1	224/225/626 Premium Thick Footrest Pad
002-10495-XXX1	224/225/626 Premium UltraFree Thick Footrest Pad

¹ - See color palette on page 19

² - 9A582001 receiver extensions required for use with 32" wide upholstery top

Welch Allyn is a registered trademark of Welch Allyn, Inc.

BACK	TO	100	NTE	NTS

Medical Guide Book Effective January 1, 2022



Midmark[®] 640

Pediatric Examination Table


- Pebble grey base and panel color only
- An upholstered end panel and beveled sides designed to help prevent falls
- Flat sides with no sharp edges
- Accurate and easy-to-use, built-in scale
- Ample storage and a hidden paper roll holder
- 42" W x 51.5" H x 23"D
- Table weight: 170 lb
- Shipping weight: 210 lb
- Working Height: 35.8" (90.8 cm)
- Table Top Length: 41.9" (106.4 cm)
- Table Top Depth: 20.8" (52.7 cm)

Model	Description
640-001-XXX ^{1,2}	Pediatric Examination Table w/ digital scale
640-001-999	Pediatric Examination Table w/ digital scale - special color
640-002-XXX ^{1,2}	Pediatric Examination Table w/ door and drawer lock with digital scale, plant installed
640-002-999	Pediatric Examination Table w/ door and drawer lock with digital scale, plant installed - special color

¹ Additional charge for UltraFree Upholstery/Accessory Color

² See color palette on page 19

Upholstery Color Palette



Example: 630-010 Chair 002-10108-843 Top 630 HUMANFORM® Procedure Chair - Oasis, Premium Upholstery Top Example: 276-001-814

276-001-814 276 Physician Seating, Stone Upholstery

Standard Upholstery Colors

Color is an important part of establishing your distinct brand. It's more than visual appeal. Color can transform your space into something you won't just see, but you will feel, too. We offer many upholstery color selections to complement the look, feel and intended function of your space.



UltraFree Upholstery Colors

There is nothing quite like the rich feel of leather, but it's too difficult to clean and maintain to be practical in the clinical environment. With UltraFree upholstery, you can have the look and feel of leather with PVC-free material. And now with seamless application, it is easier to clean and disinfect.

UltraFree Stone-865UltraFree Spice-863UltraFree Restful Path-868UltraFree Cranberry-869UltraFree Obsidian-867UltraFree Dark Linen-866UltraFree Harbor-861UltraFree Latte-870UltraFree Wheatgrass-862UltraFree Branch-864

19



Procedure Chair

- Power height, back, foot and tilt
- 450 lb patient weight capacity
- 22.5" to 40" seat height adjustment (non-rotation models)
- 26" to 43.5" seat height adjustment (rotation models)
- Knee-chest configurable
- Hand control standard
- Chair weight: 455 lb
- Shipping weight: 525 lb

М	od	el
1.4.1	ou	CI

Description

230-001	Procedure Chair w/o top
002-0824-XXX⁵	230 Seamless Upholstery top
230-001	Procedure Chair w/o top
002-10133-XXX⁵	230 Seamless UltraFree 28" Upholstery top
230-002	Procedure Chair w/ receptacle, w/o top
002-0824-XXX⁵	230 Seamless Upholstery top
230-002	Procedure Chair w/ receptacle, w/o top
002-10133-XXX⁵	230 Seamless UltraFree 28" Upholstery top
230-003	Procedure Chair w/ rotation, w/o top
002-0824-XXX⁵	230 Seamless Upholstery top
230-003	Procedure Chair w/ rotation, w/o top
002-10133-XXX⁵	230 Seamless UltraFree 28" Upholstery top

Accessories

Description

Fixed armboard

Casters

9A330001^{1,2} 9A365001²

Controls 9A371003

Headrests/Footrests

Larger 21" wide featre

Casters

Articulating armboard Hand surgery armboard

Caster base, swivel wheels

9A337001-XXX^{4,5} 9A51001-XXX^{4,5}

Larger 24" wide footrest Facial pad

Foot control, linear style

Instrument Trays/Work Shelves/IV Poles

9A427002 9A267003³ 9A77001³ Swing arm instrument tray, extension top Instrument tray, double arm IV pole

Adjustable arm system, plant installed Adjustable arm system, field installed

Adjustable arm system (arms only)

Adjustable arm system (receiver only), plant installed Adjustable arm system (receiver only), field installed

Knee Crutches

9A410001Standard knee crutches (crutches only)9A410002Standard knee crutch system9A411001Articulating knee crutches (crutches only)9A411002Articulating knee crutch system9A412001Knee crutch mounting brackets (brackets only)

Miscellaneous

9A152010 9A157003 9A205003 9A205004 9A214002 9A329003² 9A329004² 9A78001-XXX^{3,4,5} 002-0824-999 002-0950-XXX⁵ 002-0950-999 002-10135-XXX⁵ 3' Power Cord Security side panels Base rails, plant installed Base rails, field installed Security straps, field installed Urology drain system, plant installed Urology drain system, field installed Vision block screen Seamless Upholstery top Special Color Flat Upholstery top, 24" width Flat Upholstery top, 24" width

¹ - Not compatible with 9A337001

² - Not available on rotation models

³ - Requires 9A20500X base rail

⁴ - Additional charge for Special, UltraFree accessory

⁵ - See color palette on page 19

21



Ritter 244

Barrier-Free® Examination Table

- Power height, back and tilt
- Reverse Trendelenburg position
- 850 lb patient weight capacity
- 18" to 34" seat height adjustment
- 32" wide seat
- Hand control standard
- Table weight: 340 lb
- Shipping weight: 400 lb

Model

244-001

Description

244-001 002-0861-XXX¹

> Barrier-Free Examination Table w/o top Seamless UltraFree 32" Upholstery top

Seamless 32" Upholstery top

Barrier-Free Examination Table w/o top

Accessories

002-10174-XXX1

9A627001Adjustable arm system, field installed9A346001Casters9A347001Foot control9A1520103' Power Cord

Miscellaneous

002-0861-999

Seamless 32" Upholstery Top Special Color

¹ - See color palette on page 19



Ritter 225

Barrier-Free® Examination Chair

- 500 lb weight capacity
- Low height of 18" (uncompressed) to 37" high
- Power base and back
- Wired foot control
- Active Sensing Technology®
- Mist chair base color
- Chair weight: 409 lb
- Shipping weight: 465 lb

Model	Description
225-002	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top
002-2009-XXX ¹	224/225 Seamless 28" Upholstery top
225-002	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top
002-2008-XXX ¹	224/225/626 Flat 28" Upholstery top
225-002	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top
002-2013-XXX ¹	224/225 Seamless UltraFree 28" Upholstery top
225-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, w/o top
002-2009-XXX ¹	224/225 Seamless 28" Upholstery top
225-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, w/o top
002-2008-XXX ¹	224/225/626 Flat 28" Upholstery top
225-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, w/o top
002-2013-XXX ¹	224/225 Seamless UltraFree 28" Upholstery top
Arm Systems	
9A600002 9A600001	Patient Support Rails+, Field Installed Patient Support Rails, Field Installed
9A615001	Assist Arms
Knee Crutches	
9A643001	Standard Knee Crutch Kit w/ Acc. Receiver, Field Install
9A645001 9A641001	Articulating Knee Crutch Kit w/ Acc. Receiver, Field Install Mounting Brackets, Hardware w/ Acc. Receiver, Field Install
Miscellaneous	
9A610001	Stainless Steel Treatment Pan
9A638001	Supply Assistant
9A681001 002-10361-XXX1	225 Hand Control Accessory 224/225/626 Flat UltraFree 28" Upholstery top
002-10494-XXX ¹	224/225/626 Premium Thick Footrest Pad
002-10495-XXX1	224/225/626 Premium UltraFree Thick Footrest Pad

¹ - See color palette on page 19



Ritter 224

Barrier-Free® Examination Chair

- 500 lb weight capacity
- Low height of 18" (uncompressed) to 37" high
- Power base
- Wired foot control
- Active Sensing Technology® (-002 and -003 models)
- Chair weight: 409 lb
- Shipping weight: 465 lb

	Shipping weight. 403 ib
Model	Description
224-001 002-2009-XXX1	Barrier-Free Examination Chair - excludes drawer and debris tray w/o top 224/225 Seamless 28" Upholstery top
224-001 002-2008-XXX ¹	Barrier-Free Examination Chair - excludes drawer and debris tray w/o top 224/225/626 Flat 28" Upholstery top
224-001 002-2013-XXX ¹	Barrier-Free Examination Chair - excludes drawer and debris tray w/o top 224/225 Seamless UltraFree 28" Upholstery top
224-002 002-2009-XXX ¹	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top 224/225 Seamless 28" Upholstery top
224-002 002-2008-XXX ¹	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top 224/225/626 Flat 28" Upholstery top
224-002 002-2013-XXX ¹	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top 224/225 Seamless UltraFree 28" Upholstery top
224-003 002-2009-XXX ¹	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt w/o top 224/225 Seamless 28" Upholstery top
224-003 002-2008-XXX ¹	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt w/o top 224/225/626 Flat 28" Upholstery top
224-003 002-2013-XXX ¹	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt w/o top 224/225 Seamless UltraFree 28" Upholstery top

Model	Description
Arm Systems 9A600002 9A600001 9A615001	Patient Support Rails+, field installed Patient Support Rails, field installed Assist arms
Knee Crutches	

9A643001	Standard knee crutch kit w/ acc. receiver, field install
9A645001	Articulating knee crutch kit w/ acc. receiver, field install
9A641001	Mounting brackets, hardware w/ acc. receiver, field install

Miscellaneous

9A610001	Stainless steel treatment pan
9A638001	Supply assistant
002-10361-XXX1	224/225/626 flat UltraFree 28" upholstery top
002-10494-XXX1	224/225/626 Premium thick Footrest Pad
002-10495-XXX1	224/225/626 Premium UltraFree thick Footrest Pad

¹ - See color palette on page 19

25



Manual Examination Table

- 500 lb patient weight capacity
- Electrical receptacle standard
- Styled Seamless or Soft Touch stitched upholstery
- Pass-through drawer design, no need to specify RH or LH
- Manual back with air-spring assist
- Please specify upholstery/drawer panel color
- Mist table base color
- Table weight: 250 lb
- Shipping weight: 300 lb

Model Description 204-011 Manual examination table w/o top 002-10137-XXX² 204 Soft Touch 28" upholstery top 204-011 Manual examination table w/o top 002-10138-XXX² 204 Upholstery top - Styled, Seamless 204-011 Manual examination table w/o top 002-10143-XXX² 204 Seamless UltraFree 28" upholstery top 204-012 Manual examination table w/ drawer heater, pelvic tilt w/o top 002-10137-XXX² 204 Soft Touch 28" upholstery top 204-012 Manual examination table w/ drawer heater, pelvic tilt w/o top 002-10138-XXX² 204 Upholstery top - Styled, Seamless 204-012 Manual examination table w/ drawer heater, pelvic tilt w/o top 002-10143-XXX² 204 Seamless UltraFree 28" upholstery top

Accessories Description

Arm/Handle Systems

9A627001 ¹ 9A632001 ¹ 9A632002 9A632003	Adjustable arm system (arms only) Assist handles (handles only), field installed Accessory receivers only, plant installed Accessory receivers only, field installed
Knee Crutches 9A410001	Standard knee crutch only field installed
9A410004	Standard knee crutch kit with accessory receiver field installed
9A411001	Articulating knee crutch only field installed
9A411004	Articulating knee crutch kit with accessory receiver field installed
9A412003	Mounting brackets, hardware with accessory receiver field installed
Miscellaneous	
9A372001 9A70000 002-10141-XXX ² 002-10145-XXX ² 002-10193-00-XXX ²	Pass-through drawer stop (set of 3) Stainless steel treatment pan Flat Stitched Upholstery w/o pillow Flat UltraFree Upholstery top w/o pillow Thick Footrest Pad kit (for Soft-Touch and Flat uph)
002 10170 007000	iniekt ootioset aa kie (iot oote toden and that aprij

Thick Footrest Pad kit (for Seamless/UltraFree uph)

¹ - 9A632002 or 9A632003 accessory receiver required

² - See color palette on page 19

002-10194-00-XXX²

27



Ritter® 203

Treatment Table

- 400 lb patient weight capacity
- 27" W x 74" L x 31" H
- Please specify upholstery/panel color
- Panels match upholstery color
- Table weight: 135 lb
- Shipping weight: 170 lb

Model

Description

203-011-XXX^{1,2} 203-012-XXX^{1,2} **Miscellaneous** 203-011-999 203-012-999 Manual treatment table standard Manual treatment table w/drawers

Manual treatment table standard - special color Manual treatment table w/drawers - special color

¹ - Additional charge for UltraFree Upholstery/Accessory

² - See color palette on page 19



Model

95-001-XXX^{1,2} 95-001-999

Accessories 9A364001-XXX^{1,2} 9A364001-999

Ritter[®] 95

- 350 lb patient weight capacity
- 27" W x 72" L x 30" H
- Please specify upholstery color
- Frame available in a dark oak finish only
- Paper roll holder included
- Table weight: 80 lb
- Shipping weight: 100 lb

Description

Manual treatment table w/o shelf Manual treatment table w/o shelf - special color

Description

Shelf (field installed only) Shelf special color

¹ - Additional charge for UltraFree Upholstery/Accessory

² - See color palette on page 19

Medical Guide Book Effective January 1, 2022



Midmark[®] 425

Physician Seating (Hand Operated)

- Weight capacity: 225 lb •
- Stool weight: 20 lb ٠
- Shipping weight: 25 lb •

Model

425-001 002-0273-XXX^{1,2} Basic stool w/o seat cushion

Description

Seat cushion

Accessories

002-0274-XXX^{1,2} 002-0275-XXX^{1,2} 002-1546-00

Back kit Arm kit Lab height w/ black foot ring (425 only)

¹ - See color palette on page 19

² - Additional charge for Special or UltraFree Upholstery/Accessory Color



Midmark® 427 Physician Seating (Foot Operated)

- Weight capacity: 225 lb
- Stool weight: 20 lb ٠
- Shipping weight: 25 lb

Description

427-001 002-0273-XXX^{1,2}

Model

Basic stool w/o seat cushion Seat cushion

Back kit

¹ - See color palette on page 19

² - Additional charge for Special or UltraFree Upholstery/Accessory color

Midmark Chairs + Tables

BACK TO CONTENTS

29

Accessories

002-0274-XXX^{1,2} 002-0275-XXX^{1,2}

Arm kit



Physician Seating (Hand Operated)

- Chrome five leg base
- Weight capacity: 225 lb ٠
- Stool weight: 30 lb
- Shipping weight: 35 lb

Description

Model

279-001-XXX^{1,2}

Basic stool w/hand release - pneumatically adjustable



Ritter[®] 278 Physician Seating (Foot Operated)

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 30 lb
- Shipping weight: 35 lb

Model

278-001-XXX^{1,2}

Basic stool w/foot release - pneumatically adjustable



Model

277-001-XXX^{1,2} 277-201-XXX^{1,2}

Ritter[®] 277 Physician Seating with Backrest

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Basic stool w/back - pneumatically adjustable, hand release Basic stool w/back - pneumatically adjustable, foot release

hand release

Model 276-001-XXX^{1,2} 276-201-XXX^{1,2}

Ritter[®] 276 Physician Seating

- Chrome five leg base
- Weight capacity: 350 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description

Basic stool - pneumatically adjustable, hand release Basic stool - pneumatically adjustable, foot release



foot release

¹ - See color palette on page 19

² - Additional charge for Special or UltraFree Upholstery color



Description





Model 275-001-XXX^{1,2}

Physician Seating with Backrest

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Description

Basic stool w/back - manually adjustable

Ritter[®] 274

Physician Seating Chrome five leg base Weight capacity: 350 lb Stool weight: 15 lb

Shipping weight: 20 lb

Basic stool - manually adjustable

Ritter[®] 273

Weight capacity: 225 lb Stool weight: 20 lb Shipping weight: 25 lb

Black composite five leg base



Model 274-001-XXX^{1,2}



Model 273-001-XXX^{1,2}

Model

272-001-XXX1,2



Description Basic stool w/back - pneumatically adjustable

Physician Seating with Backrest

Ritter[®] 272 Physician Seating

- Black composite five leg base •
- Weight capacity: 300 lb •
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description

Basic stool - pneumatically adjustable

- ¹ See color palette on page 19
- ²- Additional charge for Special or UltraFree Upholstery color

Midmark Chairs + Tables

31



Physician Seating with Backrest

- Black composite five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

Description

Basic stool w/back - manually adjustable

Ritter[®] 270

Physician Seating

- Black composite five leg base
- Weight capacity: 300 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

Description

Basic Stool - manually adjustable

Model	Description
Miscellaneous	
9A366001-XXX ^{1,2}	Back rest kit (weight capacity, 225 lb)
27X-002-XXX ^{1,2}	Set of soft rubber casters
27X-003-XXX ^{1,2}	Set of glides (not available on foot release models)
27X-004-XXX ^{1,2}	Set of locking casters
27X-005-XXX ^{1,2}	Auto locking casters
9A367001	Foot ring, field installed (not available on foot
	release models)
9A368001	6" lab height ext. kit (not available on foot release models)
	Telease models/

¹ - See color palette on page 19

² - Additional charge for Special or UltraFree Upholstery/Accessory color



Model 270-001-XXX^{1,2}



- 400 lb patient weight capacity
- Wall-saver, kickback rear leg design
- Leveling feet
- Chair weight: 25 lb
- Shipping weight: 35 lb

Description

Basic chair (without arms)

Basic chair (with arms)

280-001-XXX¹ 280-001-XXX¹ 280-002-XXX¹ 280-002-XXX¹

Model

Miscellaneous

280-001-999 280-002-999 Basic chair Special Color (without arms)

Basic chair Special Color (with arms)

Basic chair (with arms) UltraFree upholstery

Basic chair (without arms) UltraFree upholstery

¹ - See color palette on page 19



Model

Description

282-001-XXX¹ 282-001-XXX¹ Basic chair Basic chair UltraFree upholstery

Miscellaneous 282-001-999

Basic chair Special Color

¹ - See color palette on page 19

Ritter[®] 282

- 650 lb patient weight capacity
- Wall-saver, kickback rear leg design
- Leveling feet
- Chair weight: 37 lb
- Shipping weight: 47 lb

33



Ritter® 281 Blood Draw Chair

- 400 lb patient weight capacity
- Optional storage drawer
- Chair weight: 60 lb
- Shipping weight: 75 lb



Model

Description

281-011-XXX ¹	Basic chair (without storage drawer)
281-011-XXX ¹	Basic chair (without storage drawer) UltraFree upholstery
281-012-XXX ¹	Basic chair (with storage drawer)
281-012-XXX ¹	Basic chair (with storage drawer) UltraFree upholstery
Miscellaneous	

281-011-999 281-012-999

Basic chair (without storage drawer) Basic chair (with storage drawer)

¹ - See color palette on page 19

		 4,400 K color temperature Color Rendering Index (CRI) greater than 95 7,500 fc intensity (80,700 lux) at 1 meter (39") focal length Effortless positioning Weight: 65 lb (118 lb for dual light) Shipping weight: 120 lb (190 lb for dual light) 	
Mobile Procedure Light	Single Procedure Light	Dual Procedure Light	Wall Mount
Model	Description		
255-001 255-002 255-003 255-004 255-005	Single mount, 8' ceiling Single mount, 9' ceiling Dual mount, 8' ceiling Dual mount, 9' ceiling Mobile Procedure Light		
255-006	Wall Mount		
Accessories 002-1865-00 ^{2, 3} 9A584001 ¹	10' ceiling mount adapter Conversion kit w/cover (si		
 To convert an existing sing Converts an 8' or 9' light t One required per light arr 	o a 10' ceiling height	ersion kit must be ordered in additior	n to a second 8' or 9' light
The following model numbe bottom of the ceiling surface		ne following ceiling heights (measure	d from the
Less than 8' 8' to 9'	255-005 255-001 255-003	Mobile Procedure Light 8' Single or Dual Procedure Light	

Ritter[®] 255 LED Procedure Light

Les 8' to 9' 9' to 10' Greater than 10'

Medical Guide Book Effective January 1, 2022

> 255-001, 255-003 255-002, 255-004 255-001, 255-003, 255-005

8' Single or Dual Procedure Light 9' Single or Dual Procedure Light 8' Single or Dual Light with 10' ceiling mount adapter kit(s) or Mobile Procedure Light

BACK TO CONTENTS



- LED technology for up to 50,000 hours of useful life
- Designed to consume 45-60% less energy than other lights
- 4,500 K color temperature
- 1,550 fc intensity (16,500 lux) at 16" focal length
- Sleek, ergonomic and easy-to-clean design
- Mist color
- Weight of light unit: 6 lb (253-002 light kit including bracket)
- Shipping weight: 8 lb (253-002 light kit including bracket)



Exam Chair Kit

Mobile Caster Stand Kit

Model	Description		
253-011	LED Exam Light - light only - excludes mounting hardware		
253-012 ³	LED Exam Light - exam chair kit - models 204/22X*/62X* - includes light, splines and chair mounting hardware, field installed		
253-013 ²	LED Exam Light - exam chair kit - model 626 for 32" upholstery and Patient Support Rails - includes light, splines and chair mounting hardware, field installed		
253-015	LED Exam Light - procedure chair kit - models 230/641/646/647 - includes light, base rail and chair mounting hardware, field installed		
253-016	LED Exam Light - mobile caster stand kit - includes light and caster stand		
253-017	LED Exam Light - wall mount kit - includes light, wall mounting bracket and hardware, field installed		
Miscellaneous			
9A619001	LED exam chair mounting hardware (only) - models 204/22X*/604/62X* - field installed		
9A620001	LED exam chair mounting hardware (only) - model 626 w/ 32" upholstery - field installed		
9A622001	LED procedure chair mounting hardware (only) - models 230/641/646/647 - field installed		
9A624001	LED 630 procedure chair spline mounting hardware (only), 28" wide upholstery, field installed		
9A625001	LED 630 procedure chair spline mounting hardware (only), 32" wide upholstery, field installed		
9A6260011	LED 630 procedure chair rail mount hardware (only), field installed		
9A459001	Disposable light sheaths; 5 boxes of 25 per carton		
9A628001	LED light mobile caster stand (only)		
9A623001	LED wall mount bracket (only)		
* Mounta to abairs with a 2004 or newer manufacture data			

 * Mounts to chairs with a 2006 or newer manufacture date

¹ - Requires (1) accessory rail 9A557001 or 9A557002

² - Use 253-013 with Patient Support Rails

³ - Slight interference with Patient Support Rails, use 253-013 to minimize **36**



Ritter[®] 250 LED Exam Light

- LED technology for up to 25,000 hours of useful life
- Designed to consume 45-60% less energy than other lights
- 4,000 K color temperature
- 1,800 fc intensity (19,375 lux) at 18" focal length
- Mist color
- Exam light weight: 20 lb
- Shipping weight: 25 lb



Caster Base Accessory

Model

250-003 9A621001

Description

LED exam light Caster base accessory



Midmark[®] Diagnostic Medical Products (ECG, Spirometer and Vital Signs) Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will repair or replace devices manufactured by Midmark (except for those not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of applicable devices. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective. This warranty extends only to the original retail purchaser of a product, and is not transferable or assignable.

Midmark warrants to the original retail purchaser that it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications. Midmark does not warrant that the software: (1) is error free; (2) can be used without problems or interruptions; or (3) is free from vulnerability to intrusion or attack by viruses or other methods. Midmark also agrees to provide the original retail purchaser with software updates, if any, that become generally available to the public within the warranty period stated, for the specific product owned by the purchaser.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original retail purchaser, shall be as follows:

- (a) One (1) year for software and diagnostic devices (unless otherwise noted below and except items listed under "Exclusions"),
- (b) Two (2) years for vital signs devices,
- (c) One (1) year for reusable blood pressure cuffs and SpO_2 sensors,
- (d) Ninety (90) days for all other accessories (sold in a product kit or purchased separately),

OBTAINING WARRANTY SERVICE

Warranty service must be obtained through Midmark. Contact Midmark Technical Service for warranty service inquiries or issues by phone at 1.844.856.1230, option 2, or via email at: techsupport@midmark.com

It is the original retail purchaser's obligation to arrange for delivery of a product to Midmark for warranty service or exchange. It is also the original retail purchaser's obligation to comply with the warranty service instructions provided by Midmark.

EXCLUSIONS

- This warranty does not cover and Midmark shall not be liable for the following:
- (1) any product that has been repaired, modified, or tampered with by any person or entity other than Midmark;
- (2) any product for which any label or tag attached to the product has been removed or tampered with;
- (3) defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering, or failure to seek and obtain repair or replacement in a timely manner;
- (4) products which are not installed, used, and properly cleaned and maintained as stated in the Midmark Installation, User Guide and/or Operation Manual for the applicable product;
- (5) products considered to be of a consumable nature;
- (6) accessories or parts covered by the respective original manufacturers' warranty, including but not limited to treadmills, computers, printers, scales, carts, mounting options and carrying cases, after the initial 90 days;
- (7) charges by anyone for adjustments, repairs, replacement parts, installation, or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (8) costs and expenses of routine maintenance and cleaning;
- (9) representations and warranties made by any person or entity other than Midmark; and
- (10) items provided through sales promotions.

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF DATA, LOSS OF USE, DOWNTIME, COVER, AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS, AND BENEFITS.

NO AUTHORIZATION

No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

WARRANTY DISCLAIMER

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

THIS WARRANTY WILL BE VOID IF ANY COVERED PRODUCT IS REPAIRED, MODIFIED, OR TAMPERED WITH BY ANY PERSON OR ENTITY OTHER THAN MIDMARK, OR IF ANY LABEL OR TAG ATTACHED TO A COVERED PRODUCT IS REMOVED OR TEMPERED WITH.

STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period. In the event multiple warranty periods exist with respect to a product, the ninety (90) day period provided for herein shall begin to run from expiration of the warranty period for the component to which the claim relates.

38

39

Midmark[®] Digital Diagnostic Devices Part Number Ordering Information

Note: Midmark recommends that all Midmark devices drop-ship directly to the customer.

TO ORDER:

Midmark devices require either an electronic medical record (EMR) that has integrated Midmark software or Midmark IQmanager[®] Software to operate. Devices include, but are not limited to, the part numbers shown below:

4-000-0027 4-000-0028 4-000-0724 4-000-0722 4-000-0720 4-000-0624 4-000-0622 4-000-0620

If the customer *currently has* an EMR software package that has integrated Midmark software, no additional software is required.

Place an order for the appropriate device(s) and please provide:

- EMR name and version
- Customer facility and contact name, address and phone
- Customer computer environment (Thin Client, Fat Client)

If the customer *does not have* an EMR, or **is** using Allscripts TouchWorks[™] EHR, Epic[®] EMR, AthenaHealth[®], AthenaPractice/ AthenaFlow or Cerner[®], please contact your local Midmark Representative to order the appropriate Midmark software solution. Additional charges or license fees will apply.

For ordering questions regarding Midmark devices, please contact your local Midmark Representative or call 844.856.1230, option 2.

Midmark HEINE® Diagnostic Instruments Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace diagnostic instruments manufactured by HEINE and sold by Midmark (except for those not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of invoice to the original retail purchaser of the product shall be as follows: (1) Ten (10) years for diagnostic instruments (other than GAMMA G5 Sphygmomanometers) manufactured by HEINE and sold by Midmark; (2) Five (5) years for GAMMA G5 Sphygmomanometers manufactured by HEINE and sold by Midmark; and (3) Twelve (12) months on parts and labor.

OBTAINING WARRANTY SERVICE

Warranty service must be obtained through Midmark. Midmark may be contacted for warranty service inquiries or issues via email at midmark.com, by phone at 1.844.856.1230, option 3 or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380. It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

EXCLUSIONS:

This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/
- or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Normal wear and tear
- (4) Fire or theft
- (5) Products considered to be single-use items or consumables (e.g. bulbs, cuffs and batteries)
- (6) Rechargeable batteries and batteries
- (7) Use of non-original HEINE spare parts
- (8) Exposed and accessible control and supply cables with corresponding plug connection
- (9) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (10) Costs and expenses of routine maintenance and cleaning;
- (11) Representations and warranties made by any person or entity other than Midmark.
- (12) Custom manufactured products;
- (13) Alterations or modifications to the product by any person or entity other than Midmark; and
- (14) Products that would otherwise by covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market

EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN.

WARRANTY DISCLAIMER

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

NO AUTHORIZATION

No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

*Warranty information included in product manuals

HEINE is a registered trademark of HEINE Optotechnik GmbH + Co. KG.

40

Synthesis® Cabinetry

Midmark® Medical Diagnostic Products Support Information

For technical assistance with a Midmark medical diagnostic product, please contact Technical Service at 844.856.1230.

Select: Option 2 – for Digital Vital Signs Device, IQvitals® Zone[™], Digital Spirometer, IQecg[®], Software, Supplies and Midmark IQcart[®] Diagnostic Equipment Cart. Available Monday-Friday, 8:00 am – 8:00 pm, ET.

Select: Option 3 – for Diagnostic Instrument Wall/Desktop Kits, Instrument Replacement Parts and Supplies (Otoscope/Ophthalmoscope) and Sphygmomanometer Kits. Available 8:00 am – 7:00 pm Monday-Thursday and 8:00 am – 5:00 pm Friday, ET.

For troubleshooting tips and product use information, refer to the operation manual. Additional help topics, self-guided user tools and documentation can be accessed at any time through our support section and knowledge base at midmark.com.

Return Material Authorization

To return any product for repair, an Return Material Authorization (RMA) number must be obtained from a Midmark Technical Service Specialist. The RMA number should be referenced on the package(s) containing the items to be returned.

Shipping

Before shipping any product to Midmark, be certain that an RMA number has been issued and that all guidelines regarding this authorization are followed. We highly recommend that you follow all guidelines for the shipment of medical products set forth by the shipping company you choose to use. If a question should arise regarding the appropriate method of shipment, please feel free to ask when calling for your RMA number. It is the responsibility of the customer when shipping a product to ensure that all packages and their contents get to Midmark safely.

Midmark will not assume responsibility for damage incurred by improper packaging, shipment or product use. Such actions will void all applicable warranties.

41



Shown with IQvitals[®] Zone[™] Mobile Cart Option

Mobile cart not included in pricing

Midmark IQvitals® Zone™ Vital Signs Wireless Monitor

with Touchscreen, Serenity[™] BP and Exergen TemporalScanner[™] - Includes Spot, Averaging, Continuous, Interval, SPRINT and Custom BP Modes

- 1 IQvitals Zone Wireless Monitor
- 1 Reusable Dual Tube BP Cuff Adult (26-35cm)
- 1 Reusable Dual Tube BP Cuff Large Adult (32-42cm)
- 1 Dual Tube BP Hose 6.5'
- 1 Exergen TemporalScanner Thermometer
- 1 Exergen Disposable Caps (10 ea)
- 1 IQvitals Zone Wireless USB Dongle
- 1 IQvitals Zone USB Cable 10'
- 1 IQvitals Zone Laptop Tray 15"
- 1 Front Storage Bin Kit
- 1 Device to Temperature Module Cable
- 1 Three-Arm Knobs (Set of 2)
- 1 AC Power Supply
- 1 Setup Guide
- 1 Quick Reference Guide
- 1 Operation Manual CD

Model	Description
4-000-0724	IQvitals Zone, Wireless Monitor, with Serenity BP, Exergen TemporalScanner and Masimo $^{\circledast}\text{SpO}_2$
4-000-0722	IQvitals Zone, Wireless Monitor, with Serenity BP, Exergen TemporalScanner and Nellcor [™] SpO_2
4-000-0720	IQvitals Zone, Wireless Monitor, with Serenity BP and Exergen TemporalScanner

Exergen TemporalScanner (AHB™) Thermometer is a trademark of Exergen.

Masimo is a registered trademark of Masimo Corporation.

Nellcor is a trademark of Covidien AG.



Shown with IQvitals[®] Zone[™] Mobile Cart Option

Mobile cart not included in pricing

Midmark IQvitals® Zone™ Non-Wireless Vital Signs Device

with Touchscreen, Serenity[™] BP and Exergen TemporalScanner[™] - Includes Spot, Averaging, SPRINT and Custom BP Modes

- 1 IQvitals Zone Vital Signs Device
- 1 Reusable Dual Tube BP Cuff Adult (26-35cm)
- 1 Reusable Dual Tube BP Cuff Large Adult (32-42cm)
- 1 Dual Tube BP Hose 6.5'
- 1 Exergen TemporalScanner Thermometer
- 1 Exergen Disposable Caps (10 ea)
- 1 IQvitals Zone USB Cable 10'
- 1 IQvitals Zone Laptop Tray 15"
- 1 Front Storage Bin Kit
- 1 Device to Temperature Module Cable
- 1 Three-Arm Knobs (Set of 2)
- 1 AC Power Supply
- 1 Setup Guide
- 1 Quick Reference Guide
- 1 Operation Manual CD

Model	Description
4-000-0624	IQvitals Zone, non-Wireless, with Serenity BP, Exergen TemporalScanner and Masimo® SpO ₂

- 4-000-0620 IQvitals Zone, non-Wireless, with Serenity BP and Exergen TemporalScanner

Exergen TemporalScanner (AHB™) Thermometer is a trademark of Exergen.

Masimo is a registered trademark of Masimo Corporation.

Nellcor is a trademark of Covidien AG.

43

IQvitals[®] Zone[™] Mounting Options

Model Description

6216-001-813 Vital Signs Monitor Mobile Cart 2-200-0082 IQvitals Zone Wall Mount Articulating Arm

IQvitals Zone Accessories + Supplies

Model Description

1-100-1603	Fairbanks® TeleWeigh™ Digital Scale
3-200-0030	Blood Pressure Cuff, Dual Tube, Reusable, Child (13-20cm)
3-200-0031	Blood Pressure Cuff, Dual Tube, Reusable, Small Adult (18-26cm)
3-200-0032	Blood Pressure Cuff, Dual Tube, Reusable, Adult (26-35cm)
3-200-0033	Blood Pressure Cuff, Dual Tube, Reusable, Adult Long (29-38cm)
3-200-0034	Blood Pressure Cuff, Dual Tube, Reusable, Large Adult (32-42cm)
3-200-0035	Blood Pressure Cuff, Dual Tube, Reusable, Large Adult Long (35-44cm)
3-200-0036	Blood Pressure Cuff, Dual Tube, Reusable, Thigh (42-50cm)
2-200-0020	Blood Pressure Cuff, Dual Tube, Disposable, Child (13-20cm), 5/pk
2-200-0021	Blood Pressure Cuff, Dual Tube, Disposable, Small Adult (18-26cm), 5/pk
2-200-0022	Blood Pressure Cuff, Dual Tube, Disposable, Adult (26-35cm), 5/pk
2-200-0023	Blood Pressure Cuff, Dual Tube, Disposable, Adult Long (29-38cm) , 5/pk
2-200-0024	Blood Pressure Cuff, Dual Tube, Disposable, Large Adult (32-42cm), 5/pk
2-200-0025	Blood Pressure Cuff, Dual Tube, Disposable, Large Adult Long (35-44cm), 5/pk
2-200-0026	Blood Pressure Cuff, Dual Tube, Disposable, Thigh (42-50cm), 5/pk
2-200-0041	SpO ₂ Sensor, Reusable, Masimo [®] , Pediatric (RD Set DCI-P)
2-200-0040	SpO_2^{-} Sensor, Reusable, Masimo, Adult (RD SET DCI)
2-200-0042	SpO ₂ Patient Cable, Masimo (RD SET D20-05) - 4'
3-200-0007	SpO ₂ Sensor Kit, Reusable, Nellcor™, Pediatric (D-YS and D-YSPD)
2-200-0233	SpO_2 Sensor, Reusable, Nellcor, Adult (DS100A)
2-200-0234	SpO ₂ Interface Cable, Nellcor (DOC4) - 4'
2-200-0115	Disposable Caps, Exergen Thermometer, 1000pc/box
2-200-0116	Disposable Caps, Exergen Thermometer, 100pc/box
2-200-0074	Storage Bin Kit, Dual, Cart Mounted
2-200-0084	Storage Bin Kit, Wall Mounted
3-200-0009	Tray, IQvitals Zone - 12"
3-200-0012	USB Dongle, IQvitals Zone Wireless
9A654001	Mobile Cart Scale Holder

Connectivity Supplies for Midmark® 626 Exam Chair

Model Description

9A478001	Digital Scale serial cable for IQvitals [®] , 6' length, coiled
9A478002	Digital Scale serial cable for IQvitals, 15' length, straight
9A478003	Digital Scale serial cable for IQvitals, 30' length, straight
9A478004	Digital Scale serial cable for IQvitals, 50' length, straight
002-10064-00	626 Chair Connectivity Kit

Fairbanks and TeleWeigh are trademarks of Fairbanks Scales, Inc. Masimo is a registered trademark of Masimo Corporation. Nellcor is a trademark of Covidien AG.

44



Midmark[®] Digital Vital Signs Device

with Touchscreen, Blood Pressure + Temperature

- 1 Digital Vital Signs Device w/ Blood Pressure and Temperature
- 1 AC Power Supply
- 1 Reusable Blood Pressure Cuff Adult
- 1 Reusable Blood Pressure Cuff Large Adult
- 1 Blood Pressure Hose 6.5'
- 1 Oral/Axillary Temperature Probe
- 1 Temperature Probe Covers (box of 20)
- 1 USB Cable 10'
- 1 Operation Manual CD
- 1 Quick Reference Guide

Model Description

4-000-0500*	Digital Vital Signs Device with Blood Pressure
	and Alaris [®] Turbo Temp [®] Thermometer

Digital Vital Signs Device with Touchscreen, Blood Pressure, Temperature + SpO₂

- 1 Digital Vital Signs Device w/Blood Pressure, Temperature and SpO₂
- 1 AC Power Supply
- 1 Reusable Blood Pressure Cuff Adult
- 1 Reusable Blood Pressure Cuff Large Adult
- 1 Blood Pressure Hose 6.5'
- 1 Reusable SpO₂ Adult Finger Sensor
- 1 SpO₂ Extender Cable 4'
- 1 Oral/Axillary Temperature Probe
- 1 Temperature Probe Covers (box of 20)
- 1 USB Cable 10'
- 1 Operation Manual CD
- 1 Quick Reference Guide

Model Description

4-000-0510* Digital Vital Signs Device with Blood Pressure, Alaris Turbo Temp Thermometer and SpO₂

Alaris and Turbo Temp are registered trademarks of Cardinal Health, Inc.

45

Midmark[®] Digital Vital Signs Device Accessories + Supplies

Model

Description

1-100-1603	Fairbanks® TeleWeigh™ Digital Scale
1-100-1605	Thermal Printer for Digital Vital Signs Device
3-004-2004	Mobile Cart Printer Bracket for Digital Vital Signs Device
3-004-2010	Mobile Cart Scale Mount for Digital Vital Signs Device
3-009-0010	AC Power Supply
3-009-0012	Power Supply Adapter Kit
3-009-0014	Replacement Lithium Ion Battery for Digital Vital Signs Device
3-009-0016	USB Cable - 10'
3-009-0022	Blood Pressure Hose - 6.5'
3-009-0100	Blood Pressure Hose - 10'
3-009-0021	Reusable SpO ₂ Small Finger Sensor, Envitec
3-009-0020	Reusable SpO ₂ Adult Finger Sensor, Envitec
3-009-0026	SpO ₂ Extender Cable, 4', Envitec
3-009-0024	Alaris [®] Turbo Temp [®] Oral/Axillary Temperature Probe
3-009-0058	Turbo Temp Probe Covers, Carton of 10 boxes
	(20 covers per box)
3-009-0068	Reusable Blood Pressure Cuff, Infant (8-14 cm)
3-009-0070	Reusable Blood Pressure Cuff, Child (13-20 cm)
3-009-0062	Reusable Blood Pressure Cuff, Small Adult (18-26 cm)
3-009-0064	Reusable Blood Pressure Cuff, Adult, Midmark (26-35 cm)
3-009-0066	Reusable Blood Pressure Cuff, Large Adult, Midmark (32-42 cm)
3-009-0072	Reusable Blood Pressure Cuff, Adult Long (29-38 cm)
3-009-0074	Reusable Blood Pressure Cuff, Large Adult Long (35-44 cm)
3-009-0076	Reusable Blood Pressure Cuff, Thigh (42-50 cm)
3-004-2000	Mobile Cart for Digital Vital Signs Device
3-009-0001	Countertop Mount for Digital Vital Signs Device
3-009-0003	Wall Mount with Basket for Digital Vital Signs Device
002-10592-00	Equipment Pole Mount for Digital Vital Signs Device

Fairbanks and TeleWeigh are trademarks of Fairbanks Scales, Inc.

Alaris and Turbo Temp are registered trademarks of Cardinal Health, Inc.

Connectivity Supplies for the Midmark[®] 626 Exam Chair

Model	Description
9A478001	Digital Scale serial cable for Digital Vital Signs Device, 6' length, coiled
9A478002	Digital Scale serial cable for Digital Vital Signs Device, 15' length, straight
9A478003	Digital Scale serial cable for Digital Vital Signs Device, 30' length, straight
9A478004	Digital Scale serial cable for Digital Vital Signs Device, 50' length, straight
002-10064-00	626 Chair Connectivity Kit*

*Not recommended for connecting Midmark® Digital Vital Signs Device to PC running Windows® 10 64-bit OS

Windows is a registered trademark of Microsoft Corporation.

Ritter Chairs + Tables

47



Midmark[®] Digital Spirometer

- 1 Digital Spirometer
- 10 Disposable Spirometer Mouthpieces
- 1 10/pack Disposable Nose Clips
- 1 Mouse Pad
- 1 Operation Manual CD
- 1 Quick Reference Guide
- 1 Carrying Case

Model

Description

4-000-0028	Digital Spirometer - USB version with syringe
4-000-0027	Digital Spirometer - USB version

Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

Midmark Digital Spirometer Accessories + Supplies

Model

- Description
- 1-100-0009 3-Liter Spirometer Calibration Syringe 2-100-1205 Disposable Spirometer Mouthpieces (box of 25) 2-100-1206
 - Disposable Spirometer Mouthpieces (case 4 boxes)
- 2-100-1207 Box of 500 Mouthpieces, individually bagged, bulk
- 2-244-0001 Disposable Nose Clips (box of 25)



Midmark IQecg® Digital ECG with Lead Management (Interpretive)

1 - IQeca Acquisition Module

- 1 10-Lead Patient Cable with Lead Management
 - 1 10/pack Clear ECG Clips
 - 1 100/pack Disposable ECG Electrodes
 - 1 Mouse Pad
 - 1 Operation Manual CD
 - 1 Quick Reference Guide
 - 1 Carrying Case

Model Description

4-000-0062

IQecg device - with Lead Management

Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

IQecg Supplies

Model Description

- 3-047-0001* Universal ECG Clips (3 mm, 4 mm + Snap) 10/pack
- 3-047-0005 Clear ECG Clips (4 mm + Snap) 10/pack
- 3-100-0203 Patient Cable, with Lead Management, IQecg
- 2-100-0208 Disposable ECG Electrodes (box of 1,000)
- 2-100-0209 Disposable ECG Electrodes (case of 4 boxes - 1,000/box)

* Available to the US market only

49



Midmark HEINE[®] Wall System Kits

Model	Description
A-095-12-930-166	EN200 Wall Transformer with BETA® 400 LED Otoscope and BETA® 200 LED Ophthalmoscope + GAMMA XXL Sphygmomanometer + Tip Dispenser
A-095-12-931-166	EN200 Wall Transformer with BETA 400 LED Otoscope and BETA 200 LED Ophthalmoscope + Tip Dispenser
A-095-12-932-166	EN200 Wall Transformer with K180 Otoscope and K180 Ophthalmoscope + GAMMA XXL Sphygmomanometer + Tip Dispenser
A-095-12-933-166	EN200 Wall Transformer with K180 Otoscope and K180 Ophthalmoscope + Tip Dispenser

Midmark HEINE Wall System Kit Accessories

Model	Description
X-095-12-006-166	Wallboard (for device combination of EN200 + Tip Dispenser – no instruments)
X-095-12-007-166	Wallboard (for device combination of EN200 + GAMMA XXL Sphygmomanometer + Tip Dispenser – no instruments)

HEINE and BETA are registered trademarks of HEINE Optotechnik GmbH + Co. KG.

Available to the US market only



Midmark HEINE® Diagnostic Instrument Wall Kits

- 1 EN200 Wall Transformer
- 1 Otoscope with polycarbonate, 11.5' (3.5 m) corded handle
- 1 Ophthalmoscope with polycarbonate, 11.5' (3.5 m) corded handle

Model	Description
A-095-12-208-166	EN200 Wall Transformer with BETA® 400 LED Otoscope and BETA® 200 LED Ophthalmoscope
A-095-12-209-166	EN200 Wall Transformer with BETA 200 LED Otoscope and BETA 200 LED Ophthalmoscope
A-095-12-204-166	EN200 Wall Transformer with BETA 400 XHL Otoscope and BETA 200 XHL Ophthalmoscope
A-095-12-205-166	EN200 Wall Transformer with BETA 200 XHL Otoscope and BETA 200 XHL Ophthalmoscope
A-095-12-210-166	EN200 Wall Transformer with K180 XHL Otoscope and K180 XHL Ophthalmoscope

Midmark HEINE Diagnostic Instrument Wall Kit Accessories

Model	Description
X-095-12-006-166	Wallboard (for device combination of EN200 + Tip Dispenser - no instruments)
X-095-12-007-166	Wallboard (for device combination of EN200 + Tip Dispenser + GAMMA XXL Sphygmomanometer - no instruments)
B-000-11-149-166	AllSpec Tip Dispenser - white

Midmark Chairs + Tables

51

HEINE and BETA are registered trademarks of HEINE Optotechnik GmbH + Co. KG.



Midmark HEINE[®] Diagnostic Instrument Desktop Kits

- 1 NT4 Desktop Charger
- 1 Otoscope with BETA®4 NT, Lithium-ion rechargeable handle
- 1 Ophthalmoscope with BETA4 NT, Lithium-ion rechargeable handle

Model	Description
A-853-24-420-166	NT4 Desktop Charger with BETA® 400 LED Otoscope and BETA® 200 LED Ophthalmoscope
A-832-24-420-166	NT4 Desktop Charger with BETA 200 LED Otoscope and BETA 200 LED Ophthalmoscope
A-853-23-420-166	NT4 Desktop Charger with BETA 400 XHL Otoscope and BETA 200 XHL Ophthalmoscope
A-832-23-420-166	NT4 Desktop Charger with BETA 200 XHL Otoscope and BETA 200 XHL Ophthalmoscope
A-879-23-420-166	NT4 Desktop Charger with K180 XHL Otoscope and K180 XHL Ophthalmoscope

Midmark HEINE Diagnostic Instrument Desktop Kit Accessories

Model	Description
1110 GGI	Description

B-000-11-149-166 AllSpec Tip Dispenser - white

HEINE and BETA are registered trademarks of HEINE Optotechnik GmbH + Co. KG.

52



Instrument Panel

- Instruments not included
- Shipping weight: 17 lb

Product Information
Step 1. Configuration
I1 Instrument Panel - 30" x 20"
I2 Instrument Panel - 34" x 20"
Step 2. Option Package
X No Lock
Step 3. Cabinetry Finish*
XXX Solid Color Cabinet Finish (Excludes PVC-Free)
XXX Woodgrain, Metallic + PVC-Free Cabinet Finish

Configuration Example: I2X.803

*See page 73 for Cabinetry Finish

53



Midmark HEINE® Diagnostic Instrument Replacement Parts

Model

Description

B-008-11-400-166	BETA® 400 LED Otoscope (instrument head only)
B-002-11-400-166	BETA 400 XHL Otoscope (instrument head only)
B-008-11-500-166	BETA 200 LED Otoscope (instrument head only)
B-002-11-500-166	BETA 200 XHL Otoscope (instrument head only)
B-002-11-550-166	K180 XHL Otoscope (instrument head only)
C-008-30-100-166	BETA 200 LED Ophthalmoscope (instrument head only)
C-002-30-100-166	BETA 200 XHL Ophthalmoscope (instrument head only)
C-002-30-207-166	K180 XHL Ophthalmoscope (instrument head only)
X-002-99-494-166	NT4 Desktop Charger (without instruments)
X-007-99-396-166	BETA4 NT Rechargeable Handle (with battery and handle bottom cap)
X-002-99-394-166	BETA 4 NT Handle Bottom Cap
X-095-12-250-166	EN200 Wall Transformer (without instruments)



Midmark HEINE Diagnostic Instrument Supplies

Model

Description

B-000-11-127-166	AllSpec Disposable Tips 4 mm (box of 1,000)
B-000-11-137-166	AllSpec Disposable Tips 4 mm (case of 10 boxes, 1,000/box)
B-000-11-128-166	AllSpec Disposable Tips 2.5 mm (box of 1,000)
B-000-11-138-166	AllSpec Disposable Tips 2.5 mm (case of 10 boxes, 1,000/box)
X-002-88-078-166	XHL Xenon Halogen 3.5 V Bulb for BETA 400, BETA 200 and K180 Otoscopes (package of 6)
X-002-88-070-166	XHL Xenon Halogen 3.5 V Bulb for BETA 200 Ophthalmoscope (package of 6)
X-002-88-086-166	XHL Xenon Halogen 3.5 V Bulb for K180 Ophthalmoscope (package of 6)
X-007-99-383-166	Rechargeable Lithium-ion Battery

54

Available to the US market only



Midmark HEINE® GAMMA XXL LF Sphygmomanometer Kit

- 1 GAMMA XXL LF Professional Sphygmomanometer
- 1 Twin-tube adult cuff (29–41 cm circumference)
- 1 Twin-tube adult large cuff (35–47 cm circumference)
- 1 Universal mounting plate

Model	Description
M-000-09-323-166	GAMMA XXL LF Sphygmomanometer with universal mounting plate and two twin-tube cuffs (adult and adult large)

Midmark HEINE GAMMA XXL Accessories

Model Description M-000-09-701-166 Twin-tube LF blood pressure cuff, Infant (6 x 28 cm; 8 - 13 cm circumference) Twin-tube LF blood pressure cuff, Child M-000-09-703-166 (10 x 34 cm; 13 - 20 cm circumference) M-000-09-705-166 Twin-tube LF blood pressure cuff, Adult Small (14 x 53 cm; 20 - 29 cm circumference) Twin-tube LF blood pressure cuff, Adult M-000-09-629-166 (14 x 58 cm; 29 - 41 cm circumference) M-000-09-631-166 Twin-tube LF blood pressure cuff, Adult Large (17 x 70 cm; 35 - 47 cm circumference) M-000-09-625-166 Twin-tube LF blood pressure cuff, Thigh $(20 \times 86 \text{ cm}; > 41 \text{ cm circumference})$

55

Chairs + Tables

Seating

Lighting

Diagnostics

Ritter



Midmark HEINE® GAMMA G5 Sphygmomanometer Kit

- 1 GAMMA G5 Standard Sphygmomanometer
- 1 Single-tube adult cuff (29 41 cm circ.)
- 1 Single-tube adult large cuff (35 47 cm circ.)

Model	Description
M-000-09-230-166	GAMMA G5 Sphygmomanometer with two single-tube

cuffs (adult and adult large)

Midmark HEINE GAMMA G5 Accessories

Model Description

M-000-09-700-166	Single-tube LF blood pressure cuff, Infant (6 x 28 cm; 8 - 13 cm circumference)
M-000-09-702-166	Single-tube LF blood pressure cuff, Child (10 x 34 cm; 13 - 20 cm circumference)
M-000-09-704-166	Single-tube LF blood pressure cuff, Adult Small (14 x 53 cm; 20 - 29 cm circumference)
M-000-09-619-166	Single-tube LF blood pressure cuff, Adult (14 x 58 cm; 29 - 41 cm circumference)
M-000-09-706-166	Single-tube LF blood pressure cuff, Adult Large (17 x 70 cm; 35 - 47 cm circumference)
M-000-09-615-166	Single-tube LF blood pressure cuff, Thigh (20 x 86 cm; > 41 cm circumference)

HEINE is a registered trademark of HEINE Optotechnik GmbH + Co. KG.

56



Midmark IQcart® Diagnostic Equipment Cart

(Some assembly required)

- 1 IQcart with Arm
- 1 Storage Basket
- 1 IQecg® Arm
- 1 Digital Spirometer Holster
- 1 Set of Locking Disks



Description

3-004-1000

Midmark IQcart with Arm

IQcart Accessories + Supplies

Model

Description

002-10592-00Equipment Pole Mount for Digital Vital Signs Device3-004-6028IQcart Basket (8" x 6" x 8")

57

Midmark[®] Workstations: Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation warrants to the original purchaser its Midmark Workstation products manufactured by Midmark to be free from defects in material and workmanship under normal interior use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of the parts or the products the defects of which are reported to Midmark within the applicable warranty period and which upon examination by Midmark prove to be defective.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows:

- A. Five years for all Midmark Workstation products and components (except for the shorter periods in B), including drawer fronts and casters
- B. One year for Power system components and coiled power cords

EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for: (1) Failure to install or use according to the Midmark user and installation manuals (2) Normal wear and tear, repairs and replacements required because of misuse, abuse, negligence, alteration, accident, freight damage or tampering (including, without limitation, use of product in unsuitable environments or conditions); (3) matching of color, grain or texture except to commercially acceptable standards; (4) changes in color caused by natural or artificial light; (5) products which are not installed, used or properly cleaned as required in Midmark's written Installation Manual; (6) specially manufactured products; (7) products considered to be of a consumable nature such as light bulbs and surge suppression products; (8) accessories or parts not manufactured by Midmark or use of non-authorized Midmark components; (9) charges by anyone (including Midmark's authorized dealers) for adjustments, repairs, replacement parts, installation or other work preformed upon or in connection with such products which is not expressly authorized in writing in advance by Midmark. 10) Midmark Workstation products and accessories (pgs 56-67) follow warranty on page 5.

EXCLUSIVE REMEDY

Midmarks' only obligation under this warranty is the repair or replacement of defective parts. Midmark shall not be liable for any direct, special, indirect, incidental, exemplary or consequential damages or delay.

NO AUTHORIZATION

No person or firm is authorized to create for Midmark any other obligation or liability in connection with the products.

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.



Midmark[®] 6211 Secure Laptop Workstation

- 26" of height adjustment from 23" to 49"
- Keyed laptop security system
- 23" D x 26" W worksurface
- 11" D x 20.5" W secondary worksurface

Secure Laptop Workstation

Description

- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- 19.5" D x 20.5" W base dimension
- Weight capacity: 20 lb
- Unit weight: 58 lb
- Shipping weight: 103 lb

Model 6211-001-xxx¹

Accessories		
Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/ Steel Cable Mgmnt-LH
PMSR	9A655001	Power Mount w/ Steel Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

¹ - See color code chart on page 73

² - Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.

Lighting

59

Chairs + Tables

Seating

Ritter



Model 6212-001-xxx¹

Accessories

Midmark[®] 6212 Flat Panel PC Workstation

- 26" of height adjustment rom 23" to 49"
- 23" D x 26" W worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
 - Description

Flat Panel PC Workstation

- Monitor arm w/ 9.5" independent vertical adj. and up to 40° positive tilt
- Monitor arm capacity range 5 lb min 20 lb max
- Monitor arm independent rotation of 180°
- 22.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 55 lb
- Shipping weight: 100 lb

Accessories		
Plant Installed	Field Installed	
AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/ Steel Cable Mgmnt-LH
PMSR	9A655001	Power Mount w/ Steel Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Digital Spirometer Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

60

BACK TO CONTENTS

¹ - See color code chart on page 73
 ² - Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.

Seating

Workstations

Medical Guide Book Effective January 1, 2022



Model 6213-001-xxx¹

Accessories		
Plant Installed	Field Installed	
AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PCML	9A669001	MG Coiled Power Cord for Plastic
		Bin/Cord Storage -LH
PCMR	9A669001	MG Coiled Power Cord for Plastic
		Bin/Cord Storage -RH
PMBL	9A655002	Power Mount w/ Bin Cable Mgmnt-LH
PMBR	9A655002	Power Mount w/ Bin Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHEP	3-004-6030	Digital Spirometer Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

Midmark® 6213

• 26" of height adjustment from 23" to 49"

• 11" D x 20.5" W secondary worksurface

(two non-locking casters in the rear)

Flat Panel Secure PC Workstation

• Monitor arm w/ 75 mm and 100 mm

• Two total-locking casters in the front lock

both the wheel and swivel simultaneously

• Keyed PC security system

mounting patterns

Description

• 23" D x 26" W worksurface

Flat Panel Secure PC Workstation

¹ - See color code chart on page 73

² - Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.

- Monitor arm w/ 9.5" independent vertical adj. and up to 40° positive tilt
- Monitor arm capacity range 5 lb min 20 lb max
- Monitor arm independent rotation of 180°
- 22.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 64 lb
- Shipping weight: 109 lb

61

Medical Guide Book Effective January 1, 2022



Midmark® 6214 Procedure Workstation

- 26" of height adjustment from 23" to 49"
- 21" D x 26" W worksurface

• Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)

Procedure Workstation

- Includes: rear basket, BP cuff holder, integrated spirometer holsters and equipment pole
 19.5" D x 20.5" W(base dimension
- 19.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 56 lb
- Shipping weight: 101 lb

Model 6214-001-xxx¹

Accessories

/ (0000001100		
Plant Installed	Field Installed	
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
IQEP	002-10592-00	Equipment Pole Mount for Digital Vital Signs Device
PMSL	9A655001	Power Mount w/ Steel Cable Mgmnt-LH
PMSR	9A655001	Power Mount w/ Steel Cable Mgmnt-RH
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHEP	3-004-6030	Digital Spirometer Holder Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf

Description

¹ - See color code chart on page 73



Model 6215-001-xxx¹

Accessories

Accessories		
Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/ Steel Cable Mgmnt-LH
PMSR	9A655001	Power Mount w/ Steel Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Digital Spirometer Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

¹ - See color code chart on page 73

² - Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.

- Weight capacity: 30 lb
- Unit weight: 48 lb
- Shipping weight: 93 lb

Description

Compact Workstation

Compact Workstation

19.5" D x 20.5" W base dimension
Two total-locking casters in the front lock both the wheel and swivel simultaneously

(two non-locking casters in the rear)

• 21" D x 26" W worksurface

• 26" of height adjustment from 23" to 49"

Midmark® 6215-001

Midmark Chairs + Tables

63



26" of height adjustment from 23" to 49"
17.5" D x 23.5" W worksurface

Midmark® 6215-002 Compact Rectangular Workstation

- 19.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock

both the wheel and swivel simultaneously (two non-locking casters in the rear)

- Weight capacity: 20 lb
- Unit weight: 48 lb
- Shipping weight: 93 lb

Model 6215-002-xxx¹

Accessories

Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/ Steel Cable Mgmnt-LH
PMSR	9A655001	Power Mount w/ Steel Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Digital Spirometer Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

Description

Compact Workstation

¹ - See color code chart on page 73

² - Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.



Model 6215-003-xxx¹

Accessories

Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf

¹ - See color code chart on page 73 Requires 9A651001 or 9A658001 for field installed option.

Midmark® 6215-003

Compact Kidney Workstation

- 26" of height adjustment from 23" to 49"
- 17" D x 25" W worksurface

Description

Compact Workstation

- 19.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Weight capacity: 20 lb
- Unit weight: 48 lb
- Shipping weight: 93 lb

Midmark Chairs + Tables

65



Midmark® 6217 Secure Laptop Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- Keyed laptop security system
- 23" D x 22" W worksurface w/ slideout mouse surface

Description

• 11" D x 20.5" W secondary worksurface

Secure Laptop Workstation + Mouse Tray

- 19.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Weight capacity: 30 lb
- Unit weight: 60 lb
- Shipping weight: 105 lb

Model 6217-001-xxx¹

Accessories

Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/ Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/ Steel Cable Mgmnt-LH
PMSR	9A655001	Power Mount w/ Steel Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Digital Spirometer Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount

¹ - See color code chart on page 73

² - Requires EPWS or BPCH for plant installed option.

Requires 9A651001 or 9A658001 for field installed option.



Field Installed

9A664001

9A659001

N/A

Model 6218-001-xxx1

Accessories Plant Installed

AIOM

BCHD

BCHL

Midmark[®] 6218

Flat Panel PC Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- Monitor arm capacity range 5 lb min 20 lb max
- Monitor arm independent rotation of 180°
- 23" D x 22" W worksurface w/ slideout mouse surface
- 22.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)

Description

Monitor over 10 lb

Flat Panel PC Workstation + Mouse Tray

- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. and up to 40° positive
- Weight capacity: 30 lb
- Unit weight: 55 lb
- Shipping weight: 100 lb

Bar Code Holder Mount
Bar Code Holder-LH
Bar Code Holder-RH
Blood Pressure Cuff Holder
Cart Counterweight
CPU Holder-Back
CPU Holder-Front
Double Basket
Dual Monitor Mount
Drawer w/ Lock-Mouse Trk Mnt
Drawer-Mouse Track Mnt
Equipment Pole
Equipment Pole Mount for Digital Vital Signs Device
Laptop Display Holder
Mobile Printer Shelf
Power Mount w/ Steel Cable Mgmnt-LH
Power Mount w/ Steel Cable Mgmnt-RH
Rear Basket
Single Basket-Back
Single Basket-Front
Steel Brkt + MG Coiled Cord-LH
Steel Brkt + MG Coiled Cord-RH
Digital Spirometer Holder Equipment Pole Mount
Scale Holder-LH
Scale Holder-RH
Tablet Display Holder
Adjustable Utility Shelf
VESA Mounting Plate
VOIP Shelf-LH Mount

Midmark Chairs + Tables

67

9A659001 Bar Code Holder-**BCHR** Blood Pressure C BPCH 9A658001 CCWT 9A657001 Cart Counterweic CPUB 9A663001 CPU Holder-Back CPUF 9A663001 **CPU Holder-Front** DBCM 9A653001 Double Basket DMMT 9A667001 Dual Monitor Mo DWLM 9A647001 Drawer w/ Lock-N DWNM 9A646001 Drawer-Mouse Tra **EPWS** 9A651001 Equipment Pole 002-10592-00² IQEP² Equipment Pole N LPDH 9A662001 Laptop Display H MPSF 9A648001 Mobile Printer Sh Power Mount w/ S PMSL 9A655001 9A655001 Power Mount w/ S **PMSR** RBWS 9A650001 Rear Basket SBCB 9A652001 Single Basket-Bac SBCF 9A652001 Single Basket-Fro SCML 9A668001 Steel Brkt + MG C SCMR 9A668001 Steel Brkt + MG C SHEP 3-004-6030 **Digital Spiromete** Scale Holder-LH SHLH 9A654001 SHRH 9A654001 Scale Holder-RH TBDH 9A661001 Tablet Display Ho USCB Adjustable Utility 9A656001 VESA 9A660001 **VESA Mounting P** VIPL 9A649001 VOIP Shelf-LH Mo VIPR **VOIP Shelf-RH Mount** 9A649002 VPBL 9A695001 VOIP Steel Bracket-LH Mount

VOIP Steel Bracket-RH Mount

¹ - See color code chart on page 73

VPBR

² - Requires EPWS or BPCH for plant installed option.

9A695002

Requires 9A651001 or 9A658001 for field installed option.



Midmark® 6219 Flat Panel Secure PC Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- Monitor arm capacity range 5 lb min 20 lb max
- Keyed PC security system
- Monitor arm independent rotation of 180°
- 23" D x 22" W worksurface w/ slideout mouse surface
- 22.5" D x 20.5" W base dimension
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two nonlocking casters in the rear)
 - Description

Flat Panel Secure PC Workstation + Mouse Tray

Accessories

6219-001-xxx¹

Model

710003301103		
Plant Installed	Field Installed	
AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLM	9A647001	Drawer w/ Lock-Mouse Trk Mnt
DWNM	9A646001	Drawer-Mouse Track Mnt
EPWS	9A651001	Equipment Pole
IQEP ²	002-10592-00 ²	Equipment Pole Mount for Digital Vital Signs Device
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PCML	9A669001	MG Coiled Power Cord for Plastic
		Bin/Cord Storage-LH
PCMR	9A669001	MG Coiled Power Cord for Plastic
		Bin/Cord Storage-RH
PMSL	9A655001	Power Mount w/ Steel Cable Mgmnt-LH
PMSR	9A655001	Power Mount w/ Steel Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHEP	3-004-6030	Digital Spirometer Holder Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount
VPBL	9A695001	VOIP Steel Bracket-LH Mount
VPBR	9A695002	VOIP Steel Bracket-RH Mount
		1 See color code chart on page 72

68

BACK TO CONTENTS

¹ - See color code chart on page 73
 ² - Requires EPWS or BPCH for plant installed option.
 Requires 9A651001 or 9A658001 for field installed option.

- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. + up to 40° positive tilt
- Weight capacity: 30 lb
- Unit weight: 64 lbShipping weight: 109 lb

Medical Guide Book Effective January 1, 2022



Midmark[®] 6281 Workstation

- 22" height adjustment
- 19.5" W x 14" D worksurface
- Slide-out mouse surface
- Up to 45" reach off wall
- Weight capacity: 0-8 lb

Description

Laptop Wall Mounted Workstation

Accessories

Model

6281-001

Plant Installed	Field Installed	
n/a	9A588001	20" Extension Arm
n/a	9A591001	Track Mount Bar Code Holder
n/a	9A590001	Tamper Proof PC Holder
n/a	9A589001	Adjustable PC Holder

- Folds to 9" flat to wall
- Unit weight: 33 lb
- Shipping weight: 40 lb
- Work surface accommodates up to 17" laptops

Ritter Chairs + Tables

69



Midmark[®] 6282 Workstation

- 22" worksurface height adjustment
- 25" W x 10" D worksurface
- Max monitor capacity is 20 lb
- 9.5" independent monitor height adjustment

• Independent 180° monitor rotation and 0-45° tilt

- VESA mount 75 mm and 100 mm
- Folds to 9.47" flat to wall
- Weight capacity: 5 20 lb
- Product weight: 43 lb
- Shipping weight: 50 lb

Model

6282-001

Description

Standard Duty Wall Mounted Workstation

Accessories

Field Installed	
9A593001	Dual Monitor Arm
9A588001	20" Extension Arm
9A592001	VESA Mount Bar Code Holder
9A591001	Track Mount Bar Code Holder
9A590001	Tamper Proof PC Holder
9A589001	Adjustable PC Holder
	9A593001 9A588001 9A592001 9A591001 9A590001



Midmark[®] 6283 Workstation

- 22" worksurface height adjustment
- 25" W x 10" D worksurface
- Max monitor capacity is 30 lb
- 9.5" independent monitor height adjustment
- Independent 180° monitor rotation and 0-45° tilt
- VESA mount 75 mm and 100 mm
- Folds to 9.47" flat to wall
- Weight capacity: 20-35 lb
- Product weight: 43 lb
- Shipping weight: 50 lb

Model

Description

6283-001

Heavy Duty Wall Mounted Workstation

Accessories

Field Installed	
9A593001	Dual Monitor Arm
9A588001	20" Extension Arm
9A592001	VESA Mount Bar Code Holder
9A591001	Track Mount Bar Code Holder
9A590001	Tamper Proof PC Holder
9A589001	Adjustable PC Holder
	9A593001 9A588001 9A592001 9A591001 9A590001

Midmark Chairs + Tables

Ritter Chairs + Tables

Seating

Lighting

71

Synthesis[®] Cabinetry Collection by Midmark: Limited Warranty

SCOPE OF WARRANTY

Midmark Corporation warrants to the original purchaser its Midmark Modular cabinetry products (Synthesis cabinetry Collection) manufactured by Midmark to be free from defects in material and workmanship under normal interior use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of the parts or the products the defects of which are reported to Midmark within the applicable warranty period and which upon examination by Midmark prove to be defective.

APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows:

- A. Five years for all Midmark brand casework products and components (except for the shorter periods in B and C), including door and drawer fronts and casters
- B. Three years for task lights
- C. One year for countertop surfaces and accessories

EXCLUSIONS

This warranty does not cover sinks, faucets and plumbing accessories. Any warranties on these items are extended directly by the manufacturer of these items to the original purchaser. Information on these manufacturers' warranties will be enclosed with the applicable products. Also, Midmark will furnish copies of any warranties extended by any such manufacturers upon request.

This warranty does not cover and Midmark shall not be liable for: (1) repairs and replacements required because of misuse, abuse, negligence, alteration, accident, freight damage or tampering; (2) matching of color, grain or texture except to commercially acceptable standards; (3) changes in color caused by natural or artificial light; (4) products which are not installed, used or properly cleaned as required in Midmark's written Installation Manual; (5) specially manufactured products; (6) products considered to be of a consumable nature such as light bulbs and surge suppression products; (7) accessories or parts not manufactured by Midmark; (8) charges by anyone (including Midmark's authorized dealers) for adjustments, repairs, replacement parts, installation or other work preformed upon or in connection with such products which is not expressly authorized in writing in advance by Midmark.

EXCLUSIVE REMEDY

Midmarks' only obligation under this warranty is the repair or replacement of defective parts. Midmark shall not be liable for any direct, special, indirect, incidental, exemplary or consequential damages or delay.

NO AUTHORIZATION

No person or firm is authorized to create for Midmark any other obligation or liability in connection with the products.

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

Warranty information included in product manuals

Midmark Synthesis[®] Mobile Solutions is designed for the Clinical Space to be flexible and durable to stand the test of time.



Frame/Panel Colors:

Solids:

216 Pebble 490 Fawn 649 Whisper 791 Frost **Metallics:** 8

804	Venus Silver (M)
805	Pearl Essence (M)
806	Earthen Bronze (M)
807	Flat Iron (M)

Woodgrains:

	-
459	Acorn
487	Earth
794	Dune
795	Flax
796	Nest
797	Hewn
798	Path
799	Sediment
800	Radiance
801	Timber
802	Storm
803	Henna
874	Brazilwood

Others:

793 Pebble PVC-Free

73



Mobile Treatment Cabinet M2

33" H x 27" W x 18" D (2) 4" Drawers + (2) Doors

- Serenity panel style w/ flare handles
- Shipping Weight: 155 lb

* Shown with Kydex®	
contour top	

Product Information				
Step 1: Configuration				
M2				
Step 2. Depth				
1 18" Deep				
Step 3. Option Package				
A 60 mm Casters				
B 4" Quiet Casters				
C 60 mm Casters + Keyed Locks				
D 4" Quiet Casters + Keyed Locks				
Step 4. Cabinetry Finish*				
XXX Solid Color Cabinet Finish (Excludes PVC-Free)				
XXX Woodgrain, Metallic + PVC-Free Cabinet Finish				
Step 5. Countertop Type				
KCN Kydex Contour Top				
SEB Soft Edge Bumper Top				
SEH Soft Edge Handle Top (29"W)				

Configuration Example: M21A.791.KCN

*See page 73 for Cabinetry Finish

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC



Mobile Treatment Cabinet M4

33" H x 27" W x 18" D (3) 4" Drawers + (1) 8" Drawer

- Serenity panel style w/ flare handles
- Shipping Weight: 152 lb

* Shown with soft-edge	è
bumper top	

Product Information				
Step 1: Configuration				
M4				
Step 2. Depth				
1 18" Deep				
Step 3. Option Package				
A 60 mm Casters				
B 4" Quiet Casters				
C 60 mm Casters + Keyed Locks				
D 4" Quiet Casters + Keyed Locks				
Step 4. Cabinetry Finish*				
XXX Solid Color Cabinet Finish (Excludes PVC-Free)				
XXX Woodgrain, Metallic + PVC-Free Cabinet Finish				
Step 5. Countertop Type				
KCN Kydex® Contour Top				
SEB Soft Edge Bumper Top				

SEH Soft Edge Handle Top (29"W)

Configuration Example: M41A.791.SEB

*See page 73 for Cabinetry Finish

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC

Ritter Chairs + Tables

75



Mobile Treatment Cabinet M5

33" H x 27" W x 18" D (5) 4" Drawers

- Serenity panel style w/ flare handles
- Shipping Weight: 159 lb

* Shown with soft-edge
handle top

Product Information				
Step 1: Configuration				
M5				
Step 2. Depth				
1 18" Deep				
Step 3. Option Package				
A 60 mm Casters				
B 4" Quiet Casters				
C 60 mm Casters + Keyed Locks				
D 4" Quiet Casters + Keyed Locks				
Step 4. Cabinetry Finish*				
XXX Solid Color Cabinet Finish (Excludes PVC-Free)				
XXX Woodgrain, Metallic + PVC-Free Cabinet Finish				
Step 5. Countertop Type				
KCN Kydex® Contour Top				
SEB Soft Edge Bumper Top				

SEH Soft Edge Handle Top (29"W)

Configuration Example: M51A.791.SEH

*See page 73 for Cabinetry Finish

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC

76



Mobile Procedure Cart P1

- 43" H x 18" D and available in 18", 21" and 24" W
- (2) 2" drawers, (3) 4" drawers and (1) 8" drawer
- Standard features: Slide-Out Worksurface, No-Tip Base, Full Extension Drawers, 4" Casters and Base Bumper
- Various accessories and accessory mounting locations
- Cascade panel style w/ flare handles
 - Shipping Weight:
 - o 18" 187 lb
 - o 21" 199 lb
 - o 24" 205 lb

Product	Product Information			
Step 1:	Step 1: Configuration			
P1	(2) 2″ drawers, (3) 4″ drawers, (1) 8″ drawer Slide-out Worksurface, 4″ Casters, Base Bumper			
Step 2.				
1	18" wide			
2	21" wide			
3	24" wide			
Step 3.	Options			
Ν	None			
Е	Keyed Lock			
F	Keyless Lock			
Step 4.	Cabinetry Finish	1		
XXX	Solid Color Ca	abinet Finish	(Excludes PVC-Free)	
XXX	Woodgrain, N	1etallic + PVC	C-Free Cabinet Finish	
Step 5.	Countertop Typ	е		
кст	Kydex [®] Conto	ur Top		
SCT	Soft-Edge Co	ntour Top		
	Accessories			
Plant In	stalled	Field Install	ed LHHE, RHHE, LHUA, LHDA, RHUA, RHDA and RPCK are PLANT INSTALLED ONLY	
LHHE			Side Transit Handle, LH mount	
RHHE			Side Transit Handle, RH mount	
LHUA			Side Accessory Rail, LH mount-upper only	
LHDA			Side Accessory Rail, LH mount-dual	
RHUA			Side Accessory Rail, RH mount-upper only	
RHDA			Side Accessory Rail, RH mount-dual	
RBCK			Top Accessory Rack e or more of LHUA, LHDA, RHUA, RHDA or RBCK to hold them	
SBL	ICCessories Delo	9A688001	6 x 6 Storage Storage Bins/Hinged Lid	
MWC		9A689001 9A689001	Waste Basket w/Bracket	
SWH		9A69001 9A690001	Disinfectant Wipe Bracket	
MSC		9A691001	Sharps Cabinet w/Bracket	
SGH		9A691001 9A692001	Single Glove Box Holder	
301		7A072001		

77

*See page 73 for Cabinetry Finish

9A693001

Configuration Example: P12E.216.KCT (Does not include accessories)

TGH

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC

3-Box Glove Holder



Mobile Procedure Cart P2

- 43" H x 18" D and available in 18", 21" and 24" W
- (4) 4" drawers and (1) 8" drawer
- Standard features: Slide-Out Worksurface, No-Tip Base, Full Extension Drawers, 4" Casters and Base Bumper
- Various accessories and accessory mounting locations
- Cascade panel style w/ flare handles
- Shipping Weight:
 - o 18" 187 lb
 - o 21" 199 lb
 - o 24" 205 lb

* Shown in 18" Width

Produc	ct Information		
Step 1:	: Configuration		
P2	(4) 4" drawers, (1) 8" drawer Slide-out Worksurface, 4" Casters, Base Bumper		
Step 2.	. Width		
1	18" wide		
2	21" wide		
3	24" wide		
Step 3	. Options		
N	None		
Е	Keyed Lock		
F	Keyless Lock		
Step 4.	. Cabinetry Finish		
ххх	Solid Color Cabinet Fi	nish (Excludes PVC-Free)	
ххх	Woodgrain, Metallic +	PVC-Free Cabinet Finish	
Step 5.	. Countertop Type		
КСТ	Kydex [®] Contour Top		
SCT	Soft-Edge Contour Top		
	. Accessories nstalled Field In	stalled LHHE, RHHE, LHUA, LHDA, RHUA, RHDA and RPCK are PLANT INSTALLED ONLY	
LHHE		Side Transit Handle, LH mount	
RHHE		Side Transit Handle, RH mount	
LHUA		Side Accessory Rail, LH mount-upper only	
LHDA		Side Accessory Rail, LH mount-dual	
RHUA		Side Accessory Rail, RH mount-upper only	
RHDA		Side Accessory Rail, RH mount-dual	
RBCK		Top Accessory Rack	
These	accessories below require	one or more of LHUA, LHDA, RHUA, RHDA or RBCK to hold them	
SBL	9A68800	1 6 x 6 Storage Storage Bins/Hinged Lid	
MWC	9A68900		
SWH	9A69000		
MSC	9A69100		
SGH	9A69200		
TGH	9A69300		
-	juration Example: P23F.7	94.SCT (Does not include accessories)	
78		*See page 73 for Cabinetry Finish	

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC



Mobile Procedure Cart P3

- 43" H x 18" D and available in 18", 21" and 24" W
- (2) 4" drawers, (1) 23" door w/1 adjustable shelf
- Standard features: Slide-Out Worksurface, No-Tip Base, Full Extension • Drawers, 4" Casters and Base Bumper
- Various accessories and accessory mounting locations
- Cascade panel style w/ flare handles
 - Shipping Weight:
 - o 18" 187 lb
 - o 21" 199 lb o 24" – 205 lb

_			
Produc	t Information		
Step 1:	Configuration		
Р3	(2) 4" drawers, (1) 23" door w/1 adjustable shelf Slide-out Worksurface, 4" Casters, Base Bumper		
Step 2.	Width		
1	18" wide (RH single door)		
2	21" wide (RH	single door)	
3	24" wide (dou	uble doors)	
Step 3.	Options		
N	None		
Е	Keyed Lock		
F	Keyless Lock		
Step 4.	Cabinetry Finis	h	
xxx	Solid Color C	abinet Finisł	n (Excludes PVC-Free)
xxx			/C-Free Cabinet Finish
Step 5.	Countertop Typ		
кст	Kydex [®] Conto	our Top	
SCT	Soft-Edge Contour Top		
	Accessories		
	nstalled	Field Insta	Iled LHHE, RHHE, LHUA, LHDA, RHUA, RHDA and RPCK are PLANT INSTALLED ONLY
LHHE			Side Transit Handle, LH mount
RHHE			Side Transit Handle, RH mount
LHUA LHDA			Side Accessory Rail, LH mount-upper only Side Accessory Rail, LH mount-dual
RHUA			Side Accessory Rail, EH mount-dual Side Accessory Rail, RH mount-upper only
RHDA			Side Accessory Rail, RH mount-dupler only Side Accessory Rail, RH mount-dual
RBCK			Top Accessory Rack
These a	accessories belo	w require or	ne or more of LHUA, LHDA, RHUA, RHDA or RBCK to hold them
SBL		9A688001	6 x 6 Storage Storage Bins/Hinged Lid
MWC		9A689001	Waste Basket w/Bracket
SWH		9A690001	Disinfectant Wipe Bracket
MSC		9A691001	Sharps Cabinet w/Bracket
SGH		9A692001	Single Glove Box Holder
TGH		9A693001	3-Box Glove Holder

Configuration Example: P31E.791.KCT (Does not include accessories)

*See page 73 for Cabinetry Finish

Kydex is a registered trademark of SEKISUI POLYMER INNOVATIONS, LLC

79

Midmark Chairs + Tables

Workstations

Synthesis[®] Cabinetry



Mobile Waste Cart W2

- 18" W x 18" D and available in 31" or 33" H
- Designed to fit underneath Synthesis® Base Cabinets
- Front lock and easy waste removal
- No-touch front design
- Shipping Weight: o 31" – 107 lb
 - o 33" 111 lb

Product Information				
Step 1:	Configuration			
14/2				
W2				
Step 2.	Width			
1	18" wide			
Step 3.	Height			
1	31" high			
2	33" high			
Step 4.	Cabinetry Finish			
XXX	Solid Color Cabinet Finish (Excludes PVC-Free)			
XXX	Woodgrain, Metallic + PVC-Free Cabinet Finish			
Config	uration Example: W211.802			

*See page 73 for Cabinetry Finish



Mobile Linen Cart L1

- 18" W x 18" D and available in 31" or 33" H
- Designed to fit underneath Synthesis[®] Base Cabinets
- Front lock and easy waste removal
- No-touch front design
- Shipping Weight:
- o 31" 115 lb
- o 33" 119 lb

Product Information		
Step 1: Configuration		
L1		
Step 2. Width		
1	18" wide	
Step 3. Height		
1	31" high	
2	33" high	
Step 4. Cabinetry Finish		
XXX	Solid Color Cabinet Finish (Excludes PVC-Free)	
XXX	Woodgrain, Metallic + PVC-Free Cabinet Finish	
Configuration Example: L112.794		

*See page 73 for Cabinetry Finish

Ritter Chairs + Tables

81



Ritter[®] Waste Receptacle

- Includes bag retaining ring
- 14" W x 16" D x 25" H
- Shipping weight: 17 lb
- Capacity: 43 liquid qts
- Color: Pebble

Model

Description

Waste receptacle

261-001



Model

M11-042

M11-042-P

M11-042-DL

Ritter M11[®] Steam Sterilizer

- Automatic door
- Chamber size: 11" x 18"
- Includes pouch rack, 6 slot
- ASME certified
- Sterilizer weight: 99 lb (44.9 kg)
- Shipping weight: 131 lb (59.4 kg)

Description

Ritter® M11 - large tabletop automatic
sterilizer with automatic door
Ritter M11 Steam Sterilizer with Printer (115V)
Ritter M11 Steam Sterilizer with Data Logger (115V)

Pouch Rack Included



83

Ritter

Ritter M9[®] Steam Sterilizer

- Automatic door
- Chamber size: 9" x 15"
- Includes pouch rack, 5 slot
- ASME certified
- Sterilizer weight: 73 lb (33.1 kg)
- Shipping weight: 81 lb (36.7 kg)

Pouch Rack Included



Description

M9-042 M9-042-P M9-042-DL

Model

Ritter® M9 - fully automatic sterilizer with automatic door Ritter M9 Steam Sterilizer with Printer (115V) Ritter M9 Steam Sterilizer with Data Logger (115V)



Ritter[®] M9D

Manual Door

- Chamber size: 9" x 15"
- Includes pouch rack, 5 slot
- ASME certified
- Sterilizer weight: 73 lb (33.1 kg)
- Shipping weight: 81 lb (36.7 kg)



Pouch Rack Included

Model

Description

M9D-042	Ritter M9D - Steam Sterilizer with manual door
M9D-042-P	Ritter M9D Steam Sterilizer with Printer (115V)
M9D-042-DL	Ritter M9D Steam Sterilizer with Data Logger (230V)

Sterilizer Accessories

Model	Description
002-0396-00	Speed-Clean (16 oz single bottle) (A)
002-0396-05	Speed-Clean (12 - 16 oz bottles)
9A307001	M9/M11 Cool Hand Tool (Tray and Cassette Removal Tool) (B)
9A599001	M9-042 and M11-042 Printer (C)
9A682001	Sterilizer Data Logger (D)
002-10502-00	Data Logger USB Device (E)





Model

9A5860021

Double VistaCool[™] Direct-to-Drain System

• Designed to be semestible with Bitter M111® or MO® at

- \bullet Designed to be compatible with Ritter M11 $^{\ensuremath{\texttt{m}}}$ or M9 $^{\ensuremath{\texttt{m}}}$ sterilizer
- Tank dimensions w/fittings: 4" W x 13" L x 17" H (10.2 cm x 33.1 cm x 43.2 cm)
- Tank volume: 3 gal (11.40 liters)
- Tank weight w/water and fittings: 28 lb (12.7 kg)
- Approximate shipping weight: 10 lb (4.5 kg)
- Max coolant water pressure: 100 psi (6.9 bar)

Description

VistaCool double unit for connecting to any two current Ritter[®] sterilizers

¹Required adapter kits included. Ritter M9 and M11 only compatible with the double unit.

Accessories Model

Description

9A586004Air Gap Kit (needed only if required by local code)9A586005Single VistaCool adjustable tank stand9A586006Double VistaCool adjustable tank stand

VistaCool™ Direct-to-Drain System for Autoclave Wastewater is a trademark of Crosstex International, Inc., a Cantel Medical Company, Hauppauge, NY

Midmark Chairs + Tables

85



Ultrasonic Cleaner QC6/QC6R

- Tank Capacity: 6.6 gal/25L
- Tank Size: 19.9" L x 11.8" W x 7.9" H (50.5 x 30 x 20 cm)
- Overall Size: 21.6" L x 16.5" W x 14.3" H (55 x 42 x 36.5 cm)
- Tabletop Unit Weight: 29.6 lb (13.44 kg)
- Tabletop Unit Shipping Weight: 32.6 lb (14.78 kg)
- Recessed Unit Weight: 32.9 lb (14.92 kg)
- Recessed Unit Shipping Weight: 35.8 lb (16.26 kg)
- Safety basket INCLUDED

Model	Description

QC6-01	QuickClean Ultrasonic Cleaner – 6.6 gal Tabletop Unit
QC6R-01	QuickClean Ultrasonic Cleaner – 6.6 gal Recessed Unit



QuickClean®

Ultrasonic Cleaner QC3/QC3R

- Tank Capacity: 3.3 gal/12.5L
- Tank Size: 13" L x 11.8" W x 5.9" H (33 x 30 x 15 cm)
- Overall Size: 14.7 " L x 16.5 " W x 10.4 " H (37.5 x 42 x 26.5 cm)
- Tabletop Unit Weight: 18.6 lb (8.42 kg)
- Tabletop Unit Shipping Weight: 20.9 lb (9.50 kg)
- Recessed Unit Weight: 22.3 lb (10.10 kg)
- Recessed Unit Shipping Weight: 24.6 lb (11.18 kg)
- Safety basket INCLUDED

Model	Description
QC3-01	QuickClean Ultrasonic Cleaner – 3.3 gal Tabletop Unit
QC3R-01	QuickClean Ultrasonic Cleaner – 3.3 gal Recessed Unit

87

BACK TO CONTENTS



Medical Guide Book Effective January 1, 2022

QuickClean®

Ultrasonic Cleaner QC1

- Tank Capacity: 1.2 gal/4.5L
- Tank Size: 11.8" L x 5.9" W x 5.9" H (30 x 15 x 15 cm)
- Overall Size: 13.5" L x 10.2" W x 10.4" H (30 x 26 x 26.5 cm)
- Tabletop Unit Weight: 12.7 lb (5.76 kg)
- Tabletop Unit Shipping Weight: 14.6 lb (6.64 kg)
- Safety basket INCLUDED

Model	Description
QC1-01	QuickClean Ultrasonic Cleaner – 1.2 gal Tabletop Unit

Ultrasonic Parts + Accessories

9A612001	QC1 QuickClean Beaker Accessory Kit
9A613001	QC3/QC3R QuickClean Beaker Accessory Kit
9A614001	QC6/QC6R QuickClean Beaker Accessory Kit
9A298001	Midmark Enzymatic Cleaner (32 oz)
9A296001	Midmark General Purpose Cleaner (32 oz)
9A297001	Midmark Tartar and Stain Remover (32 oz)
002-10007-00	QC1 Safety Basket (1 basket comes standard with QC1 unit)
002-10008-00	QC3/QC3R Safety Basket (1 basket comes standard with QC3/QC3R unit)
002-10009-00	QC6/QC6R Safety Basket (1 basket comes standard with QC6/QC6R unit)



Designing better care®



CARB 93120.2 Phase 2 Compliant and TSCA Title VI Compliant

The color examples shown are the best representation of the original material. Actual color may vary slightly. We strongly recommend that you contact Midmark Customer Experience at 1.800.MIDMARK to request a sample before placing your order.

Midmark is an ISO 13485 and ISO 9001 Certified Company. Certain products are not included. See the complete list at: <u>midmark.com/ISO</u>

HEINE is an ISO 13485 Certified Company.

For more information, contact your Midmark dealer or call: 1.800.MIDMARK Outside the USA call: 1.937.526.3662 or visit our website: <u>midmark.com</u>

© 2022 Midmark Corporation, Miamisburg, Ohio USA Products subject to improvement changes without notice. Litho in USA.

MKT00176 (1/22)

