



# EnzyMed™

Instrument Processing



# Trusted solutions

For over 15 years, the EnzyMed™ Family of cleaning solutions has provided trusted effective chemistries to the CSD's and MDRD's across Canada. Our market leading natural enzyme-based pre-soak & cleaner formula assists in the break down & removal of organic soils at a neutral pH. Enzymed's multi-enzymatic formula offers a broader range of action and excellent in use. It is more effective at low concentrations and lower temperatures than detergents.

**Non-toxic,  
phosphate free,  
& biodegradable**

Safe for staff & an  
environmentally  
conscious choice

**Low foam &  
easy rinsing**

Facilitates faster  
processing for  
your teams

**Concentrated  
formulations  
& fast acting**

Saves you time  
and money

**Fragrance free  
and safe for  
manual cleaning**

Contact your Valued  
representative for more information,  
or visit [www.valuemed.ca](http://www.valuemed.ca)



We deliver **Better Value**

[valuemed.ca](http://valuemed.ca) | [sales@valuemed.ca](mailto:sales@valuemed.ca)







## Asset Protection

Water hardness or poor water quality can affect your chemistry selection. The **Enzymed™ Family** of instrument cleaning solutions offer you the options to achieve the balance needed to **protect your valuable instruments and washers.**



Non-Toxic



Biodegradable



Phosphate Free

## Eco-Friendly

The **Enzymed™ Family** cleaners cause less impact on the environment by reducing waste and overuse with concentrated solutions. Streamline your inventory and use with dosing systems for **clean, safe, and cost effective results.**

- Biodegradable
- Safe and Non-Toxic
- Packaged in Recyclable containers
- **Made in Canada.**

Health Canada registered ISO 9001: 2015 certified facility.



We deliver **Better Value**

[valuemed.ca](http://valuemed.ca) | [sales@valuemed.ca](mailto:sales@valuemed.ca)



## The EnzyMed™ Family

All of our EnzyMed™ solutions include a high volume of multi enzymes, and can be used as an instrument pre-soak and cleaner solution. Our formulation is tinted to make it easily visible.



- Easy Rinsing
- Neutral pH
- Non-Toxic
- Fast Acting
- Biodegradable
- Non-corrosive
- Fragrance free
- Highly visible tint



### EnzyMed™ Transport Gel

An instrument gel pre-cleaner and hydrator.

- Gel clings to instruments
- Colored for easy visibility
- Superior hydration

Item No.	Description	Pkg
MDTG001-12E	Bottle, 1 L	12/case



### EnzyMed™

Our original multi-enzymatic instrument pre-soak & cleaner formula.

- Concentrated
- Multi-Enzyme Formula
- Low Foaming

Item No.	Description	Pkg
MD0001-12E	Bottle, 1 L	12/case
MD0004-16E	Bottle, 4 L	4./case
MD0010-10E	Carboy, 10 L	1 each
MD0020-20E	Pail, 20 L	1 each
MD0056-56E *	Drum, 56 L	1 each
MD0115-115E *	Drum, 115 L	1 each

\* 56 L and 115 L are special order items

Rev. 09-23



### EnzyMed™ Plus

A super concentrated multi-enzymatic instrument pre-soak & cleaner formula.

- Super concentrated
- Multi-Enzyme Formula
- Low Foaming

Item No.	Description	Pkg
MDP0004-16E	Bottle, 4 L	4/case
MDP00010-10E	Carboy, 10 L	1 each
MDP00020-20E	Pail, 20 L	1 each
MDP00056-56E *	Drum, 56 L	1 each
MDP00115-115E *	Drum, 115 L	1 each



### EnzyMed™ Plus AW

An ultra concentrated multi-enzymatic instrument pre-soak & cleaner formula that can be used in manual decontamination and washer processes.

- Ultra concentrated
- Multi-Enzyme Formula
- Low Foaming

Item No.	Description	Pkg
MDPAW00010-10E	Carboy, 10 L	1 each
MDPAW00020-20E	Pail, 20 L	1 each
MDPAW00056-56E *	Drum, 56 L	1 each
MDPAW00115-115E *	Drum, 115 L	1 each



### EnzyMed™ Ortho P

A multi-enzymatic instrument pre-soak & cleaner formula with a grease removal additive.

- Concentrated
- Multi-Enzyme Formula
- Degreasing Formulation
- Low Foaming

Item No.	Description	Pkg
MO0001-12E	Bottle, 1 L	12/case
MO0004-16E	Bottle, 4 L	4/case
MO0010-10E	Carboy, 10 L	1 each
MO0020-20E	Pail, 20 L	1 each
MO0056-56E *	Drum, 56 L	1 each
MO0115-115E *	Drum, 115 L	1 each