



Valuemed

Professional Products Ltd.

**VALUEMED
HANDPIECE
REPAIR 2026**



We now offer handpiece repair!

Our technicians specialize in 30 of the most popular handpieces on the market. Simply call in, email, or message us for pickup and we will schedule the courier.

All pre-approved repairs for the models listed on the back will be repaired and on their way back to you within 72 hours of receiving your approval.

All of our highspeed repairs are backed by our 12 month hassle free guarantee.

Call 1-800-919-8444 or email us at sales@valuemed.ca to get started

FIRST TIME ORDER PROMOTION:

GET A \$50 REPAIR CREDIT TOWARDS YOUR SECOND REPAIR!

2026 HANDPIECE REPAIR PRICE LIST

We currently offer 3 levels of repair on air driven handpieces:

Level 1

05-VMShp-A1

Opening, cleaning, and o-ring gasket replacement **\$129.00**

Level 2

05-VMShp-A2

Includes **Level 1** + bearing and end cap replacement (as required) **\$229.00**

Level 3

05-VMShp-A3

Includes **Level 1 & level 2** + turbine replacement **\$429.99**

(prices are estimated on best value OEM parts and may be subject to an additional cost)

Popular Handpieces to Order from Valuemed

Sable Air High Speeds

SAB-2000017 SABLE ACCESS PRO-NL
 SAB-2000018 SABLE ROTAMAX PRO-NL
 SAB-2000020 SABLE ROTAMAX PRO-KL
 SAB-2000021 SABLE ACCESS Pro KL
 SAB-2000022 Sable Rotamax Pro KL
 SAB-2000023 Sable Rotamax Pro 45-KL
 SAB-2000024 Sable Rotamax Pro 45
 SAB-2000028 SABLE ACCESS K4
 SAB-2000029 SABLE ROTAMAX K4
 SAB-2000030 SABLE ACCESS K6
 SAB-2000031 SABLE ROTAMAX K6

NSK Air High Speeds

NSK M600
 NSK X600/X700
 NSK Z800/Z900

Kavo Air High Speeds

625/630/635/636/640/642/645/646/650
 6000/6500/8000
 E679/680
 4500/5000
 M8900/M9000

Flight Air High Speeds

F-DE-11
 F-140
 F-160
 F-160K
 F-161K
 F-162K
 F-163K

Flight Low Speeds

FDE41PC
 FDE41PS
 F-DE4
 F-DE-01
 F-DE-74
 F-DE11

Sable Air Low Speeds

1600214 One Piece
 1600216 One Piece
 1600218 One Piece
 1600305 Access Prophy
 1600111 20K Mini Motor



04-HPREPAIR

There is a \$62.50 diagnostic fee per handpiece. This charge is waived if the repair is approved, a replacement handpiece is purchase from us, or the handpiece is repaired under warranty. If we do not receive approval of the repair quote within 30 days, the diagnostic fee will be charged, and the handpiece returned to you.