

Hygenic® Non-Latex FR Dental Dam

4+2



- FR designates Flexible Retraction which provides higher elasticity and less retraction force on the isolated area

6"x 6"	Purple	30/pack	15482
6"x 6"	Green	30/pack	15483
6"x 6"	Blue	30/pack	15561

Hygenic® Non-Latex FR Framed Dental Dam

4+2



Purple 20/pack 15514

Hygenic® Non-Latex HR Dental Dam

4+2



- HR designates High Retraction which provides a latex-like fit around the isolated area

5"x 5"	Medium	Teal	15/pack	15552
6"x 6"	Medium	Teal	15/pack	15550
6"x 6"	Medium	Teal	75/pack	15551

Hygenic® Trio Dental Dam

4+2



- Assorted colours of pink, purple, and light blue
- Subtle tropical fruit flavor

6"x 6" Medium 36/pack 15547

Hygenic® Latex Dental Dam

4+2



- Low protein and powder free for greater patient comfort
- 360-piece bulk pack also available for high volume practices

6"x 6"	Thin	Green	36/pack	15498
6"x 6"	Medium	Green	36/pack	15499
6"x 6"	Heavy	Green	36/pack	15500
5"x 5"	Medium	Green	52/pack	15502
6"x 6"	Medium	Blue	36/pack	15530

ISODAM® POLYISOPRENE DAM

- Thirteen-time Dental Advisor top non-latex dental dam award winner
- Light blue colour helps reduce eye fatigue
- Compatible with the use of chloroform



3+1
or
7+3
Standard
Package

4+1
Economy
Package

Standard Package

5"x 5"	Medium	20/pack	ISO01800520
5"x 5"	Heavy	20/pack	ISO02300520
6"x 6"	Medium	15/pack	ISO01800615
6"x 6"	Heavy	15/pack	ISO02300615
6"x 6"	X-Heavy	15/pack	ISO02800615


Economy Package

5"x 5"	Medium	100/pack	ISO01800510
5"x 5"	Heavy	100/pack	ISO02300510
6"x 6"	Medium	75/pack	ISO01800675
6"x 6"	Heavy	75/pack	ISO02300675
6"x 6"	X-Heavy	75/pack	ISO02800675

ISODAM® HD POLYISOPRENE DAM

- Royal blue colour provides enhanced contrast to the teeth
- Compatible with the use of chloroform

5"x 5"	Medium	20/pack	ISOHDB01800520
5"x 5"	Heavy	20/pack	ISOHDB02300520
6"x 6"	Medium	15/pack	ISOHDB01800615
6"x 6"	Heavy	15/pack	ISOHDB02300615



3+1
or
7+3

BIODAM® POLYISOPRENE DAM

- Performs like natural rubber without the potential allergen of latex proteins
- Violet colour provides enhanced contrast to the teeth

6"x 6"	Medium	20/pack	LFBI0VIO01800620
6"x 6"	Heavy	20/pack	LFBI0VIO02300620



2+1

BIODAM® LATEX DAM

- Low protein latex dam
- Green and blue colours

Green	5"x 5"	Medium	52/pack	BIOGRN02412852
	6"x 6"	Medium	36/pack	BIOGRN02415236
Blue	5"x 5"	Medium	52/pack	BIOBLU02412852
	6"x 6"	Medium	36/pack	BIOBLU02415236



5+5

BIODAM®+ LATEX DAM

- Proprietary surface finish aids in the retention of the gum
- Cobalt blue colour

6"x 6"	Medium	36/pack	LABDPBLU01800636
6"x 6"	Heavy	36/pack	LABDPBLU02300636
6"x 6"	X-Heavy	36/pack	LABDPBLU02800636



5+5

Other sizes and gauges of Four D Rubber dental dams are also available