


3+1 FREE: Tofflemire Dead-Soft & Regular Bands

Slick Bands™ - When a sectional matrix just won't do

- Non-stick technology
- Improved tear resistance
- **Exclusive** Garrison shapes

Margin Elevation Band!

The deep margin elevation technique is a useful non-surgical approach to treating deep localized margins.




Right Curve™

Reaches over to the adjacent tooth more easily to help you create broader, deeper contacts.

- Extremely wide preps are simplified.
- Improved visibility for core buildups.

Tighter gingival seal is beneficial for all restorations.



Dead-soft Bands
with non-stick coating

#1 Universal for general restorations

SMT200D10

#13 Pedodontic for short and primary teeth

SMT050D10

#2 Subgingival for deep restorations

SMT300D10

Regular Bands
with non-stick coating

SMT200H10

SMT400H10

SMT050H10

SMT300H10

SMT500H10

Dead Soft Bands	Item	Bulk Qty.	Price
#1 Universal Tofflemire	SMT200D10	100	\$108
Right-Curve	SMT400D10	100	\$108
#13 Pedodontic	SMT050D10	100	\$108
#2 Subgingival	SMT300D10	100	\$108


Regular Bands	Item	Bulk Qty.	Price
#1 Universal Tofflemire	SMT200H10	100	\$108
Right-Curve	SMT400H10	100	\$108
#13 Pedodontic	SMT050H10	100	\$108
#2 Subgingival	SMT300H10	100	\$108

Regular Bands	Item	Bulk Qty.	Price
Margin elevation band	SMT500H10	100	\$141



VariStrip™ Our most popular anterior matrix


- **Curved Form:** Optimum for interproximal anatomy
- **Variable Height:** Tapered design to match tooth height
- **Blue Tint:** Contrast between matrix and tooth



Item	Description	Qty.	Price
VS01	VariStrip	100	\$190

Buy 3 Packs of
VariStrip (VS01)
Get 1 Pack
> FREE! <

FitStrip™ is Your Easy Solution



- Rotate the coloured barrel to lock the strip curved or straight.
- Easily attach a handle for extra control - less hand fatigue for you - more comfortable for your patient!

Refills	Single Sided	Double Sided	4 Pack
Super fine 15 Micron	FPSXFSS-4 0.08 mm thick	FPSXFDS-4 0.11 mm thick	Either Refill \$132
Fine 30 Micron	FPSFSS-4 0.10 mm thick	FPSFDS-4 0.15 mm thick	Either Refill \$132
Medium 46 Micron	FPSMSS-4 0.13 mm thick	FPSMDS-4 0.21 mm thick	Either Refill \$132
Coarse 90 Micron	FPSCSS-4 0.18 mm thick	FPSCDS-4 0.30 mm thick	Either Refill \$132
Serrated	FPSSAW-4 0.05 mm thick		Refill \$132
Medium 46 Micron	N/A	FPSTMDs-4 0.40 mm thick	Refill \$132
Medium 46 Micron	N/A	FPSXTMDs-4 0.50 mm thick	Refill \$132
Super fine 15 Micron	FPSGXFSS-4 0.08 mm thick	N/A	Refill \$132
Fine 30 Micron	FPSGFSS-4 0.10 mm thick	N/A	Refill \$132
Medium 46 Micron	FPSGMSS-4 0.13 mm thick	N/A	Refill \$132
FitStrip® replacement handles (2-pack) FPSHANDLE-2			\$47
IPR Gauge FPSGAUGE			\$65



IPR Kit

FPSK05 contains:

2 handles, 1 IPR gauge, 10 double sided (2 of each .11 mm, .21 mm, .30 mm, .40 mm, and .55 mm)

OFFER:

Buy FPSK05, Get 2
Saws (FPSSAW-1) **FREE!**

Buy 4 Refill Packs, Get 1
Refill Pack **FREE!**