



Fusion™ Anterior

The **Fusion™ Anterior Matrix System** from Garrison® gives you the materials and techniques to restore many types of anterior lesions with confidence.

For years, our customers have shown us how they use posterior metal matrices to restore anterior teeth. When placed vertically, what was the buccal-lingual curvature of the band creates ideal gingival-incisal curvature for an anterior tooth. We worked with these customers to design a metal matrix band that takes this technique to a whole new level of refinement.

Add to that the intelligent design of **Fusion™ Anterior Wedges** and you get a complete system. These radically curved wedges hug the bands to the tooth all the way from the facial aspect around to the lingual, helping to free up your hands. Composite placement and sculpting are simplified, allowing you to easily express your artistic nature restoring beautiful smiles for your patients.



ANK01



ANK02

ANK01 Kit Contains:

50 each Short/Tall Bands
25 each:
XS, SM, M, and LG Wedges

ANK02 Kit Contains:
25 each Short/Tall Bands
10 each:
XS, SM, M, and LG Wedges

Garrison.

order toll free at **888-437-0032** | www.garrisondental.com

PMCA1026OH


Garrison.

150 DeWitt Lane
Spring Lake, MI 49456

Restorative Confidence Starts with the Right Matrix


Keys to Success with ReelMatrix for Deep Margin Elevation

Bands that naturally slide below the height of contour



The tooth's natural height of contour prevents standard Tofflemire bands from sealing properly on deep restorations.

ReelMatrix DME bands proprietary shape allows them to produce an excellent seal and easily slide down the tooth as far as you need them to while their extra-narrow design provides excellent visibility in the deepest cases.




Retainerless design

Band stays in position

- Lightweight extra-small reels won't torque the matrix band out of position or accidentally dislodge it

Excellent access and visibility

- Visibility and access are at a premium in deep margin elevation cases. With ReelMatrix, there's nothing in your way!



DME-01

Deep Margin Elevation Kit Contains:

50 RM500H Margin Elevation band reels
1 RMH-100 Reel Matrix™ tensioning instrument
5 FXH150 3D Fusion™ Firm bands
10 FXH200 3D Fusion™ Firm bands
20 FXH300 3D Fusion™ Firm bands
20 each FXYL, FXBL, FXOR, and FXGR 3D Fusion™ wedges


**Buy DME Kit Get
1 Pack of REELS
FREE!**



Choose from: RM500H, RM100D, RM200D, RM100H, and RM200H

Cervical Matrices

Fill, press, cure. Class Vs are hassle-free when you use **Transparent Cervical Matrices**. Anatomical and flexible, our cervical matrices eliminate the air-inhibited layer and reduce your finishing time while protecting the restoration from contamination. Faster, more profitable Class Vs with outstanding aesthetics will make both you and your patients smile.



Flexible

- adapts to individual anatomy for improved contours
- use it to directly compact the composite into the restoration
- thin edges allow excess composite to flow out for improved anatomy

Convenient

- protects the restoration from contamination
- eliminates time spent hand sculpting
- eliminates air-inhibited layer
- reduces finishing time
- 6 sizes for a perfect fit

Blue Tint

- provides contrast between matrix and tooth structure
- does not compromise composite polymerization

Kit Contains:

60 each:
KCMA10, KCMA20, KCMA30, KCMA40,
KCMA60, KCMA70

1 KCMINST1 Placement Instrument



KCMK2



To redeem free goods, email a scanned copy of the invoice dated 1/10/26 – 20/11/26 to ekammers@garrisondental.com with **PMCA1026OH** in the subject line. If no free item is chosen at time of redemption, Garrison will choose the free item for you. Allow 4-6 weeks for delivery.

END OF SEASON SAVINGS ON RINGS | BANDS | WEDGES


Buy 2 packs of Garrison® rings,
Get 1 ring FREE!

NEW!



All-new asymmetrical **Quad Ring** design with specialized 'driver-tip' maximizes Quad Wedge adaptation and cervical seal.

- Adapts the split-tip Quad wedge to both sides of the embrasure simultaneously
- Seals both cervical margins
- Drives deeper for improved matrix band adaptation and tooth separation
- Soft-Face™ design for superior flash reduction!



Tall rings (orange) with both left and right-side driver tips are available for simplified stacking.

Quad Rings



Q400 – 2/pk
1 each Standard – Left and Right Driver Tips
\$423



Q500 – 2/pk
1 each Tall – Left and Right Driver Tips
\$423

Canadian MSRP

Maximum adaptation

The right and left driver tips of the Quad system allow you to select the ideal ring for each quadrant and for each restorative challenge.



garrisdental.com/quad-resource-center

Strata-G™

\$399/pk



Universal Ring
SG400
Blue 2pk



Tall Ring
SG500
Orange 2pk



Wide Prep Ring
SG600
Green 2pk

Composi-Tight® 3D Fusion™

\$370/pk



Short Ring
FX400
Blue 2pk



Tall Ring
FX500
Orange 2pk



Wide Prep Ring
FX600
Green 2pk

You can buy with complete confidence when you select Garrison products. We support all our products with expert customer service, outstanding warranties, and unrivaled satisfaction guarantees.

It's that simple!

Garrison.

To redeem free goods, email a scanned copy of the invoice dated 1/10/26 – 20/11/26 to ekammers@garrisdental.com with **PMCA1026OH** in the subject line. If no free item is chosen at time of redemption, Garrison will choose the free item for you. Allow 4-6 weeks for delivery.



Buy 4 packs of Garrison® bands,
Get a 5th pack FREE!

Free pack must be of equal value.

NTG Series

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps	Item	Qty.	Price
NTiTi™ Full contour matrix bands for ideal anatomy.									
thickness 0.0015"							NTG100	50	\$93
							NTG150	50	\$93
							NTG175	50	\$93
							NTG200	50	\$93
							NTG300	50	\$93
height >>	NA	NTG100 4.0 mm w/ext 4.4 mm	NTG150 4.1 mm w/ext 6.0 mm	NTG175 5.0 mm w/ext 5.6 mm	NTG200 5.4 mm w/ext 6.6 mm	NTG300 6.0 mm w/ext 8.7 mm			

SG Series

Strata-G™ matrix bands with optimized marginal ridge							Item	Qty.	Price
thickness 0.0016"							SG050	100	\$199
							SG100	100	\$199
							SG175	100	\$199
							SG200	100	\$199
							SG300	100	\$199
height >>		SG050 3.5 mm	SG100 4.5 mm	NA	SG175 5.5 mm	SG200 6.5 mm	SG300 7.5 mm		

FXH Series

Firm matrix bands with eZ-Place™ technology							Item	Qty.	Price
thickness 0.0015"							FXH100	100	\$187
							FXH150	100	\$187
							FXH175	100	\$187
							FXH200	100	\$187
							FXH300	100	\$187
height >>	NA	FXH100 4.3 mm	FXH150 4.3 mm w/ext 5.9 mm	FXH175 5.3 mm	FXH200 6.1 mm	FXH300 6.1 mm w/ext 7.7 mm			

FX Series

Full Curve Bands with non-stick coating and occlusal placement tab							Item	Qty.	Price
thickness 0.0016"							FX100	100	\$199
							FX150	60	\$199
							FX175	100	\$199
							FX200	100	\$199
							FX300	60	\$199
height >>	NA	FX100 4.0 mm w/ext 4.4 mm	FX150 4.1 mm w/ext 6.0 mm	FX175 5.0 mm w/ext 5.6 mm	FX200 5.4 mm w/ext 6.6 mm	FX300 6.0 mm w/ext 8.7 mm			

SM Series

Non-stick coating							Item	Qty.	Price
thickness 0.0016"							SM050	60	\$189
							SM100	100	\$189
							SM150	60	\$189
							SM175	100	\$189
							SM200	100	\$189
							SM300	60	\$189
height >>		SM050 3.2 mm w/ext 4.2 mm	SM100 4.6 mm	SM150 3.8 mm w/ext 5.4 mm	SM175 5.5 mm	SM200 6.4 mm	SM300 6.4 mm w/ext 8.6 mm		

M Series

Original bands							Item	Qty.	Price
thickness 0.0015"							M050	50	\$140
							M100	100	\$140
							M150	50	\$140
							M175	100	\$140
							M200	100	\$140
							M300	50	\$140
height >>		M050 3.2 mm w/ext 4.2 mm	M100 4.6 mm	M150 3.8 mm w/ext 5.4 mm	M175 5.5 mm	M200 6.4 mm	M300 6.4 mm w/ext 8.6 mm		

Canadian MSRP

For more sectional matrix band choices, please go to our website.




Buy 4 packs of Garrison® wedges,
Get a 5th pack FREE!

Free pack must be of equal value.

Quad Wedge

NEW!


Utilizes a revolutionary split-tip design to effectively seal two matrix bands (back-to-back Class IIs) at the same time - even with mismatched cervical margins!



Refills	Item	Qty.	Price
Large Green	QGR	100	\$102
Medium Orange	QOR	100	\$102
Small Blue	QBL	100	\$102
Medium Anatomical Purple	QPL	100	\$102

Strata-G™ Wedge

Hollow underside allows for stacking two wedges while the overall shape and soft fins adapt to concavities and irregularities.



Refills	Item	Qty.	Price
X-Small Yellow	SGYL	100	\$92
Small Blue	SGBL	100	\$92
Medium Orange	SGOR	100	\$92
Large Green	SGGR	100	\$92

Composi-Tight® 3D Fusion™ Wedge


Combining Soft-Face™ adaptive materials, a firm inner core and advanced mechanical features to produce a wedge that truly works. Every time.



Refills	Item	Qty.	Price
X-Small Yellow	FXYL	100	\$89
Small Blue	FXBL	100	\$89
Medium Orange	FXOR	100	\$89
Large Green	FXGR	100	\$89

WedgeWands®


WedgeWands® are designed to twist off the wand, but many prefer to leave it attached to help retract the cheek and for easy removal after the procedure is complete.



Refills	Item	Qty.	Price
X-Small Yellow	SWYL	100	\$80
Small Blue	SWBL	100	\$80
Medium Orange	SWOR	100	\$80
Large Green	SWGR	100	\$80

G-Wedge™

If you prefer traditional wedge placement with either a hemostat or cotton pliers, G-Wedges give you the same great shape as WedgeWands®, just without the applicator wand.



Refills	Item	Qty.	Price
X-Small Yellow	GSWYL	100	\$69
Small Blue	GSWBL	100	\$69
Medium Orange	GSWOR	100	\$69
Large Green	GSWGR	100	\$69

For more more sectional wedge choices, please go to our website.

Canadian MSRP

PMCA1026OH