

# ATS 500A Deluxe Drop Arm Commode

**Features:** Chrome plated steel frame. Drop down armrests

**Seat Width:** 18" | **Seat Depth:** 16" | **Overall Width:** 22 ¼"

**Floor Seat Height:** 19" | 21" (With 2" grow kit) | 22" (With 3" grow kit)

**Floor to Bottom of Frame:** 18"

**Weight Capacity:** 250 lbs


**Sold per box:** 1 per box

**Colours Available:** Chrome with blue upholstery



AADL: BENCHMARK

## Standard Seat Pan



## Optional Toilet Seat c/w Lid

Upgrade at no extra charge



## Optional Padded Seat



## Optional Grow Kit

Available in 2" or 3"



Measurements may vary  $\pm \frac{1}{4}$ " | All products priced individually. Price, Products, Terms & Colours are subject to change without notice.

**ATS 500A Commode**

# ATS 500A Deluxe Drop Arm Commode

## Features



### Pail c/w Lid



### Bucket Options

#### **Pail Specifications**

Pail Diameter: 13"  
(Outside to Outside Lips)  
Pail Diameter: 14"  
(Outside to Outside Tabs)  
Pail Depth: 6 1/4"

### Bed Pan Mount

Stainless Steel Bedpan



Bedpan Mount

Measurements may vary  $\pm \frac{1}{4}$ " | All products priced individually. Price, Products, Terms & Colours are subject to change without notice.

Standard Features 500A