

# ATS 500AF Adjustable Frame Deluxe Drop Arm Commode

20" to 22"

- Features: Chrome plated steel frame. Adjustable seat widths 20" or 22", drop down armrests
- Seat Width: 20"-22"
- Seat Depth: 16"
- Floor to Seat Height: 19" | 21" (With 2" grow kit) | 22" (With 3" grow kit)
- Overall Width: 23 <sup>15</sup>/<sub>16</sub>"
- Floor to Bottom of Frame: 18"
- Weight Capacity: 450 lbs
- Sold per box: 1 per box
- Colours Available: Chrome with black upholstery



AADL: BENCHMARK

**Standard**  
Toilet Seat c/w Lid



**Optional**  
Padded Seat



**Optional**  
Grow Kit



Available in  
2" or 3"

Measurements may vary  $\pm \frac{1}{4}$ " | All products priced individually. Price, Products, Terms & Colours are subject to change without notice.

ATS 500AF Commode

# ATS 500AF Adjustable Frame Deluxe Drop Arm Commode

## Features



### Pail c/w Lid



### Bucket Options

#### **Pail Specifications**

Pail Diameter: 13"  
(Outside to Outside Lips)  
Pail Diameter: 14"  
(Outside to Outside Tabs)  
Pail Depth: 6 1/4"

Measurements may vary  $\pm 1/4$ " | All products priced individually. Price, Products, Terms & Colours are subject to change without notice.

Standard Features 500AF