

ATS 600 Stainless Steel Shower Commode

16" to 22"

Features: Stainless steel frame and a plastic toilet seat. Frame is adjustable by 2".
Optional padded seat designed for individual patients.

Seat Width: 16" | 18" | 20" | 22" | 24"

Seat Depth: 16" | 18" | 20"

Floor to Seat Height: 22"

Floor to Bottom Frame: 19 3/16"

(Without Mags) Overall Width: 16"- 21 1/2" | 18"- 23 1/2" | 20"- 25 1/2" | 22"- 27 1/2" | 24"- 29 1/2"

(Mags) Overall Width: 16"- 25 1/2" | 18"- 27 1/2" | 20"- 29 1/2" | 22"- 31 1/2" | 24"- 33 1/2"

Weight Capacity: 550 lbs

Sold per box: 1 per box

Colours Available: Stainless steel with black upholstery

Adjustable Frame

16" to 18" | 18" to 20" | 20" to 22"



AADL: BENCHMARK

Standard
Toilet Seat Open Front



Optional
Padded Seat



Bucket
Rectangle Bucket



Measurements may vary ± 1/4" | All products priced individually. Price, Products, Terms & Colours are subject to change without notice.