

Masks, Oropharyngeal Airways and Resuscitation Bags

Face Masks



VBM proposes a very large range of face masks with several types of materials and shapes. It covers all patient's size from newborn to adult. There are two main product groups with reusable and single use masks.

Amongst others VBM offers masks with inflatable cushion and silicone lip.

Order information



Silicone Face Masks

reusable

Bi-Mask with silicone lip
Silicone mask with internal plastic shell (PSU)

Size	Patient	Connection	Colour code	REF	Box
2	Child	I.D. 22 mm	green	35-65-222	1
3	Adult, small	I.D. 22 mm	yellow	35-65-223	1
4	Adult, medium	I.D. 22 mm	red	35-65-224	1
5	Adult, large	I.D. 22 mm	purple	35-65-225	1
6	Adult, extra large	I.D. 22 mm	blue	35-65-226	1



Bi-Mask with inflatable cushion
Silicone mask with internal plastic shell (PSU)

Size	Patient	Connection	Colour code	REF	Box
3	Adult, small	I.D. 22 mm	yellow	35-75-223	1
4	Adult, medium	I.D. 22 mm	red	35-75-224	1
5/6	Adult, large	I.D. 22 mm	purple	35-75-225	1



Mask with plastic dome and inflatable cushion

Size	Patient	Connection	Colour code	REF	Box
2	Child	I.D. 22 mm	green	35-70-222	1
3	Adult, small	I.D. 22 mm	yellow	35-70-223	1
4	Adult, medium	I.D. 22 mm	red	35-70-224	1
5	Adult, large	I.D. 22 mm	purple	35-70-225	1
6	Adult, extra large	I.D. 22 mm	blue	35-70-226	1



Rendell Baker Mask

Size	Patient	Connection	REF	Box
0	Newborn	I.D. 22 mm	35-60-000	10
1	Infant	I.D. 22 mm	35-60-001	10
2	Child, small	I.D. 22 mm	35-60-002	10
3	Child	I.D. 22 mm	35-60-003	10



Pediatric Mask, round

Size	Patient	Connection	REF	Box
0	Newborn	I.D. 15 mm	80-11-000	10
1	Infant	I.D. 15 mm	80-11-001	10
2	Child	I.D. 22 mm	80-11-002	10



Order information



PVC Face Masks

for single use

Mask with inflatable cushion

Size	Patient	Connection	Colour code	REF	Box
1	Newborn	I.D. 15 mm	-	REF 35-40-201	10
2	Infant	I.D. 15 mm	-	REF 35-40-202	10
3	Child	I.D. 22 mm	yellow	REF 35-40-203	10
4	Adult, small	I.D. 22 mm	red	REF 35-40-204	10
5	Adult, medium	I.D. 22 mm	purple	REF 35-40-205	10
6	Adult, extra large	I.D. 22 mm	blue	REF 35-40-206	10



Masks, Oropharyngeal Airways and Resuscitation Bags

Face Masks

Accessories



for Silicone and PVC Face Masks

Hook ring, Ø 33 mm
fits masks with plastic dome,
autoclavable, Material: Metal



REF 80-11-266 Box 10



Hook ring, Ø 30-33 mm
fits Bi-Masks,
autoclavable, Material: PSU



REF 35-60-266 Box 10



Fixation tape
autoclavable, Material: Silicone



Patient	REF	Box
Child	35-70-155	10
Adult	35-70-255	10



Oropharyngeal Airways

Guedel Airway

Benefits and Features

- ▶ Atraumatic device made of soft material with efficient bite block
- ▶ ET Tube guidance and fixation thanks to the open flange and lateral clamp
- ▶ Tube-guiding channel with same lumen from proximal to distal end to allow largest possible suction catheters
- ▶ Open side with smooth inner surface facilitates suctioning and easy removal of suction catheter



Order information



Guedel Airway

for single use, sterile, Material: EVA



Size	Colour code	Length	max. ET Tube	Guedel Airway	Guedel Airway with elastic necktape	Box
000	light blue	3.5 cm	I.D. 3.0 mm	REF 40-20-007-1	REF 40-26-007-1	30
00	dark blue	5.0 cm	I.D. 4.0 mm	REF 40-20-008-1	REF 40-26-008-1	30
0	black	6.0 cm	I.D. 4.5 mm	REF 40-20-009-1	REF 40-26-009-1	30
1	white	7.0 cm	I.D. 6.0 mm	REF 40-20-010-1	REF 40-26-010-1	30
2	green	8.5 cm	I.D. 7.0 mm	REF 40-20-020-1	REF 40-26-020-1	30
3	yellow	9.0 cm	I.D. 8.0 mm	REF 40-20-030-1	REF 40-26-030-1	30
4	red	10.0 cm	I.D. 9.0 mm	REF 40-20-040-1	REF 40-26-040-1	30
5	purple	11.0 cm	I.D. 9.5 mm	REF 40-20-050-1	REF 40-26-050-1	30

autoclavable, Material: PUR



Guedel Airway	Box
REF 40-40-007	10
REF 40-40-008	10
REF 40-40-009	10
REF 40-40-010	10
REF 40-40-020	10
REF 40-40-030	10
REF 40-40-040	10
REF 40-40-050	10