

The image features a detailed view of an anesthesia circuit. A large, circular, green corrugated breathing tube is the central focus. It is connected to a white plastic Y-piece. A clear plastic tube with a blue cap is attached to the side of the Y-piece. Below the Y-piece, a blue corrugated tube is visible. The background is a light gray gradient.


LIMITLESS

POSSIBILITIES.

General Anesthesia Circuits and Masks
Selection, expertise and value to meet anesthesia needs

WHAT YOU NEED. WHEN YOU WANT IT.

Whether you choose to customize or go stock, you'll find the selection, expertise and value to meet your needs at Medline. Place your order and count on a strategically located distribution network to deliver innovative solutions on time, every time.




Build The Perfect Circuit.

Use this guide and attached order form to select the best options for each point along the circuit. No extra supplies to sift through, just the components you know you'll use. Whatever your preference, Medline meets the demands of your clinicians and facility.

Anesthesia Turnover Kits.

Reduce SKUs and improve room turnover efficiency with Medline anesthesia turnover kits. Our custom and stock solutions offer exceptional value while increasing ordering efficiency, freeing up storage space and improving product usage.



Reliability


Medline is the largest privately held manufacturer and distributor of healthcare products and #1 manufacturer of custom healthcare solutions in the world.

YOU MAY WANT TO CONSIDER STOCK.

If you're new to Medline general anesthesia supplies, our stock inventory is a great place to start. Work with your Medline representative to fill your anesthesia gear requirements. Once you try them, and find what works best for you, customize your future orders.

Efficiency

With over 30 distribution centers nationwide, and a fleet of 200+ delivery trucks, circuits and supplies are on time, every time.



UNILIMB CIRCUITS

- » Increase the temperature and humidity of inspired gasses
- » Declutter and organize the field
- » Choose from color-coded adult and pediatric circuits
- » Wide variety of accessory options for a complete system

Adult Unilimb Circuits

Item No.	Length	BV filter or HMEF	Breathing bag*	Mask	GS line	Pkg
DYNJAAF4002	40" (102 cm)	1	3 L			20/cs
DYNJAAF4300A	40" (102 cm)	2	3 L		M/M	20/cs
DYNJAAF6002	60" (152 cm)	1	3 L			20/cs
DYNJAAF6380	60" (152 cm)	2	3 L Slim		M/M	20/cs
DYNJAAF6390	60" (152 cm)	2	3 L Slim	DYNJAAMASK5	M/M	20/cs
DYNJAAF6500	60" (152 cm)	2	2 L Curved		M/M	20/cs
DYNJAAF6400	72" (183 cm)	2	3 L		M/M	20/cs
DYNJAAF6410	72" (183 cm)	2	3 L	DYNJAAMASK5	M/M	20/cs
DYNJAAF6463	72" (183 cm)	HMEF	3 L Slim	DYNJAAMASKX5	M/M	20/cs
DYNJAAF9600	96" (244 cm)	2	3 L			20/cs

*All breathing bags are not made with natural rubber latex.

Pediatric Unilimb Circuits

Item No.	Length	BV filter or HMEF	Breathing bag*	Mask	GS line	Pkg
DYNJAPF4000A	40" (102 cm)	1	1 L			20/cs
DYNJAPF6002	60" (152 cm)	1	1 L		M/M	20/cs
DYNJAPF6030	72" (183 cm)	2	1 L		M/M	20/cs

*All breathing bags are not made with natural rubber latex.

EXPANDABLE DUAL LIMB CIRCUITS

Adult Expandable Dual Limb Circuits

Item No.	Length	WYE	BV filter or HMEF	Breathing bag*	Mask	GS line	Pkg
DYNJAA0109	90" (229 cm)	Parallel	1	3 L		M/M	20/cs
DYNJAA0110	90" (229 cm)	Parallel	1	3 L	DYNJAAMASK5	M/M	20/cs
DYNJAA0111†	90" (229 cm)	Parallel	1	3 L	DYNJAAMASK5	M/M	20/cs
DYNJAA0220	90" (229 cm)	Parallel	1	3 L	DYNJAAMASK35	M/M	20/cs
DYNJAA9010A	90" (229 cm)	Parallel	1	3 L	DYNJAAMASK6	M/M	20/cs
DYNJAA9015	90" (229 cm)	Parallel		3 L			20/cs
DYNJAA9020	90" (229 cm)	Parallel		3 L	DYNJAAMASK6	M/M	20/cs
DYNJAA9027	90" (229 cm)	Parallel	2	3 L			20/cs
DYNJAA9120	90" (229 cm)	Parallel	1	3 L			20/cs
DYNJAA9300†	90" (229 cm)	Parallel with port	2	3 L	DYNJAAMASK5	M/M	20/cs
DYNJAA9450	90" (229 cm)	Swivel with port	2	3 L	DYNJAAMASK5	M/M	20/cs
DYNJAA9500	90" (229 cm)	Parallel	2	3 L	DYNJAAMASK6	M/M	20/cs
DYNJAA9980	90" (229 cm)	Parallel	HMEF	3 L		M/M	20/cs
DYNJAA10150	90" (229 cm)	Parallel with port	1	3 L		M/M	20/cs
DYNJAA10160	90" (229 cm)	Parallel with port	1	3 L	DYNJAAMASKX5	M/M	20/cs
DYNJAA10200	90" (229 cm)	Parallel	HMEF	3 L	DYNJAAMASK35	M/M	20/cs
DYNJAA10285	90" (229 cm)	Parallel	HMEF	3 L	DYNJAAMASK5	M/M	20/cs
DYNJAA10230	120" (305 cm)	Parallel	1	3 L		M/M	20/cs
DYNJAA10232	120" (305 cm)	Parallel	2	3 L		M/M	20/cs

*All breathing bags are not made with natural rubber latex. †Yankauer: Bulb Tip, Suction Tubing: 0.25" x 10 ft

Pediatric Expandable Dual Limb Circuits


Item No.	Length	WYE	BV filter or HMEF	Breathing bag*	Mask	GS line	Pkg
DYNJAP7207	72" (183 cm)	Bi-wye		1 L			20/cs
DYNJAP7209	72" (183 cm)	Bi-wye		2 L			20/cs
DYNJAP7211	72" (183 cm)	Bi-wye	1	1 L			20/cs
DYNJAP7213	72" (183 cm)	Bi-wye	1	2 L			20/cs
DYNJAP7237	90" (229 cm)	Bi-wye	1	1 L			20/cs
DYNJAP9810	90" (229 cm)	Parallel	HMEF	1 L		M/M	20/cs
DYNJAP9820	90" (229 cm)	Parallel	1	1 L	DYNJAAMASK3	M/M	20/cs
DYNJAP9830	90" (229 cm)	Parallel	1	1 L		M/M	20/cs
DYNJAP9840	90" (229 cm)	Parallel	2	1 L		M/M	20/cs
DYNJAP9850	90" (229 cm)	Parallel	2	2 L		M/M	20/cs

*All breathing bags are not made with natural rubber latex.

Neonatal Expandable Dual Limb Circuits

Item No.	Length	WYE	BV filter or HMEF	Breathing bag*	Mask	GS line	Pkg
DYNJAP9900	60" (152 cm)	Parallel	2	0.5 L		M/M	20/cs
DYNJAP9901	60" (152 cm)	Parallel	2	0.5 L		M/F	20/cs

*All breathing bags are not made with natural rubber latex.



Work with your Medline representative to choose from the many sizes and configurations to meet your facility's unique needs.

CORRUGATED DUAL LIMB CIRCUITS

Adult Corrugated Dual Limb Circuits

Item No.	Length	WYE	BV filter or HMEF	Breathing bag*	Mask	GS line	Pkg
DYNJAA0011	40" (102 cm)	Parallel	1	3 L			20/cs
DYNJAA0103	40" (102 cm)	Parallel	1	3 L		M/M	20/cs
DYNJAA0104	40" (102 cm)	Parallel	1	3 L	DYNJAAMASK5	M/M	20/cs
DYNJAA4840	40" (102 cm)	Parallel	2	3 L	DYNJAAMASK5	M/M	20/cs
DYNJAA4005	40" (102 cm)	Parallel		3 L			20/cs
DYNJAA6006	60" (152 cm)	Parallel		3 L			20/cs
DYNJAA6011	60" (152 cm)	Parallel	1	3 L			20/cs
DYNJAA6021	60" (152 cm)	Parallel	2	3 L			20/cs
DYNJAA6040	60" (152 cm)	Parallel	1	3 L	DYNJAAMASK7	M/M	20/cs


*All breathing bags are not made with natural rubber latex.

Pediatric Corrugated Dual Limb Circuits

Item No.	Length	WYE	BV filter or HMEF	Breathing bag*	Mask	GS line	Pkg
DYNJAP0010	40" (102 cm)	Parallel	1	1 L		M/M	20/cs
DYNJAP0012	40" (102 cm)	Parallel	2	1 L			20/cs
DYNJAP4003	40" (102 cm)	Bi-wye		1 L			20/cs
DYNJAP4007	40" (102 cm)	Bi-wye		2 L			20/cs
DYNJAP4012	40" (102 cm)	Parallel	1	1 L			20/cs
DYNJAP4015	40" (102 cm)	Bi-wye	1	2 L			20/cs
DYNJAP4025	40" (102 cm)	Parallel	1	2 L		M/M	20/cs
DYNJAP6109	60" (152 cm)	Parallel	1	1 L			20/cs

*All breathing bags are not made with natural rubber latex.

Choose from a variety of accessory options, including port placement, gas sampling lines, masks and HMEs.




ENGINEERED FOR THE RIGHT GRIP.

After observing and interviewing clinicians in operating rooms across the country, we discovered the secrets to making our anesthesia mask the perfect fit in your hand.


Ergonomic Anesthesia Masks

Item No.	Description	Size	Pkg
DYNJAAMASKX4	Small Adult Anesthesia Mask	4	50/cs
DYNJAAMASKX5	Regular Adult Anesthesia Mask	5	50/cs
DYNJAAMASKX6	Large Adult Anesthesia Mask	6	50/cs
DYNJAAMASKX7	Wide Adult Anesthesia Mask	7	50/cs




Flexible Anesthesia Masks

Item No.	Description	Size	Pkg
DYNJAAMASK31	Neonatal Anesthesia Mask	1	50/cs
DYNJAAMASK32	Infant Anesthesia Mask	2	50/cs
DYNJAAMASK33	Child Anesthesia Mask	3	50/cs
DYNJAAMASK34	Small Adult Anesthesia Mask	4	50/cs
DYNJAAMASK35	Regular Adult Anesthesia Mask	5	50/cs
DYNJAAMASK36	Large Adult Anesthesia Mask	6	50/cs



Tail Valve Anesthesia Masks

Item No.	Description	Size	Pkg
DYNJAAMASK1	Neonatal Anesthesia Mask	1	40/cs
DYNJAAMASK2	Infant Anesthesia Mask	2	40/cs
DYNJAAMASK3	Child Anesthesia Mask	3	40/cs
DYNJAAMASK4	Small Adult Anesthesia Mask	4	40/cs
DYNJAAMASK5	Regular Adult Anesthesia Mask	5	20/cs
DYNJAAMASK6	Large Adult Anesthesia Mask	6	20/cs
DYNJAAMASK7	Wide Adult Anesthesia Mask	7	20/cs



Clear dome

Helps visualize patient breathing.



Standard Anesthesia Masks, Pre-Inflated

Item No.	Description	Size	Pkg
DYNJAAMASK9	Small Premie Anesthesia Mask	00	40/cs
DYNJAAMASK10	Premie Anesthesia Mask	0	40/cs
DYNJAAMASK11	Neonatal Anesthesia Mask	1	40/cs
DYNJAAMASK12	Infant Anesthesia Mask	2	40/cs
DYNJAAMASK13	Child Anesthesia Mask	3	40/cs
DYNJAAMASK14	Small Adult Anesthesia Mask	4	40/cs
DYNJAAMASK15	Regular Adult Anesthesia Mask	5	20/cs
DYNJAAMASK16	Large Adult Anesthesia Mask	6	20/cs

Scented Anesthesia Masks

Item No.	Description	Scent	Size	Pkg
DYNJAAMASK31B	Neonatal Anesthesia Mask	Bubble Gum	1	50/cs
DYNJAAMASK31C	Neonatal Anesthesia Mask	Cherry	1	50/cs
DYNJAAMASK31S	Neonatal Anesthesia Mask	Strawberry	1	50/cs
DYNJAAMASK32B	Infant Anesthesia Mask	Bubble Gum	2	50/cs
DYNJAAMASK32C	Infant Anesthesia Mask	Cherry	2	50/cs
DYNJAAMASK32S	Infant Anesthesia Mask	Strawberry	2	50/cs
DYNJAAMASK33B	Child Anesthesia Mask	Bubble Gum	3	50/cs
DYNJAAMASK33C	Child Anesthesia Mask	Cherry	3	50/cs
DYNJAAMASK33S	Child Anesthesia Mask	Strawberry	3	50/cs



HEAT/MOISTURE EXCHANGERS (HME) AND FILTERS

Medline's broad portfolio of HMEs and HMEFs meet your patient's passive humidification needs and provides effective airway humidification.

- » Wide range of tidal volumes
- » Includes non-filtered, electrostatic-filtered, and HEPA filtered options
- » Choose from both angled and straight designs
- » Gas Sampling Ports available with tethered cap

High-Efficiency HMEs

Item No.	Description	Pkg
DYNJAAHME10	HME with Port, Straight	50/cs
DYNJAAHME12	HME No Port, Straight	50/cs
DYNJAPHME18	HME No Port, Straight Neonatal	50/cs
DYNJAPHME20	HME with Port, Straight, Pediatric	50/cs

High-Efficiency Filtered HMEs

- » Products with the BV filter included have an electrostatic hydrophobic membrane that helps protect against cross-contamination
- » Filtered HMEs provide a BFE and VFE of 99.999%
- » HEPA filtration protection available
- » Pleated hydrophobic mechanical filters available

Item No.	Description	Pkg
DYNJAAHME1	HMEF w/Port, Straight, LG	50/cs
DYNJAAHME5	HMEF No Port, Straight, LG	50/cs
DYNJAAHME9	HMEF w/Port, Straight	50/cs
DYNJAAHME11	HMEF No Port, Straight	50/cs
DYNJAAHME13	HMEF w/ Port, Angled	50/cs
DYNJAPHME17	HMEF w/ Port, Straight, Pediatric	50/cs
DYNJAAHMEHEPA	HMEF HEPA, Adult	50/cs

High-Efficiency Filters

- » HEPA filter provides a BFE and VFE of 99.9999%
- » BV filter provides a BFE and VFE of 99.97%

Item No.	Description	Pkg
DYNJAAHEPA	Filter, HEPA	50/cs
DYNJAABV1	Filter, Bacterial Viral	50/cs

*A breathing circuit bacterial viral filter is a device that is intended to remove microbiological and particulate matter from the gasses in the breathing circuit.



DYNJAAHME10



DYNJAPHME18



DYNJAAHME9



DYNJAAHME13

Medline also provides combination HME and bacterial viral protection* in a complete offering of HMEFs.

ANESTHESIA ACCESSORIES

Breathing Bags*

- » Accommodates all patient needs and physician preferences
- » Offered in a variety of sizes, from 0.5L to 3L
- » Available in 2 pleats or a slim 3-pleat style that inflates in a narrow shape to accommodate various hand sizes

Soft Bushing

Item No.	Description	Pkg
DYNJAA005S	Breathing Bag, 0.5L LF	20/cs
DYNJAA01S	Breathing Bag, 1L LF	20/cs
DYNJAA02S	Breathing Bag, 2L LF	20/cs
DYNJAA03S	Breathing Bag, 3L LF	20/cs
DYNJAACRV02S	Breathing Bag, 2L, Curved LF	20/cs

*All breathing bags are not made with natural rubber latex.

Hard Bushing

Item No.	Description	Pkg
DYNJAA005	Breathing Bag, 0.5L LF	20/cs
DYNJAA01	Breathing Bag, 1L LF	20/cs
DYNJAA02	Breathing Bag, 2L LF	20/cs
DYNJAA03	Breathing Bag, 3L LF	20/cs
DYNJAACRV02	Breathing Bag, 2L, Curved LF	20/cs

*All breathing bags are not made with natural rubber latex.



Soft bushing

New soft bushing is more flexible to allow for a more secure grip on the anesthesia machine.



DYNJAA005, DYNJAA01,
DYNJAA02, DYNJAA03

All Medline breathing bags are

ISO 5362 compliant.

Each bag is inspected to ensure the highest levels of quality and delivers uniform inflation to provide a consistent feel.

Gas Sampling Lines

Each connector meets ISO standards and ensures a secure fit to both anesthesia circuit and machine.

- » Choose from coextruded and standard gas sampling lines
- » Both M/M and M/F connections available

Item No.	Description	Pkg
DYNJAA04	Gas Sampling Line (0.06 ID), 10 ft., Male/Male	50/cs
DYNJAA05	Gas Sampling Line (0.06 ID), 10 ft., Male/Male w/ Filter, 0.8 Micron	50/cs
DYNJAA06	Gas Sampling Line (0.06 ID), 10 ft., Male/Female w/ Filter, 0.8 Micron	50/cs
DYNJAA07	Gas Sampling Line Filter, 0.8 Micron	50/cs
DYNJAA11	Gas Sampling Line (0.06 ID), 10 ft., Male/Female	50/cs
DYNJAA15	Gas Sampling Line (0.047 ID), 10 ft., Male/Male, Coextruded	50/cs
DYNJAA18	Gas Sampling Line (0.05 ID), 10 ft., Male/Male	50/cs
DYNJAA22	Gas Sampling Line (0.06 ID), 15 ft., Male/Male	50/cs



DYNJAA04

Anesthesia Accessories

- » Wide variety of connectors with or without gas sampling ports available in bulk or packaged with a circuit
- » Extension sets available to add length to your current anesthesia circuits in both pediatric and adult sizes

Item No.	Description	Pkg
DYNJAA10	Adult Expandable Extension Set, 30"-90" (76-229 cm), 2 tubes	20/cs
DYNJAA20	Pediatric Expandable Extension Set, 30"-90" (76-229 cm), 2 tubes	20/cs
DYNJAA08	Straight Connector w/ Gas Sampling Port	50/cs
DYNJAA14	Flexible Catheter Mount	50/cs
DYNJAATREE1	Steel Anesthesia Tubing Holder, Tree	50/cs
DYNJAASTRAP	Anesthesia Head Strap	30/cs



DYNJAA14



DYNJAA08

With customized circuits you pay for only the products you know you'll use.

Visit medline.com/go/customkits for more information.



Custom Anesthesia Circuit Ordering Guide

Account number/name:


Custom description: (27 characters max)

Monthly usage (each):


Custom Circuit Options. Most popular highlighted in blue.

Tubing


Expandable



Unlimb



Corrugated



BV Filters


BV Filters
HEPA Filter

Connections


Tube to Wye



Tube to Bacterial/Viral Filter




Elbow




Elbow Cap




Additional Options



Mask




Breathing Bag



Gas Sampling Line



Heat Moisture Exchanger





Syringes/Needles

Suction Tubing

Suction Canister

Pulse Oximeter

Monitor Electrodes

Patient Positioners

Airways

Temperature Monitoring

Esophageal Stethoscope

Multi Use Temp Probe

Tympanic Temp Probe

Skin Temp Probe

Turnover Supplies

Other

Manuf./Item #:					
Manuf./Item #:					
Manuf./Item #:					

Airways

Dimensions	Berman	Colored Berman	Hard Guedel	Soft Guedel
30 mm				
40 mm				
50 mm				
60 mm				
70 mm				
80 mm				
90 mm				
100 mm				
110 mm				
120 mm				

Manuf/Item #: _____

Laryngeal Mask Airway

Disposable

Size 1	
Size 1.5	
Size 2	
Size 2.5	
Size 3	
Size 4	
Size 5	

Manuf/Item #: _____

[illegible]

Other Items

NOTES

[illegible]

THE SOURCE FOR YOUR ANESTHESIA PRODUCT NEEDS.

Medline is your full-service anesthesia partner. Our broad portfolio includes supplies for general anesthesia, monitored anesthesia care, pain management and regional anesthesia.

Contact your Medline Representative to help you design a customized product package.



Medline Industries, Inc. Three Lakes Drive, Northfield, IL 60093

Medline United States
1-800-MEDLINE (633-5463)
medline.com | info@medline.com

Medline Canada
1-800-396-6996
medline.ca | canada@medline.com

Medline México
01-800-831-0898
medlinemexico.com | mexico@medline.com

FOLLOW US

