

Why should you upgrade to CaviWipes™ 2.0?



While both CaviWipes™ and CaviWipes 2.0 offer effective disinfection and great material compatibility, here is why you might want to **switch to CaviWipes 2.0!**



CaviWipes™



CaviWipes™ 2.0

Clinical Efficacy	16 Pathogens	55 Pathogens
Steps Required for Disinfection	2-Steps	1-Step ¹
Quaternary Ammonia (Active Ingredient)	1 st Gen	5 th Gen ²
Disinfection Time	3 Minutes	2 Minutes
Formula and pH	Non-Buffered Formula	Buffered Formula to Stabilize pH (Neutral pH Formula)

1. Precleaning is required for visible soil, bloodborne pathogens and *Candida auris*.
2. Fifth generation Quaternary ammonia compounds have outstanding germicidal performance, are active under more hostile conditions and are safer to use.

To Learn More
Scan the QR code
or visit **Metrex.com**

