



Limit uncertainty with quality instrumentation
from the brand you trust

Tip-It™ Instrument Protectors help protect instrument tips and edges from damage during handling and sterilization. Integra Miltex offers a line of Tip-It Instrument Protectors in four styles: Clear Vented, Tinted Vented, Solid Vented, and Solid Non-Vented style.

Tip-It Instrument Protectors feature:

- Four styles in nine colors and sizes for a variety of choices
- Available in vented styles to help reduce vacuum and moisture build up, allowing for easy removal from the instrument
- Steam and gas permeable to allow for sterilization of instruments
- Available in packages of 50
- Non-sterile, single-use only



** Images not to scale

Reference Clear, Vented	Reference Tinted, Vented	Reference Solid, Vented	Reference Solid, Non-Vented	Dimensions	Size Code	Color
3-2501C	3-2501T	3-2501V	3-2501	1/16" X 3/4" (1.6mm X 19mm)	1	White*
3-2502C	3-2502T	3-2502V	3-2502	5/64" X 3/4" (2mm X 19mm)	2	Blue*
3-2503C	3-2503T	3-2503V	3-2503	7/64" X 3/4" (2.8mm X 19mm)	3	Green*
3-2504C	3-2504T	3-2504V	3-2504	1/8" X 1" (3.2mm X 25.4mm)	4	Red*
3-2505C	3-2505T	3-2505V	3-2505	3/16" X 1" (4.8mm X 25.4mm)	5	Yellow*
3-2506C	3-2506T	3-2506V	3-2506	13/32" X 3/4" (10.3mm X 19mm)	6	Maroon*
3-2507C	3-2507T	3-2507V	3-2507	1/16" X 3/8" X 1" (1.6mm X 9.5mm X 25.4mm)	7	Orange*
3-2508C	3-2508T	3-2508V	3-2508	5/64" X 5/8" X 1" (2mm X 16mm X 25.4mm)	8	Brown*
3-2509C	3-2509T	3-2509V	3-2509	1/8" X 1" X 1" (3.2mm X 25.4mm X 25.4mm)	9	Black*
3-2510C	3-2510T	3-2510V	3-2510	Assorted, equal number of sizes plus 5	-	-

*Color does not apply to clear vented pattern

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

Integra ■ 589 Davies Drive, York, PA 17402
 866-854-8300 USA ■ 717-840-2763 outside USA
 866-854-8400 fax
integralife.com/integra-miltex

Manufactured for:

Integra York PA, Inc.
 589 Davies Drive
 York, PA 17402 ■ USA