

ASTM LEVEL 3 MASK

PROCEDURE EAR LOOPS | INDIGO

REF PG4-1273

ABOUT THE MASK

Based on ASTM F2100 standards, this mask provides the highest level of protection during medical procedures. Advanced fluid-resistant and filtration capabilities of the mask provide all-in-one protection for wound care, line changes and dental work. This mask is recommended for procedures with a high risk of fluid splash, when handling trauma or caring for ER and ICU patients.

PROTECTION STANDARDS



BFE (%):	≥ 98
PFE (%):	≥ 98
ΔP (mmH ₂ O/cm ²):	< 6.0
Fluid Resistance:	160mm Hg
Flame Spread:	Class 1

Variation in lot results exist.
All results meet or exceed
ASTM F2100 requirements.

FEATURES & BENEFITS

ULTRASONICALLY WELDED SEAMS

Ensures mask layers remain secure and rest comfortably against the clinician's skin.

WIDE NOSEPIECE

Constructed with high-quality aluminum to keep the mask in place across the nose and cheekbones.

SOFT MASK AND EAR LOOP MATERIAL

Carefully engineered and clinically preferred for extended periods of use without discomfort.



QUANTITIES

MASKS/BOX	BOXES/CASE	MASKS/CASE
50	10	500



PRIMED Medical Products Inc.
200, 2003-91 Street SW, Edmonton, AB, Canada, T6X 0W8
TF +1.877.877.4633 | F +1.780.497.7670 | W PRIMED.com | E info@PRIMED.com

NA-REVB1