

# ASTM LEVEL 3 VISOR MASK

ANTI-FOG | ANTI-GLARE | PROCEDURE EAR LOOPS | INDIGO

REF PG4-1573

## ABOUT THE MASK

Based on ASTM F2100 standards, this mask provides the highest level of protection during medical procedures. The attached visor and advanced fluid-resistant and filtration capabilities of the mask provide all-in-one protection for wound care, line changes and dental work. This mask is recommended for procedures with a high risk of fluid splash, when handling trauma or caring for ER and ICU patients.

## PROTECTION STANDARDS



BFE (%):	≥ 98
PFE (%):	≥ 98
$\Delta P$ (mmH <sub>2</sub> O/cm <sup>2</sup> ):	< 6.0
Fluid Resistance:	160mm Hg
Flame Spread:	Class 1

Variation in lot results exist.  
All results meet or exceed  
ASTM F2100 requirements.

## FEATURES & BENEFITS

### ANTI-GLARE VISOR

Designed with wrap-around wings, anti-glare film and anti-fog properties, providing all-in-one splash protection with an unobstructed view.

### ULTRASONICALLY WELDED SEAMS

Ensures mask layers remain secure and rest comfortably against the clinician's skin.

### WIDE NOSEPIECE WITH ANTI-FOG FOAM

Constructed with high-quality aluminum, anti-fog foam and a black anti-glare strip, providing enhanced visibility and comfort.

### SOFT MASK AND EAR LOOP MATERIAL

Carefully engineered and clinically preferred for extended periods of use without discomfort.



## QUANTITIES

MASKS/BOX	BOXES/CASE	MASKS/CASE
25	4	100



PRIMED Medical Products Inc.  
200, 2003-91 Street SW, Edmonton, AB, Canada, T6X 0W8  
TF +1.877.877.4633 | F +1.780.497.7670 | W PRIMED.com | E info@PRIMED.com

NA-REV B1