

ASTM LEVEL 3 VISOR MASK

ANTI-FOG | ANTI-GLARE | SURGICAL TIES | INDIGO

REF PG4-5073

ABOUT THE MASK

Based on ASTM F2100 standards, this mask provides the highest level of protection during surgical procedures. The attached visor and advanced fluid-resistant and filtration capabilities of the mask provide all-in-one protection for procedures with a high risk of fluid splash or surgical debris, such as orthopedic surgery or trauma.

PROTECTION STANDARDS



BFE (%): ≥ 98
PFE (%): ≥ 98
 ΔP (mmH₂O/cm²): < 6.0
Fluid Resistance: 160mm Hg
Flame Spread: Class 1

Variation in lot results exist.
All results meet or exceed
ASTM F2100 requirements.

FEATURES & BENEFITS

LONG SURGICAL TIES

Long ties allow for easy donning and doffing, with flexibility to adjust the tightness and fit of the mask over a surgical or bouffant cap.

ANTI-GLARE VISOR

Designed with wrap-around wings, anti-glare film and anti-fog properties, providing all-in-one splash protection with an unobstructed view.

WIDE NOSEPIECE WITH ANTI-FOG FOAM

Constructed with high-quality aluminum, anti-fog foam and a black anti-glare strip, providing enhanced visibility and comfort.

ULTRASONICALLY WELDED SEAMS

Ensures mask layers remain secure and rest comfortably against the clinician's skin.



QUANTITIES

MASKS/BOX	BOXES/CASE	MASKS/CASE
25	4	100



PRIMED Medical Products Inc.
200, 2003-91 Street SW, Edmonton, AB, Canada, T6X 0W8
TF +1.877.877.4633 | F +1.780.497.7670 | W PRIMED.com | E info@PRIMED.com

NA-REVB1