


enbio®




Enbio[®]
World's fastest autoclave.
Vacuum+drying.

Designed  in Switzerland

Capacity



x8



Load space
11x7x1.3 in

Max. packages amount (circa) **8**
Max. recommended batch weight **1.1 lb**
Process chamber volume **2.7 l**



Power supply

Enbio^s works with any typical electrical outlet.



Data archiving

All process data are automatically written to a USB disk.



Water supply

Enbio uses distilled water from external tank.



Process reports

Enbio Data Viewer (included) allows generating and printing detailed process reports from the data on USB.

Dimensions



Pressure & temperature

Working pressure **2.1 Bar/30.5 psi**
Max. pressure **2.3 Bar/33.4 psi**
Max. process temperature **137°C/278.6°F**

Electrical

Power supply **110-120V/60Hz**
Installed power **3.23 kW**
Max. power consumption **15 A**

Misc

Ingress protection **IP20**
Data archiving **USB**
Water conductivity* **< 15 µS/cm**

* - demineralized water only



NOW AVAILABLE IN CANADA
For more information, contact
Valuemed Professional Products Ltd.
sales@valuemed.ca | 1-800-919-8444
Visit us at www.valuemed.ca



Enbio Group AG
Eichengasse 3, 4702 Oensingen
Switzerland

FDA Establishment Registration Number: 301677357