




Valuemed
Professional Products Ltd.

ENBIO®



Enbio^S
World's fastest autoclave.
Vacuum+drying.

Designed 
in Switzerland

Bone dry pouches in no time.



15 min





Enbio can sterilize a load of wrapped instruments in just 15 minutes.

Full drying means that pouches are bone dry. Not somewhat dry.

In case of emergency, when you need an instrument sterilized for immediate use, FAST process does the job in just 7 minutes.

7 min



4 full cycles with drying per hour.



Enbio S is 6 times smaller and 3 times lighter than a typical vacuum autoclave.


But when it comes to work efficiency, it keeps up the pace.

	Pouches / h	Instruments / h
Enbio ^S	32	50



The new definition of sterilization speed.

Time is a relative value. It flows differently under different circumstances.
But when you need to have your instruments ready, you need them, well, FAST.
While other autoclaves might test your patience, Enbio bends the laws
of physics to cut the waiting to the invisible.





Inaudibly silent.





Enbio autoclaves work unbelievably silent.
Wherever you use it, nobody will notice it's working.

Enbio^S. Speedmaster.

Enbio S is compact to fit even a confined space. Unbelievably short process times will save your procedure in emergency, or allow sterilization of the instruments as you go, between patients.

Unwrapped
FAST / 134°C
7 min



Unit weight

33 lb

Capacity

Max. packages amount (circa)	8
Max. recommended batch weight	1,1 lb
Process chamber volume	2,7 l

Pressure & temperature

Working pressure	2,1 Bar
Max. pressure	2,3 Bar
Max. process temperature	137°C





Vacuum

Pouches
+drying

FDA Cleared

Wrapped
Vacuum / 134°C

15 min

Wrapped
Vacuum / 121°C

45 min



reddot design award
winner



Dimensions

Width **9.8 in**
Height **6.3 in**
Depth **22 in**

Electrical

Power supply **110V/60Hz**
Installed power **3,23 kW**
Max. power consumption **15 A**

Misc

Ingress protection level **IP20**
Data archiving **USB stick**
Water conductivity* **< 15 µS/cm**

* - demineralized water only

enbio

www.valuemed.ca | sales@valuemed.ca | 1-800-919-8444



Enbio user's manual




1. Plug in
2. Choose program
3. Start


Installation and operation of Enbio autoclaves
does not require any special training.



Enbio Care


Perfect technology meets perfect care.
All Enbio devices are covered with world's unique
warranty and care program, with backup device
delivered to you for the time of repair.






20.000 units in 4 years.

In just 4 years, Enbio's worldwide community has grown to 20,000 users, who sterilize instruments in our autoclaves, for the complete safety of their customers and patients.









Designed  in Switzerland

Worldwide presence.

Enbio autoclaves are manufactured in the EU and distributed globally through logistics centers covering Europe, North America and the Middle East.



enbio®

Designed 
in Switzerland



For more information, contact
Valuemed Professional Products Ltd.
sales@valuemed.ca | 1-800-919-8444
Visit us at www.valuemed.ca



Enbio Group AG
Eichengasse 3, 4702 Oensingen
Switzerland

**FDA Establishment Registration
Number: 3016777357**