

ALPEN®

ShapeGuard

Strong in Polishing.
Gentle on Surfaces.




 COLTENE
www.coltene.com

 COLTENE

ADVANTAGES AT A GLANCE

A perfect complement to existing ALPEN® diamond impregnated polishing systems.



- **Adapts** to any surface
- **Uniform pressure distribution** due to flexible lamellas
→ Morphology of tooth is preserved
- **Fast, easy and intuitive** finishing & polishing
- **Optimized** polishing systems for composite, ceramic, and zirconia
- Outstanding **high luster** polishing results

ORDER INFORMATION

ShapeGuard Composite Plus

Adds a new head shape to our well established 2-step diamond impregnated Composite Polishing System **Composite Plus Step ①** and **Composite Plus Step ②**.



RPM: 10,000 – 12,000

Step ①

Step ②



Dia. in 1/10 mm	110	140	110	140
Head L in mm	1.60	1.60	1.60	1.60
Polisher type	Spiral Wheel			

RA/ 204

Order No.	60021904	60021905	60021906	60021907
Content	5	5	5	5

60021910 ShapeGuard Composite Plus Assorted Pack

Content	1	1	1	1
---------	---	---	---	---

ShapeGuard Ceramic Plus

ShapeGuard is also available as an expansion to the broadly used 3-step diamond impregnated Ceramic Polishing System **Ceramic Plus Step ①**, **Ceramic Plus Step ②** and **Ceramic Plus Step ③**.



RPM: 10,000 – 12,000

Step ①

Step ②

Step ③



Dia. in 1/10 mm	140	140	140
Head L in mm	1.60	1.60	1.60
Polisher type	Spiral Wheel		

RA/ 204

Order No.	60021911	60021912	60021913
Content	5	5	5

60021914 ShapeGuard Ceramic Plus Assorted Pack

Content	1	1	1
---------	---	---	---

ShapeGuard Zirconia Plus

Is now available as an enhancement to our 2-step diamond impregnated Zirconia Polishing System **Zirconia Plus Step ①** and **Zirconia Plus Step ②**.



RPM: 10,000 – 12,000

Step ①

Step ②



Dia. in 1/10 mm	140	140
Head L in mm	1.60	1.60
Polisher type	Spiral Wheel	

RA/ 204

Order No.	60022205	60022206
Content	5	5

60022207 ShapeGuard Zirconia Plus Assorted Pack

Content	2	2
---------	---	---

ORDER INFORMATION

COMPOSITE PLUS

Step ① Shaping & Smoothing



Shape	Minipoint	Point	Cup	Small Lens	Spiral Wheel	Spiral Wheel
Dia. in 1/10 mm	030	040	060	100	110	140
Head L in mm	7.50	10.00	7.50	2.50	1.60	1.60

RA/ 204

5 pcs	Order No.	RP230291RA	RP230293RA	RP230295RA	RP230297RA	60021904	60021905
-------	-----------	------------	------------	------------	------------	----------	----------

Step ② High Luster Polishing



Shape	Minipoint	Point	Cup	Small Lens	Spiral Wheel	Spiral Wheel
Dia. in 1/10 mm	030	040	060	100	110	140
Head L in mm	7.50	10.00	7.50	2.50	1.60	1.60

RA/ 204

5 pcs	Order No.	RP230299RA	RP230301RA	RP230303RA	RP230305RA	60021906	60021907
-------	-----------	------------	------------	------------	------------	----------	----------

CERAMIC PLUS

Step ① Removing & Contouring



Shape	Minipoint	Flame	Cup	Small Lens	Wheel	Spiral Wheel
Dia. in 1/10 mm	030	040	060	100	110	140
Head L in mm	7.50	10.00	7.50	2.50	2.00	1.60

RA/ 204

5 pcs	Order No.	RP230313RA	RP230317RA	RP230319RA	RP230321RA	RP230325RA	60021911
-------	-----------	------------	------------	------------	------------	------------	----------

Step ② Smoothing



Shape	Minipoint	Flame	Cup	Small Lens	Wheel	Spiral Wheel
Dia. in 1/10 mm	030	040	060	100	110	140
Head L in mm	7.50	10.00	7.50	2.50	2.00	1.60

RA/ 204

5 pcs	Order No.	RP230331RA	RP230335RA	RP230337RA	RP230339RA	RP230343RA	60021912
-------	-----------	------------	------------	------------	------------	------------	----------

Step ③ High Luster Polishing



Shape	Minipoint	Flame	Cup	Small Lens	Wheel	Spiral Wheel
Dia. in 1/10 mm	030	040	060	100	110	140
Head L in mm	7.50	10.00	7.50	2.50	2.00	1.60

RA/ 204

5 pcs	Order No.	RP230349RA	RP230353RA	RP230355RA	RP230357RA	RP230361RA	60021913
-------	-----------	------------	------------	------------	------------	------------	----------

ZIRCONIA PLUS

Step ① Smoothing



Shape	Flame	Cup	Spiral Wheel
Dia. in 1/10 mm	040	060	140
Head L in mm	10.00	7.50	1.60

RA/ 204

5 pcs	Order No.	RP230379RA	RP230381RA	60022205
-------	-----------	------------	------------	----------

Step ② High Luster Polish



Shape	Flame	Cup	Spiral Wheel
Dia. in 1/10 mm	040	060	140
Head L in mm	10.00	7.50	1.60


RA/ 204

5 pcs	Order No.	RP230391RA	RP230393RA	60022206
-------	-----------	------------	------------	----------

APPLICATION POSSIBILITIES


ShapeGuard facilitates intraoral finishing and polishing of the following areas:

Interproximal
Occlusal
Incisal
Labial
Buccal
Lingual



UTILIZATION

Exert minimal, constant pressure.
Use wet when polishing intraorally.





Coltène/Whaledent Inc. 
235 Ascot Parkway
Cuyahoga Falls, OH 44223 / USA
Tel. USA & Canada +1 800 221 3046
+1 330 916 8800
Fax +1 330 645 8704
info.us@coltene.com

Coltène/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten / Switzerland
Tel. +41 (0)71 757 53 00
Fax +41 (0)71 757 53 01
info.ch@coltene.com

Coltène/Whaledent GmbH + Co. KG
Raiffeisenstrasse 30
89129 Langenau / Germany
Tel. +49 (0)7345 805 0
Fax +49 (0)7345 805 201
info.de@coltene.com

 **COLTENE**

P/N 31417B