

CanalPro X-Move

CORDLESS ENDO MOTOR

The CanalPro X-Move is a cordless endodontic motor with an integrated root canal apex locator. It can be used as an endodontic motor for preparation of root canals and as an apex locator device to determine the root canal working length. The time-saving and intuitive user interface includes pre-programmed settings for HyFlex files.

- Continuous rotation
- Reciprocating motion with adjustable angles
- Integrated Apex Locator
- Wireless operation
- Super mini insulated contra angle
- Pre-programmed settings for COLTENE files



CanalPro X-Move

CORDLESS ENDO MOTOR

CanalPro X-Move Endo Motor

65002739	X-Move Motor	65002747	Measuring wire
65002742	X-Move Motor Handpiece	65002748	File Holder
65002743	Charging base	65002749	Lip Clip
65002744	Contra angle	65002750	Touch Probe
65002745	Spray Nozzle	65002752	Battery
65002746	Power Adapter		

