Coltolux Comfort

Coltolux" Comfort

SCOLTENE

Simple. Efficient. Reliable.



Coltolux Comfort





Simple. Efficient. Reliable.

Successor to the reliable Coltolux LED, the Coltolux Comfort is a simple to use, one button, single-mode cordless curing light. Its powerful 1200mW/cm2 irradiance is compatible with all campherquinone-

The Coltolux Comfort was designed with the doctor and patient in mind. Lightweight and ergonomic, it can reduce hand arm fatigue, while the flat tip 90 angle provides easy access to the posterior, offering more comfort for the patient. The newly designed intraoral light shield allows placement of the lens to be flush with the curing plane blocking reflected blue light from the tooth surface.

The convenient user replaceable battery provides a minimum of 120 full power 20 second curing cycles with a fully charged battery. With an 11mm curing diameter, the homogenous beam profile provides efficient and reliable curing, while the durable material construction makes the Coltolux Comfort easy to disinfect.

O ····· Non-contact inductive charging is durable and protects against corrosion

Coltolux Comfort LED Curing Ligtht

Technical Data		
Irradiance	1200 mW/cm2 (+/-10%)	
Wavelength range	380 nm to 515 nm	
Peak wavelength	449 nm	
Handpiece Power Supply		
Lithium ion battery	Minimum of 120, twenty-second cures with a new, fully charged battery	
Charger Base Power Supply		
Operation volatage	100-240v	
Frequency	50/60 Hz	
Operation		
Operation	One-button, single-mode	
Curing Cycle	20-second cycle timer with audible indication after 10 seconds	

Order Information

REF	Description	
60034390	Coltolux Comfort Kit, International Contains: Handpiece, battery, charger base, USB-C power sup- ply with plug adapters, 3 intraoral light shields, 3 extraoral light shields, 20 lens caps, 100 barrier sleeves	
60034391	Spare battery	
60034392	Universal power supply with plug adapters	
60034393	Barrier sleeves (pack of 1,000)	
60034394	Intraoral light shield (pack of 5)	
60034395	Extraoral light shield (pack of 5)	
60034396	Lens caps (pack of 20)	
60034446	Composite test black	

Handpiece Dimensions	
Length	8.4 in / 21.5 cm
Diameter	.78 in / 2.0 cm
Weight	2.7 oz / 77g
Light output diameter	11 mm

Warranty	
Curing Light	2 years
Battery	1 year



P/N 31673A

Coltene International Dental Group Dent4You AG Bahnhofstrasse 2 CH-9435 Heerbrugg US Coltène/Whaledent Inc. 235 Ascot Parkway Cuyahoga Falls, OH 44223 Service Center service@coltene.com

www.coltene.com © 2023

