



Sterile HyFlex® CM Files

- › Sterile instruments eliminate the risk of cross contamination between patients, while reducing exposure of the cleaning and disinfection personnel
- › Save processing and handling time
- › Open only the files needed chairside
- › Can be autoclaved for multiple uses

HyFlex CM NiTi files have been manufactured utilizing a unique process that controls the material's memory, making the files extremely flexible, with no re-bound. This enhancement in NiTi file production technology gives the file the ability to more closely follow the anatomy of the canal, and reduces the risk of file separation, binding and complications such as ledging, transportation, and perforation.

Sterile HyFlex[®] CM Files

Catalog Number				Description	Package Size
19 mm	21 mm	25 mm	31 mm		
SBH8190825				.08/25	6 files
	SBH8210415	SBH8250415	SBH8310415	.04/15	6 files
	SBH8210420	SBH8250420	SBH8310420	.04/20	6 files
	SBH8210425	SBH8250425	SBH8310425	.04/25	6 files
	SBH8210430	SBH8250430	SBH8310430	.04/30	6 files
	SBH8210435	SBH8250435	SBH8310435	.04/35	6 files
	SBH8210440	SBH8250440	SBH8310440	.04/40	6 files
	SBH8210445	SBH8250445	SBH8310445	.04/45	6 files
	SBH8210450	SBH8250450	SBH8310450	.04/50	6 files
	SBH8210460	SBH8250460	SBH8310460	.04/60	6 files
	SBH8210620	SBH8250620	SBH8310620	.06/20	6 files
	SBH8210625	SBH8250625	SBH8310625	.06/25	6 files
	SBH8210630	SBH8250630	SBH8310630	.06/30	6 files
	SBH8210635	SBH8250635	SBH8310635	.06/35	6 files
	SBH8210640	SBH8250640	SBH8310640	.06/40	6 files
	SBH821ASST	SBH825ASST	SB60018460	.08/25, .04/20, .04/25, .06/20, .04/30, .04/40	1 of each
	SBH821SML	SBH825SML	SBH831SML	.08/25, .04/35, .04/30, .04/25, .04/20, .04/15	1 of each
	SBH821MED	SBH825MED	SBH831MED	.08/25, .04/60, .04/50, .04/45, .04/40, .04/35	1 of each
	SBH821LRG	SBH825LRG	SBH831LRG	.08/25, .06/40, .06/35, .06/30, .06/25, .06/20	1 of each



Conversion Chart

Catalog Number		Product Description
Sterilized	Current	
SBH8190825	H8190825	HyFlex CM File 08 taper, size 25, 19 mm
SBH8210415	H8210415	HyFlex CM File 04 taper, size 15, 21 mm
SBH8210420	H8210420	HyFlex CM File 04 taper, size 20, 21 mm
SBH8210425	H8210425	HyFlex CM File 04 taper, size 25, 21 mm
SBH8210430	H8210430	HyFlex CM File 04 taper, size 30, 21 mm
SBH8210435	H8210435	HyFlex CM File 04 taper, size 35, 21 mm
SBH8210440	H8210440	HyFlex CM File 04 taper, size 40, 21 mm
SBH8210445	H8210445	HyFlex CM File 04 taper, size 45, 21 mm
SBH8210450	H8210450	HyFlex CM File 04 taper, size 50, 21 mm
SBH8210460	H8210460	HyFlex CM File 04 taper, size 60, 21 mm
SBH8210620	H8210620	HyFlex CM File 06 taper, size 20, 21 mm
SBH8210625	H8210625	HyFlex CM File 06 taper, size 25, 21 mm
SBH8210630	H8210630	HyFlex CM File 06 taper, size 30, 21 mm
SBH8210635	H8210635	HyFlex CM File 06 taper, size 35, 21 mm
SBH8210640	H8210640	HyFlex CM File 06 taper, size 40, 21 mm
SBH821ASST	H821ASST	HyFlex CM File Assorted, 21mm
SBH821LRG	H821LRG	HyFlex CM File Asst. Large , 21mm
SBH821MED	H821MED	HyFlex CM File Asst. Medium , 21mm
SBH821SML	H821SML	HyFlex CM File Asst. Small, 21mm
SBH8250415	H8250415	HyFlex CM File 04 taper, size 15, 25 mm
SBH8250420	H8250420	HyFlex CM File 04 taper, size 20, 25 mm
SBH8250425	H8250425	HyFlex CM File 04 taper, size 25, 25 mm
SBH8250430	H8250430	HyFlex CM File 04 taper, size 30, 25 mm
SBH8250435	H8250435	HyFlex CM File 04 taper, size 35, 25 mm
SBH8250440	H8250440	HyFlex CM File 04 taper, size 40, 25 mm
SBH8250445	H8250445	HyFlex CM File 04 taper, size 45, 25 mm
SBH8250450	H8250450	HyFlex CM File 04 taper, size 50, 25 mm
SBH8250460	H8250460	HyFlex CM File 04 taper, size 60, 25 mm

Catalog Number		Product Description
Sterilized	Current	
SBH8250620	H8250620	HyFlex CM File 06 taper, size 20, 25 mm
SBH8250625	H8250625	HyFlex CM File 06 taper, size 25, 25 mm
SBH8250630	H8250630	HyFlex CM File 06 taper, size 30, 25 mm
SBH8250635	H8250635	HyFlex CM File 06 taper, size 35, 25 mm
SBH8250640	H8250640	HyFlex CM File 06 taper, size 40, 25 mm
SBH825ASST	H825ASST	HyFlex CM File Assorted, 25mm
SBH825LRG	H825LRG	HyFlexCM File Asst. Large , 25mm
SBH825MED	H825MED	HyFlex CM File Asst. Medium , 25mm
SBH825SML	H825SML	HyFlex CM File Asst. Small, 25mm
SBH8310415	H8310415	HyFlex CM File 04 taper, size 15, 31 mm
SBH8310420	H8310420	HyFlex CM File 04 taper, size 20, 31 mm
SBH8310425	H8310425	HyFlex CM File 04 taper, size 25, 31 mm
SBH8310430	H8310430	HyFlex CM File 04 taper, size 30, 31 mm
SBH8310435	H8310435	HyFlex CM File 04 taper, size 35, 31 mm
SBH8310440	H8310440	HyFlex CM File 04 taper, size 40, 31 mm
SBH8310445	H8310445	HyFlex CM File 04 taper, size 45, 31 mm
SBH8310450	H8310450	HyFlex CM File 04 taper, size 50, 31 mm
SBH8310460	H8310460	HyFlex CM File 04 taper, size 60, 31 mm
SBH8310620	H8310620	HyFlex CM File 06 taper, size 20, 31 mm
SBH8310625	H8310625	HyFlex CM File 06 taper, size 25, 31 mm
SBH8310630	H8310630	HyFlex CM File 06 taper, size 30, 31 mm
SBH8310635	H8310635	HyFlex CM File 06 taper, size 35, 31 mm
SBH8310640	H8310640	HyFlex CM File 06 taper, size 40, 31 mm
SB60018460	60018460	HyFlex CM File Assorted, 31mm
SBH831LRG	H831LRG	HyFlex CM File Asst. Large, 31mm
SBH831MED	H831MED	HyFlex CM File Asst. Medium, 31mm
SBH831SML	H831SML	HyFlex CM File Asst. Small 31mm

HyFlex CM Paper Points & Guttapercha Points

HyFlex CM		
Size / Taper	Adhesive-free Paper points 100 pcs	Latex-free Guttapercha points 60 pcs
15/.04	60022671	-
20/.04	60022672	60022689
25/.04	60022673	60022690
30/.04	60022674	60022691
35/.04	60022675	60022692
40/.04	60022676	60022693
45/.04	60022677	60022694
50/.04	60022678	60022695
60/.04	60022679	60022696
Assortment 20-45 /.04	60022680	60022703
20/.06	60022682	60022697
25/.06	60022683	60022698
30/.06	60022684	60022699
35/.06	60022685	60022700
40/.06	60022686	60022701
45/.06	60022687	60022702
Assortment 20-45 /.06	60022688	60022704



Coltene/Wahledent Inc.
235 Ascot Parkway
Cuyahoga Falls, OH 44223 / USA
Tel. USA & Canada + 1 800 221 3046
+ 1 330 916 8800
Fax +1 330 645 8704
info.us@coltene.com

Coltene/Wahledent AG
Feldwiesenstrasse 20
9450 Altstätten / Switzerland
Tel. +41 (0)71 757 53 00
Fax +41 (0)71 757 53 01
info.ch@coltene.com

Coltene/Wahledent GmbH + Co. KG
Raiffeisenstrasse 30
89129 Langenau / Germany
Tel. +49 (0)7345 805 0
Fax +49 (0)7345 805 201
info.de@coltene.com

coltene.com

COLTENE

P/N 31665A