

# BioSonic® UC150

Ultrasonic Cleaning System



Your Patients Trust You. You Can Trust BioSonic.

With its sleek, Swiss design, the **BioSonic UC150** is built for versatility, efficiency and convenience. The 1.5 gallon stainless steel tank offers a small footprint but generous capacity, holding up to 4 exam cassettes. The UC150 can be used as a countertop or recessed unit, allowing flexibility in any sterilization center setting.



### Versatility

1.5 Gallon/5.7L stainless steel tank.

Holds up to 4 exam instrument cassettes.

Can be used as a countertop or recessed unit, allowing for flexibility in any sterilization centre setting.



### Convenience

Quiet operation (63 dB).

Glove-friendly touch screen.

Quick draining with externally accessible drain.



### Efficiency

Data logging capability.

Fits multiple cassettes at once.

Multiple cleaning configurations with basket, beakers, and cassettes.



 **COLTENE**

## Technical Specifications

Tank Capacity	6.6L (1.75 US gal)
Use Capacity	5.7L (1.5 US gal)
Overall Dimensions	34.5 × 27.4 × 35.5 cm (13.6" × 10.8" × 14")
Tank Dimensions	28.7 × 16.0 × 17.8 cm (11.3" × 6.3" × 7")
Operating frequency	48 - 50.5 kHz
Warranty	2 Years

### Maximized Counter Space


The BioSonic UC150 can sit directly on the countertop or be recessed offering ultimate versatility in your work space.



## Accessories



**OPTIM® 28**  
Ultrasonic Cleaning Solution  
SCI2800



**Beaker**  
600 ml beaker with cover & positioning ring.  
UC53



**Beaker holder**  
60022041

Distributed by:  
**SciCan Ltd.**  
1440 Don Mills Rd.,  
Toronto, ON M3B 3P9 / Canada  
Tel: 416 445 1600  
Toll free: 1 800 667 7733  
customerservice@scican.com

Manufactured by:  
**Coltène/Whaledent Inc.**   
235 Ascot Parkway  
Cuyahoga Falls, OH 44223 / USA  
Tel. USA & Canada +1 800 221 3046  
+1 330 916 8800  
Fax +1 330 645 8704  
info.us@coltene.com

[www.scican.com](http://www.scican.com)

 **COLTENE**