

# FLIGHT ULTRASONIC CLEANER




### Simple to Operate

The Flight SonicArc is easy to use right out of the box, so you can get your staff up and running with minimal training time. Features like a digital LCD control panel that displays prompts, options for preset or customized settings, intelligent software that remembers the last cycle settings, and a fluid sensor that ensures the bath always has sufficient water levels help simplify your processes and increase productivity.

- LCD Digital Display
- Multi-function Integration
- High-Precision Temp Control
- Dynamic Light Bar

### SAFETY

Optimize your cleaning process with the Flight Ultrasonic Cleaner, available in four convenient sizes: 3.2L, 6.5L, 10L, and 22L. Designed for dental, medical, veterinary, laboratory, industrial, tattoo, and eyewear applications, this intelligent cleaner features an LCD digital display, high-precision temperature control, and a dynamic light bar for simple operation. Its ultra-powerful ultrasonic technology ensures thorough cleaning while minimizing worker exposure to contaminants and sharps. With a sleek, modern design and advanced functionality, the Flight Ultrasonic Cleaner is a smart, safe, and stylish addition to any workspace.



### Range of sizes and design options

The SonicArc comes in four table top sizes (3.2L, 6.5L, 10L, 22L), So you choose the option that best fits your configuration and workflow. And the stainless steel finish and attractive design coordinate well with any decor.



### ACCESSORIES AND SERVICE PARTS



Feature / Model	UC-6032	UC-6065	UC-6100	UC-6220
Tank Capacity gal/L	0.85Gal/3.2L	1.7Gal/6.5L	2.64Gal/10L	5.8Gal/22L
Ultrasonic Frequency	40 kHz	40 kHz	40 kHz	40 kHz
Heater Option	Yes	Yes	Yes	Yes
Timer Control	Yes	Yes	Yes	Yes
Drainage Valve	No	Yes	Yes	Yes
Warranty	1 Year	1 Year	1 Year	1 Year
Heating Power	200W	300W	400W	600W
Ultrasonic Power	120/60W	180/90W	240/120W	480/240W
Tank Size Dimensions (W x D x H)	9.5 x 5.3 x 3.9"	11.8 x 5.9 x 5.9"	11.8 x 11.4 x 11"	19.7 x 11.8 x 5.9"
Product Size (W x D x H)	10.6 x 7.3 x 9"	14.8 x 7.9 x 10.8"	15 x 11.4 x 11"	23 x 13 x 11"
Weight	7.9lbs (3.6kg)	14.1lbs (6.4kg)	16.8lbs (7.6kg)	29.3(13.3kg)
Noise Level	<60db	<60db	<60db	<60db