


BRILLIANT Bulk Fill Flow



One Multi-Shade Composite that blends brilliantly. BRILLIANT Bulk Fill Flow is the ideal filling material for fast and easy restorations.



No top layer needed

Thanks to high abrasion resistant properties, quality restorations are achieved in one efficient step since there is no need for a covering layer.

4 mm reliable curing depth

Increments up to 4 mm thickness are reliably cured in 20 sec ($\geq 1000 \text{ mW/cm}^2$). This feature combined with the surface wettability and flowability of the material offers time saving working steps.

One adaptive Multi Shade material

The chameleon effect of BRILLIANT Bulk Fill Flow allows the Multi Shade material to adapt to a wide range of tooth shades.




Courtesy: Dr. Lndon Albrecht

60033280 BRILLIANT Bulk Fill Flow Refill
1 syringe 2.3 g + 11 application needles

High flexural strength, low shrinkage and low abrasion allow for fillings in the posterior region even without an additional covering layer.


FLEXURAL STRENGTH / MPA



Among the tested bulk composites, **BRILLIANT** Bulk Fill Flow is one of the materials with high flexural strength.


BRILLIANT Bulk Fill Flow exhibits comparatively low shrinkage compared to competitor products.

SHRINKAGE / %



BRILLIANT Bulk Fill Flow shows the lowest abrasion-related material loss in the investigated material group.

ABRASION / μm



Clinical procedure →



Coltène/Whaledent Inc.
 235 Ascot Parkway
 Cuyahoga Falls, Ohio 44223 / USA
 T +1 330 916 8800
 F +1 330 645 8704
 info.us@coltene.com

Coltène/Whaledent AG
 Feldwiesenstr. 20
 9450 Altstätten / Switzerland
 T +41 71 757 53 00
 F +41 71 757 53 01
 info.ch@coltene.com

