

V10

CELEBRATING

25

YEARS OF  
EXCELLENCE

# SABLE PRODUCT CATALOGUE



**SABLE**  
INDUSTRIES INC.

[WWW.SABLEINDUSTRIESINC.COM](http://WWW.SABLEINDUSTRIESINC.COM)

### HANDPIECES

High Speed Handpieces and Couplers..... 1-2  
Sable NSK & STAR Cross Reference Chart.....3  
Slow Speed Handpieces & Components ..... 4-7  
MIDWEST Type Slow Speed Components.....8  
U-Type Prophy Heads & Contra Angles.....9  
Sable KAVO & MIDWEST Cross Reference Chart .....10  
KAVO Type Premium Slow Speed Handpieces .....11  
Handpiece Repair Tools and Parts .....12  
Sable EZ Lube Spray .....13

### HANDPIECE REPLACEMENT PARTS

Sable 1:5 Electric Handpiece Parts  
KAVO 25LHA, KAVO LPA/LPR, NSK Electric Parts .....14-16  
Replacement Turbines by Brand.....17-27  
Bearing Table (Steel and Ceramic) ..... 28  
NSK & OEM Canister Chuck Type Turbines ..... 29  
NSK & OEM Canister PB Type Turbines ..... 30  
NSK Type Economy Canister Turbines..... 31-33  
NSK & OEM Back Caps ..... 33  
Handpiece Replacement Bulbs ..... 34  
LED Fiber Optic Handpiece Replacement Bulbs..... 35

### EQUIPMENT

Light Control..... 36  
Handpiece Tubing ISO Fiber Optic, EURO, ADEC.....37  
Buck-Boost Transformers .....38-39  
Autoscaler & Parts ..... 40  
Ultrasonic Scaling Tips ..... 40  
Statim Compressors & Autoclave Gaskets..... 40

### BIO-PURE

Bio-Pure Evacuation System Cleaner..... 41-42  
Bio-Pure EPA Compliance Plus Software (USA Only) .... 42  
Ultrasonic 3-in-1 Presoak & Cleaner..... 43  
Product Warranty List ..... 44  
Product Returns, Pricing, Conditions, Terms ..... 45

The following are registered trademarks of the indicated companies:  
A-dec - A-dec Inc., Bien Air - Bien Air Dental SA, Kavo Kaltenbach & Voight GmbH & Co., Sirona - Sirona Dental System GmbH, W&H - W&H Dentalwerk Burmose Inc., NSK - NSK Nakanishi Inc., Star - Dental EZ Group, DENTSPLY International Inc (Midwest, Cavitron, RDH, Bobcat, SPS), Mk-Dent - Mk-Dent GmbH, Faro - Faro S.p.A., Statim 2000 - SciCan Ltd., Whaledent - Coltene/Whaledent Inc., Hu-Friedy- Hu Friedy Mfg. Co., LLC, Premier- Premier Products Co., Impact Air - Palisades Dental, LLC, Siemens - Siemens AG, Vector - Vector Research and Development Inc., Safo - Safo Dental Supply Inc., Dabi Atlante Industrias S/A, Pelton & Crane - Pelton & Crane. Sable Industries Inc. is in no way affiliated with the above companies.

### SABLE PRO SERIES FIBER OPTIC HANDPIECES

- 2 year warranty
- Quadruple Water Spray
- Back cap design to remove exhaust air more efficiently

ACCESS™ PRO-KL Mini Head **2000021**

ROTAMAX™ PRO-KL Torque Head **2000022**



Refer to page 2 for Compatible Couplers

### ROTAMAX 45° DEGREE HANDPIECES

ROTAMAX™ 45-KL Torque Head **2000023**

ROTAMAX™ 45 Torque Head **2000024**

Features a Fixed Backend



### NON-FIBER OPTIC HIGH SPEED HANDPIECES

SABLE ACCESS™ II **2000003**  
PB Miniature Head without coupler  
12 Month Warranty

SABLE ROTAMAX™ II **2000004**  
PB Torque Head without coupler  
12 Month Warranty

4 HOLE / 360° **2000101**  
Quick Connect Coupler



4 Hole / 360° / KAVO MULTIflex Type Coupler  
**2000114**

KAVO MODEL# 457 OEM P/N: 553-1210



5 Hole / 360° / KAVO MULTIflex Type Coupler  
**2000105**

KAVO MODEL#1390 P/N 0553.1390



6 Hole / 360° / KAVO MULTIflex LUX Type Coupler  
**2000116**

KAVO MODEL# 465LRN P/N: 1001.7599



6 Hole / 360° / KAVO MULTIflex LUX Type LED Coupler  
**2000119**

KAVO MODEL# 460LED P/N: 1007.3201



6 Hole NSK Type NL400 Fiber Optic Coupler  
**2000107**

NSK MODEL# PTL-CL-4HV



6 Hole NSK Type NL400 LED Fiber Optic Coupler  
**2000108**



**NSK Type** - Custom Options Available

Please Specify:

6 Hole NSK Type NL400 FIBER OPTIC Coupler **2000107**

6 Hole NSK Type NL400 LED Coupler **2000108**

4 Hole NSK Type NL400 Non-Optic Coupler **Specify Custom 4 hole**

**WARRANTY:**

All SABLE Couplers carry a 1-Year Warranty from date of Installation.

Product Description	Sable Part #	NSK Part
---------------------	--------------	----------

**DOCTORS PACKAGE**

20K 4-Hole Air Motor, 1:1 Straight Nosecone, 1:1 Contra Angle, Standard Swing Latch Head	<b>1600008</b>	EX-203 Y100-104 SET
20K Air Motor, 1:1 Nosecone, 1:1 Contra Angle, Standard Latch Head	<b>1600009</b>	EX-203 Y100-104 SET

**AIR MOTORS – 4 HOLE, ISO E-TYPE**

Accept all brands of ISO E-type Attachments

20K Air Motor	<b>1600110</b>	EX-203 M138
5K Air Motor	<b>1600102</b>	EX-103 M147
20K Air Mini Motor	<b>1600111</b>	FX-204

**ATTACHMENTS**

1:1 Contra Angle (Blue Band)	<b>1600210</b>	EC C100
4:1 Contra Angle (Green Band)	<b>1600202</b>	E-4R C220
16:1 Contra Angle	<b>1600203</b>	E-16R C219
1:1 Straight Nosecone	<b>1600310</b>	EX-6 M116
4:1 Straight Nosecone	<b>1600303</b>	EX-5 M115

**HEADS - ISO E-TYPE, 12 TOOTH (NSK)**

Snap-on Prophy Head	<b>1600401</b>	AR-YK C188
Screw-in Prophy Head	<b>1600402</b>	AR-YS C187
Standard Latch Head	<b>1600403</b>	NAC-Y C032
Ball Bearing Latch Head	<b>1600407</b>	BB-YL C035
Miniature Latch Head	<b>1600405</b>	NM-Y C247
Push Button Ball Bearing Latch Head	<b>1600421</b>	FPB-Y C875 NPB-Y C808

**SLOW SPEED ONE PIECE HEADS**

E-Type 1:1 BB PB Latch Contra- Angle	<b>1600214</b>	FX 25 C1052
E-Type 4:1 Screw in Prophy Contra- Angle	<b>1600216</b>	FX 57 C1054
E-Type 4:1 BB PB Latch Contra- Angle	<b>1600218</b>	FX15 C1053

**SABLE / STAR CROSS REFERENCE CHART**

Product Description	Sable Part #	NSK Part
Star Titan 5 Star Snap-On Head	<b>1600408</b>	263176
Star Titan 5 Star Screw-In Head	<b>1600409</b>	263176
Star Titan Standard Latch Head	<b>1600406</b>	n/a
Star Titan Ball Bearing Latch Head	<b>1600410</b>	259005

**DOCTOR'S PACKAGES**

**Standard Swing Latch Head Sable NSK Style 1600008**

NSK MODEL# EX-203 Y100-104 SET 1600403 Standard Swing Latch Head  
1600210 EC Type 1:1 Contra Angle  
1600310 EX-6 Type 1:1 Straight Nosecone  
1600110 EX-203 Type 20K 4-Hole Air Motor



**Push Button BB Latch Head Sable NSK Style 1600009**

1600421 Push Button BB Latch Head  
1600210 EC Type 1:1 Contra Angle  
1600310 EX-6 Type 1:1 Straight Nosecone  
1600110 EX-203 Type 20K 4-Hole Air Motor



**Push Button BB Latch Head Sable KaVo Style 1600010**

**\*with Internal Air & Water** 1610810 KAVO Type 1:1 Ball Bearing PB Latch Head  
1610210\* 1:1 Contra Angle  
1610310\* 1:1 Straight Nosecone  
1610110\* IS-205 Type 20K 4-Hole Motor



**HYGIENIST PROPHY PACKAGES**

**Hygienist Package #1 - Snap-On E-Type 1600011**

1600110 20K Air Motor, 1600202 4:1 Reduction Contra Angle,  
1600401 Snap-on Sealed Propphy Head



**Hygienist Package #2 - Screw-In E-Type 1600012**

1600110 20K Air Motor, 1600202 4:1 Reduction Contra Angle,  
1600402 Screw-in Sealed Propphy Head



**Hygienist Package #3 - Standard Latch E-Type 1600013**

1600110 20K Air Motor, 1600202 4:1 Reduction Contra Angle,  
1600403 Standard Latch Head



**Hygienist Package #4 - Ball Bearing Swing Latch E-Type 1600014**

1600110 20K Air Motor, 1600202 4:1 Reduction Contra Angle,  
1600407 Ball Bearing Swing Latch Head



Package **1600023**  
(1600214 + 1600110)  
1:1 BB PB Latch



Package **1600024**  
(1600215 + 1600110)  
BB PB Friction Grip



Package **1600025**  
(1600216 + 1600110)  
4:1 Screw-In Propphy



Package **1600026**  
(1600218 + 1600110)  
4:1 BB PB Latch Propphy



**ACCESS™ PROPHY CLEAN**

Sable 5K U-Type Hygienist Propphy Swivel Handpiece  
High torque. 360° Swivel nosecone. Accepts all hygienist  
disposable and metal prophy angles.  
Friction grip style chuck.

Weights 2.8 oz.  
Fully autoclavable.



**1600005**

**AIR MOTORS - 4 Holes, ISO E-Type.**

Accepts all brands of ISO E-Type attachments.

20K Air Motor NSK MODEL# EX-203 M138

**1600110**



20K Mini Air Motor NSK MODEL# FX-204

**1600111**



5K Air Motor NSK MODEL# EX-103 M147

**1600102**



**NOSECONES**

1:1 Straight Nosecone NSK MODEL# EX-6 M116

**1600310**



4:1 Straight Reduction Nosecone NSK MODEL# EX-5 M115

**1600303**



**ATTACHMENTS** - ISO E-Type 12 Tooth NSK Type,  
Couple to all brands of ISO E-Type motors.

1:1 Contra Angle  
**1600210**  
NSK MODEL# EC C100



4:1 Reduction Contra Angle  
**1600202**  
NSK MODEL# E-4R C220



Sable E-Type 1:1 BB PB  
Latch Contra Angle  
One Piece Design  
**1600214**  
FX 25 2.35mm Burs



Sable E-Type 1:1 PB  
Friction Grip Contra Angle  
One Piece Design  
**1600215**  
1.6mm Burs



Sable E-Type 4:1 Screw-in  
Prophy Contra Angle  
One Piece Design  
**1600216**  
FX 57



Sable E- Type 4:1 Snap-on  
Prophy Contra Angle  
One Piece Design  
**1600217**



Sable E-Type 4:1 PB Latch  
Prophy Contra Angle  
One Piece Design  
**1600218**  
FX 15 2.35mm Burs



**HEADS** - ISO E-Type 12 Tooth (SABLE NSK & STAR TYPE)

Snap-On Prophy Head 5,000 RPM  
**1600401** NSK MODEL# AR-YK C188



Screw-In Prophy Head 5,000 RPM  
**1600402** NSK MODEL# AR-YS C187



Standard Swing Latch Head 5,000 RPM  
**1600403** NSK MODEL# NAC-Y C032



Ball Bearing Swing Latch Head 40,000 RPM  
**1600407** NSK MODEL# BB-YL C035



Miniature Latch Head 5,000 RPM  
**1600405** NSK MODEL# NM-Y C247



Ball Bearing Push Latch Head 40,000 RPM  
**1600404**



Ball Bearing Push Button Latch Head  
40,000 RPM  
**1600421** NSK MODEL# FPB-Y C873, NPB-Y C808



**STAR TYPE HEADS W/ ELBOW**

STAR Type Ball Bearing Push Button  
Latch Head W/ Elbow 20,000 RPM  
**1600412\*** (\*1600421 with Elbow)



STAR Type Standard Latch W/ Elbow  
5,000 RPM  
**1600406**



STAR Type Snap-On  
Prophy Head W/ Elbow 5,000 RPM  
**1600408**



STAR Type Screw-In  
Prophy Head W/ Elbow 5,000 RPM  
**1600409** STAR MODEL# 263176



STAR Type Ball Bearing Swing Latch  
Head W/ Elbow 20,000 RPM  
**1600410** STAR MODEL# 259005



STAR Type 12-Tooth Elbow Only  
**1604004**



## MIDWEST TYPE SLOW SPEED COMPONENTS

1:1 MIDWEST Shorty & Rhino Type Contra Angle  
**1600204**



MIDWEST MODEL# 710074

1:1 MIDWEST Shorty & Rhino Type Nosecone  
**1600302**



MIDWEST MODEL# 720005

MIDWEST Type 14 Tooth Snap-On Prophy Head  
**1600501**



MIDWEST Type 14 Tooth Screw-In Prophy Head  
**1600502**



MIDWEST Type 14 Tooth Standard Latch Head  
**1600503**



MIDWEST Type 14 Tooth Ball Bearing Push Button Latch Head  
**1600504**



MIDWEST Type 14 Tooth Ball Bearing Swing Latch Head  
**1600507**



MIDWEST MODEL# 720406

WARRANTY: 6 Months on Attachments, 3 Months on Heads.

## U-TYPE PROPHY HEADS & CONTRA ANGLES

### PROPHY HEADS

U-Type Snap-On Prophy Head, 5,000 RPM  
**1300103**



U-Type Screw-In Prophy Head, 5,000 RPM  
**1300104**



### CONTRA ANGLES W/ HEAD

U-Type Standard Latch, 5,000 RPM Non-Sectional  
**1300301**



U-Type Standard Latch, 5,000 RPM Sectional  
**1300302**



U-Type Ball Bearing Push Latch, 20,000 RPM Sectional  
**1300303**



U-Type Friction Grip W/ Bur Tool, 20,000 RPM Sectional  
**1300304**



WARRANTY: 6 Months on Contra Angles, 3 Months on Heads.



## SABLE / KAVO CROSS REFERENCE CHART

Product Description	Sable Part #	Interchangeable All Features	Interchangeable Most Features
<b>SLOW SPEED</b>			
20K Air Motor w/ Internal Air & Water Non-Fiber Optic	<b>1610110</b>	KAVO 181H OEM P/N: 535-5200	
1:1 Contra Angle w/ Internal Air & Water Non-Fiber Optic	<b>1610210</b>	KAVO 20C 20 Series OEM P/N: 534-5220	KAVO 20, 20A, 0B, 20CH, 20CHC, 20CN
1:1 Nosecone w/ Internal Air & Water Non-Fiber Optic	<b>1610310</b>	KAVO 10C 10 Series OEM P/N: 524-1070	KAVO 10, 10A, 10B, 10CH, 10CN, 10CNC
1:1 Ball Bearing Push Button Latch Head	<b>1610810</b>	KAVO 68LH 68 Series OEM P/N: 540-2040	68LDN/68LU
1:1 Screw-In Propy Head	<b>1610820</b>	KAVO 31LR 31 Series OEM P/N: 549-338	
<b>HIGH SPEED</b>			
Access PRO-KL Mini Head	<b>2000021</b>	8700	642, 635, 645, 4500B, 4500BR, 8700
Rotamax PRO-KL Torque Head	<b>2000022</b>	8900	640, 647, 649, 659, 6500, 6500BR, 8900
KAVO type MULTiflex Coupler (Non-Fiber Optic)	<b>2000114</b>	457 OEM P/N: 553-1210	
KAVO type MULTiflex LUX Coupler (Fiber Optic)	<b>2000116/2000119</b>	465LRN OEM P/N: 553-1550	
KAVO type Bulb (3-Pack)	<b>1000802</b>	OEM P/N: 553-3881	
<b>FIBER OPTIC SYSTEMS</b>			
Light Control Board	<b>2300101</b>	OEM P/N: 905-8400	
6 Hole 7' Grey F/O Tubing	<b>2100103</b>	OEM P/N: KLO35009	
6 Hole 10' Grey F/O Tubing	<b>2100105</b>	OEM P/N: KLO35010	

## SABLE / MIDWEST CROSS REFERENCE CHART

Product Description	Sable Part #	Midwest Part #
Sable MIDWEST Rhino/Shorty 1:1 Type Contra Angle	<b>1600204</b>	710074
Sable MIDWEST Rhino/Shorty 1:1 Nosecone	<b>1600302</b>	720005
Sable MIDWEST Standard Latch Head	<b>1600503</b>	n/a
Sable MIDWEST BB Push Latch Head	<b>1600504</b>	n/a
Sable MIDWEST BB Swing Latch Head	<b>1600507</b>	720406

## KAVO TYPE PREMIUM SLOW SPEED HANDPIECES

**AIR MOTORS** - ISO E-TYPE. Couple to all brands of ISO E-Type attachments (with or without air & water).

20K Air Motor W/ Internal Air & Water 20,000 Rpm. Accepts All ISO E-Type Attachments

**1610110**

KAVO MODEL#  
KAVO 181H OEM P/N: 535-5200



**ATTACHMENTS** - ISO E-TYPE. Couple to all brands of ISO E-Type air or electric motors (with or without air & water).

1:1 Ball Bearing Contra Angle W/ Internal Air & Water 40,000 RPM Accepts All KAVO & Sable Premium Type Heads

**1610210**

KAVO MODEL#  
KAVO 20C OEM P/N: 534-5220



1:1 Ball Bearing Straight Nosecone W/ Internal Air & Water 40,000 RPM

**1610310**

KAVO MODEL#  
KAVO 10C OEM P/N: 524-1070



**HEADS** - Attach to SABLE and KAVO Contra Angles.

Ball Bearing 1:1 Push Button Latch Head

**1610810**

KAVO MODEL#  
KAVO 68LH OEM P/N: 540-2040



1:1 Screw-In Propy Head  
Maximum Speed 5,000 RPM  
(Available only in Canada)

**1610820**



WARRANTY: 1 Year on all KAVO Type Premium Handpieces.

SABLE ROTAMAX™ ACCESS™ I & II  
Lubricating Nozzle



**2001220**

Lubrication Nozzle

(KAVO Type Contra-Angles & Heads)



**1611201**

KAVO MULITflex Lubricating Nozzle



**2001228**

NSK Lubricating Nozzle

**2001227**

NSK Coupler O-Ring Set

**2001232**

KAVO MULITflex Coupler O-Ring Set



**2001231**

Coupler O-Ring Set For 2000101 & 2000102

**2001217**

Coupler O-Ring For 1600363 (Complete Set)

**2001223**

E-Type Motor O-Ring 3 pack

**2001233**

Coupler Rear Gasket (4-Hole)



**2001218**

6-Hole Coupler Rear Gasket



**2001224**

SABLE ACCESS™ I & II Cap Key



**2001215**

SABLE ROTAMAX™ I & II Cap Key

**2001216**

SABLE ACCESS™ S4500KL /

ROTAMAX™ S6500KL Type Cap Key



**2001234**

Clean Out Wire

**2001230**

**SPRAY, DROPPER & REFILL  
HANDPIECE LUBRICANT AND CLEANER**



- Cleans & Maintains as it Lubricates.
- Synthetic Food Grade Aerospace Industry Level Lubricant.
- Multi-Use Dental Handpiece Lubricant, will not breakdown in Sterilizer Temperatures

Sable EZ Lube SPRAY 500ml (sold only as case of 6) **2200103**

Sable EZ Lube REFILL 500ml (sold only as case of 6) **2200104**

Sable EZ Lube DROPPER 60ml (sold only as case of 6) **2200105**

E-TYPE Contra Angles, Straight Attachments  
for KAVO Heads Lubricating Nozzle **1611201**



NSK High Speed Lubricating Nozzle **2001227**



KAVO Multiflex Lubricating Nozzle **2001228**





**KAVO 25LHA**

KAVO Type 25LHA Gear Insert Cartridge  
**1608051**  
OEM Part # 0.553.7341



KAVO Type 25LHA Intermediate Shaft  
**1608055**  
OEM Part # 0.553.7201



KAVO Type 25LHA  
Cartridge & Shaft Kit  
**1608056**  
OEM Part # 0.553.7341 + 0.553.7201



KAVO Type 25LHA  
Cartridge, Shaft & Idler Gear Kit  
**1608057**  
OEM Part # 0.553.7341 +  
0.553.7201 + 0.553.7251



KAVO Type 25LHA  
5-Pack Front Gear Intermediate Shaft  
**1608059**

KAVO Type 25LHA  
5-Pack Rear Gear Intermediate Shaft  
**1608060**

KAVO Type M25L  
Gear Insert, Intermediate Shaft & Idler Gear Kit  
**1608079**

**KAVO LPA/LPR**

KAVO Type 25LPA/LPR Gear Insert Cartridge  
**1608052**  
OEM Part # 1.001.3520



KAVO Type 25LPA/LPR  
Intermediate Shaft  
**1608053**  
OEM Part # 1.001.3515



KAVO Type 25LPA/LPR  
Cartridge & Shaft Kit  
**1608054**  
OEM Part # 1.0002.6135 + 1.001.3515



KAVO Type 25LPA/LPR  
Cartridge, Shaft & Idler Gear Kit  
**1608058**  
OEM Part # 1.0002.6135 +  
1.001.3515 + 1.001.3519 or 1.006.2965



KAVO Type 25LPA/LPR Idler Shaft  
**1608064**

KAVO Type 25LPA/LPR  
5-Pack Front Gear Intermediate Shaft  
**1608065**

KAVO Type 25LPA/LPR  
5-Pack Rear Gear Intermediate Shaft  
**1608066**

KAVO Type 25LPA/LPR  
Gear Insert Cart. + Front Gear Int. Shaft  
**1608067**

KAVO Type E25L  
Intermediate Shaft (OEM)  
**1608072**



KAVO Type E25L  
Gear Insert Cartridge + Intermediate Shaft (OEM)  
**1608073**



KAVO Type E25L  
Gear Insert Cartridge  
**1608077**



KAVO Type M25L  
Gear Insert Cartridge  
**1608078**

WARRANTY: 1 Year.

**NSK ELECTRIC PARTS**

NSK Type Ti-Max X95L  
Gear Insert Cartridge  
**1604050**



NSK Type Ti-Max X95L  
Intermediate Shaft (OEM)  
**1604051**



NSK Type Ti-Max X95L  
Gear Insert Cartridge + Int. Shaft (OEM)  
**1604052**



NSK Type Ti-Max Z95L Head Cartridge  
**1604060**

W&H Synea WA-99LT 3 Port Spray Head Cartridge  
**1607060**

W&H Synea WA-99LT 3 Port Spray Intermediate Shaft  
**1607061**

W&H Synea WA-99LT 3 Port Spray Cartridge & Shaft Kit  
**1607062**

W&H Synea WA-99LT 5 Port Spray Head Cartridge  
**1607070**

W&H Synea WA-99LT 5 Port Spray Intermediate Shaft  
**1607071**

W&H Synea WA-99LT 5 Port Spray Cartridge & Shaft Kit  
**1607072**

**KAVO 632,  
633, 639,  
642, 643B,  
645 Type**  
Turbine Steel 1000101  
Ceramic 1000137  
Conversion Kit 1001101  
Rebuild Kit 1000701  
Spindle/Impeller 1000501  
Spindle 1000201  
Steel Bearing 1000601M  
Ceramic Bearing 1000699M  
O-Ring(s) 1000921  
Back Cap 1000401  
Back Cap (643B & 645) 1000423  
Cap Tool 1001205



**KAVO 625C  
630/40  
Type**  
Turbine Steel 1000102  
Ceramic 1000138  
Conversion Kit Steel 1001102  
Ceramic 1001103  
Rebuild Kit 1000701  
Spindle/Impeller 1000502  
Spindle 1000202  
Steel Bearing 1000601M  
Ceramic Bearing 1000699M  
O-Ring(s) 1000921  
Back Cap 1000402-640  
1000403-625/630  
Cap Tool 1001205



**KAVO  
635B  
Type**

Turbine Steel 1000115  
Ceramic 1000139  
Rebuild Kit 1000715  
Spindle/Impeller 1000517  
Spindle 1000217  
Steel Bearing 1000601M (Front)  
1000613M (Rear)  
1000699M (Front)  
1000695M (Rear)  
Ceramic Bearing  
O-Ring(s) 1000921  
Air Dam 1001211  
Back Cap 1000415  
Cap Tool 1001217



**KAVO  
650  
Type**

Turbine Steel 1000116  
Ceramic 1000117  
Rebuild Kit 1000701  
Spindle/Impeller 1000516  
Spindle 1000218  
Steel Bearing 1000601M  
Ceramic Bearing 1000699M  
O-Ring(s) 1000921  
Air Dam 1001213  
Back Cap 1000422  
Cap Tool 1001214



**KAVO  
635BR  
Type**

Turbine Ceramic 1000188  
Spindle/Impeller 1000517  
Spindle 1000217  
Ceramic Bearing 1000699M (Front)  
1000695M (Rear)  
O-Ring(s) 1000921  
Air Dam 1001219  
Back Cap 1000415  
Cap Tool 1001217



**KAVO  
6500BR  
Type**

Turbine Ceramic 1000189  
Spindle/Impeller 1000563  
Spindle 1000234  
Ceramic Bearing 1000697M (Front)  
1000698M (Rear)  
O-Ring(s) 1000921  
Air Dam 1001236  
Spring Washer(s) 1001223  
Back Cap 1000442  
Cap Tool 1001214



# REPLACEMENT TURBINES BY BRAND

## KAVO 647B, 649B Type



Turbine Steel	1000122
Ceramic	1000140
Rebuild Kit	1000701
Spindle/Impeller	1000522
Spindle	1000219
Steel Bearing	1000601M
Ceramic Bearing	1000699M
O-Ring(s)	1000921
Air Dam	1001212
Back Cap	1000422
Cap Tool	1001214

## KAVO 6000 Type



Turbine Steel	1000125
Ceramic	1000161
Rebuild Kit	1000733
Spindle/Impeller	1000531
Spindle	1000231
Steel Bearing	1000608M (Front)
	1000605M (Rear)
	1000693M (Front)
	1000698M (Rear)
Ceramic Bearing	1000921
O-Ring(s)	1001223
Spring Washer	1001224
Back Cap	1000425
Cap Tool	1001214

## KAVO 659B, 659BR Type



Turbine Steel	1000175
Ceramic	1000170
Rebuild Kit Steel	1000736
Ceramic	1000731
Spindle/Impeller	1000561
Spindle	1000236
Steel Bearing	1000613M (Front)
	1000617M (Rear)
Ceramic Bearing	1000695M (Front)
	1000690M (Rear)
O-Ring(s)	1000921
Spring Washer(s)	1001223
Back Cap	1000440
Cap Tool	1001214

## KAVO 7000 Type



Turbine Ceramic	1000126
Rebuild Kit Ceramic	1000762
Spindle/Impeller	1000532
Spindle	1000233
Steel Bearing	1000605M
Ceramic Bearing	1000698M
O-Ring(s)	1000921
Spring Washer(s)	1001223
Back Cap	1000425
Cap Tool	1001214

## KAVO 655, 660B Type



Turbine Steel	1000131
Ceramic	1000132
Rebuild Kit Ceramic	1000760
Spindle/Impeller	1000533
Spindle	1000232
Steel Bearing	1000611M
Ceramic Bearing	1000697M
O-Ring(s)	1000921
Spring Washer(s)	1001223
Back Cap	1000425
Cap Tool	1001214

## KAVO 634 Type



Turbine Steel	1000160
Ceramic	1000141
Rebuild Kit	1000701
Spindle/Impeller	1000501
Spindle	1000201
Steel Bearing	1000601M
Ceramic Bearing	1000699M
O-Ring(s)	1000921
Back Cap	1000401
Cap Tool	1001205

## KAVO 637B Ceramic Type



Turbine Ceramic	1000165
Rebuild Kit Ceramic	1000761
Spindle/Impeller	1000517
Spindle	1000217
Steel Bearing	1000601M (Front)
	1000613M (Rear)
Ceramic Bearing	1000699M (Front)
	1000695M (Rear)
O-Ring(s)	1000921
Air Dam	1001219
Back Cap	1000415
Cap Tool	1001217

## KAVO 8000 Ceramic Type



Turbine Ceramic	1000166
Rebuild Kit Ceramic	1000763
Spindle/Impeller	1000563
Steel Bearing	1000611M (Front)
	1000605M (Rear)
Ceramic Bearing	1000697M (Front)
	1000698M (Rear)
O-Ring(s)	1000921
Spring Washer(s)	1001223
Back Cap	1000442
Cap Tool	1001214

## KAVO 6500B Type



Turbine Steel	1000173
Ceramic	1000178
Rebuild Kit Steel	1000735
Ceramic	1000738
Spindle/Impeller	1000563
Spindle	1000234
Steel Bearing	1000608M (Front)
	1000605M (Rear)
Ceramic Bearing	1000693M (Front)
	1000698M (Rear)
O-Ring(s)	1000921
Air Dam	1001236
Spring Washer	1001223 (Front)
	1001224 (Rear)
Back Cap	1000442
Cap Tool	1001214

## KAVO 625CD 636, 646 Type



Turbine Steel	1000174
Ceramic	1000179
Rebuild Kit Steel	1000701
Ceramic	1000739
Spindle/Impeller	1000564
Spindle	1000235
Steel Bearing	1000601M
Ceramic Bearing	1000699M
O-Ring(s)	1000921
Back Cap	1000443
Cap Tool	1001214

## KAVO 4500B/BR, M4500L Ceramic Type



Turbine Ceramic	1000171
Rebuild Kit Ceramic	1000732
Spindle/Impeller	1000562
Steel Bearing	1000613M
Ceramic Bearing	1000695M
O-Ring(s)	1000921
Spring Washer(s)	1001224
Back Cap	1000441
Cap Tool	1001234

## KAVO SMARTtorque 609C & S619L



Turbine Steel	1000191
Steel Bearing	1000611M
Ceramic Bearing	1000697M
O-Ring(s)	1000921
Spring Washer	1001223 (Front)
	1001223 (Rear)
Back Cap	1000425
Cap Tool	1001214

# REPLACEMENT TURBINES BY BRAND

## KAVO SMARTtorque Mini S605C & S615L



Turbine Steel	1000192
Steel Bearing	1000617M
Ceramic Bearing	1000690M
O-Ring(s)	1000921
Spring Washer(s)	1001223
Cap Tool	1001234

## KAVO MASTERtorque LUX M8700L Type



Turbine Ceramic	1000198
Ceramic Bearing	1000695M (Front)
	1000697M (Rear)
Spindle/Impeller	1000588
O-Ring (s)	1000921 (Front)
	1000921 (Rear)
Spring Washer (s)	1001224 (Front)
	1001223 (Rear)
Push Button Back Cap	1000485
Push Button Cap Tool	1001234

## KAVO E679L, E680L Turbine



Ceramic Turbine	1000194
Spindle/Impeller	1000586
Ceramic Bearing	1000697M (Front)
Ceramic Bearing	1000695M (Rear)
O-Ring (s)	1000921 (x2)
Spring Washer Front	1001224
Spring Washer Rear	1001223
Push Button Back Cap	1000483
Push Button Cap Tool	1001214

## MIDWEST Tradition XGT Type



Turbine Steel	1000110
Conversion Kit	1001110
Rebuild Kit	1000717
Spindle/Impeller	1000510
Spindle	1000210
Steel Bearing	1000604M
Ceramic Bearing	1000694M
O-Ring(s)	1000903
Air Dam	1001220
Back Cap	1000410
Cap Tool	1001243

## KAVO EXPERTtorque Mini 675L, 677L



Turbine Ceramic	1000197
Ceramic Bearing	1000695M (Front)
	1000690M (Rear)
O-Ring(s)	1000921
Spring Washer	1001224 (Front)
	1001223 (Rear)
Cap Tool	1001234
Back Cap	1000441

## KAVO MASTERtorque LUX M8900L Type



Turbine Ceramic	1000199
Ceramic Bearing	1000695M (Front)
	1000697M (Rear)
Spindle/Impeller	1000587
O-Ring (s)	1000921 (Front)
	1000921 (Rear)
Spring Washer (s)	1001223 (Front)
	1001223 (Rear)
Push Button Back Cap	1000484
Push Button Cap Tool	1001214
Turbine Brake	1001247
Turbine Brake Tool	1001248

## J. Morita TWINPOWER 4-Hex & 45 Head Type



Turbine Ceramic	1000196
Spindle/impeller	1000585
Ceramic Bearing	1000698M
O-Ring	1000906
Spring Washer	1001223
Note: Canister lid, Canister and External O-Rings must be re-used	

## MIDWEST Quiet Air Type



Turbine Steel	1000112
Conversion Kit	1001111
Rebuild Kit	1000704
Spindle/Impeller	1000515
Spindle	1000216
Steel Bearing	1000604M (Front)
	1000605M (Rear)
Ceramic Bearing	1000694M (Front)
	1000698M (Rear)
O-Ring(s)	1000902 (Front)
	1000903 (Rear)
Back Cap	1000411
Cap Tool	1001243

## MIDWEST Stylus, 540S, 541S, 360S, 361S, Type



Turbine Steel	1000153
Rebuild Kit	1000724
Spindle/Impeller	1000567
Steel Bearing	1000614M
Ceramic Bearing	1000689M
O-Ring(s)	1000914
Air Dam	1001218
Back Cap	1000424
Cap Tool	1001253

## MIDWEST Stylus ATC 990 Mini Type



Turbine Steel	1000177
Ceramic	1000182
Rebuild Kit Steel	1000737
Ceramic	1000740
Spindle/Impeller	1000566
Spindle	1000238
Steel Bearing	1000611M
Ceramic Bearing	1000697M
O-Ring(s)	1000921
Spring Washer(s)	1001223
Back Cap	1000444
Cap Tool	1001253

## MIDWEST Tradition Lever Type



Turbine Steel	1000183
Rebuild Kit	1000706
Spindle/Impeller	1000583
Steel Bearing	1000604M (Front)
	1000605M (Rear)
O-Ring	1000902 (Front)
	1000903 (Rear)
Spring Washer(s)	1001223
Cap Tool	1001243

Spring Washer(s)	1001223
Cap Tool	1001243

## MIDWEST Stylus ATC 890 Standard Type



Turbine Steel	1000176
Ceramic	1000181
Rebuild Kit Steel	1000737
Ceramic	1000740
Spindle/Impeller	1000565
Spindle	1000237
Steel Bearing	1000611M
Ceramic Bearing	1000697M
O-Ring(s)	1000921
Spring Washer(s)	1001223
Back Cap	1000444
Cap Tool	1001253

## MIDWEST Quiet Air Lever Type



Turbine Steel	1000184
Rebuild Kit	1000704
Spindle/Impeller	1000584
Steel Bearing	1000604M (Front)
	1000605M (Rear)
O-Ring(s)	1000902 (Front)
	1000903 (Rear)
Cap Tool	1001243

## MIDWEST\* Tradition Screw Type



Turbine Steel	1000185
Rebuild Kit	1000704
O-Ring(s)	1000902 (Front)
	1000903 (Rear)
Bur Tool	1001239

## MIDWEST\* Quiet Air Screw Type



Turbine Steel	1000186
Rebuild Kit	1000706
O-Ring	1000902 (Front)
	1000903 (Rear)
Cap Tool	1001243
Bur Tool	1001241

12 Month warranty on Sable steel & ceramic bearing turbines.  
\*6 Month warranty on Midwest Tradition Screw Type.

12 Month warranty on Sable steel & ceramic bearing turbines.

# REPLACEMENT TURBINES BY BRAND

## MIDWEST\* Tradition Plus Type

Turbine Ceramic	1002102
Spindle/Impeller	1000591
Ceramic Bearing	1000697M
O-Ring(s)	1000921 (Front)
Bur Tool	1001253
Air Dam	1001254
Spring Washer(s)	1000123 (Front)
	1000123 (Rear)
PB Back Cap	1000488
Cap Tool	1001253

## MIDWEST Stylus Plus Phoenix PB Type

Turbine Ceramic	1002109
Back Cap	1002402
Cap Tool	1001256

## MIDWEST Tradition, Phoenix, Stylus (9.7MM Head) Ceramic PB Type Turbine

Turbine Ceramic	1002113
-----------------	---------

## MK DENT Large Head Triple Spray Type PB

Turbine Steel	1000128
Rebuild Kit	1000722
Spindle/Impeller	1000528
Spindle	1000228
Steel Bearing	1000601M
Ceramic Bearing	1000699M
O-Ring(s)	1000921
Back Cap	1000428

## NSK Ti-Max X500 / NSK NL9000M Type

Turbine Steel	1000167
Rebuild Kit	1000728
Spindle/Impeller	1000558
Steel Bearing	1000608M (Front)
	1000605M (Rear)
Ceramic Bearing	1000693M (Front)
	1000698M (Rear)
O-Ring(s)	1000902 (Front)
	1000912 (Rear)
Air Dam	1001231
Spring Washer(s)	1001223
Cap Tool	1001227

## MIDWEST Tradition PRO TLF Lever Type

Turbine Ceramic	1002107
Back Cap	1000494
Cap Tool	1001255

## MIDWEST Tradition Pro TBF PB Type Turbine

Turbine Ceramic	1002108
Back Cap	1000494
Cap Tool	1001255

## MK DENT Small Head Type

Turbine Steel	1000127
Rebuild Kit	1000722
Spindle/Impeller	1000527
Spindle	1000227
Steel Bearing	1000601M
Ceramic Bearing	1000699M
O-Ring (s)	1000921
Back Cap	1000427

## MK DENT Large Head Single Spray Type

Turbine Steel	1000151
Rebuild Kit	1000722
Spindle/Impeller	1000534
Spindle	1000228
Steel Bearing	1000601M
Ceramic Bearing	1000699M
O-Ring(s)	1000921

## NSK Ti-Max X600 / NSK, NL4500, NL9000S, X450, Sabra Power 45S Type

Turbine Steel	1000168
Rebuild Kit	1000729
Spindle/Impeller	1000559
Steel Bearing	1000605M
Ceramic Bearing	1000698M
O-Ring(s)	1000912
Air Dam	1001232
Spring Washer(s)	1001223
Cap Tool	1001228

## NSK Ti-MAX X700 / NSK NL9000T Type

Turbine Steel	1000169
Rebuild Kit	1000729
Spindle/Impeller	1000560
Steel Bearing	1000605M
Ceramic Bearing	1000698M
O-Ring(s)	1000912
Air Dam	1001233
Spring Washer(s)	1001223
Cap Tool	1001229

## NSK Ti-Max Z900L Type

Turbine Ceramic	1002100
Spindle/Impeller	1000589
O-Ring(s)	1000912 (Front)
	1000912 (Rear)
Air Dam	1001250
Spring Washer(s)	1001223 (Rear)
Back Cap	1000486
Back Cap Tool	1001249

## PALISADES IMPACT AIR 45 OEM Type For Use with OEM Back Cap

Turbine Steel	1002104
Conversion Kit	1001163
Rebuild Kit	1000716
Steel Bearing	1000604M
Ceramic Bearing	1000694M
O-Ring(s)	1000902 (Rear)
	1000906 (Front)
Back Cap	1002401

## SABLE ROTAMAX™

Turbine Steel	2009902
Rebuild Kit	1000720
Spindle/Impeller	2009502
Spindle	2009802
Steel Bearing	1000605M
Ceramic Bearing	1000698M
O-Ring(s)	1000920
Spring Washer(s)	1001223
Back Cap	2009602
Cap Tool	2001216

## NSK Ti-Max Z800L Type

Turbine Ceramic	1002101
Spindle/Impeller	1000590
O-Ring(s)	1000912 (Front)
	1000912 (Rear)
Air Dam	1001252
Spring Washer(s)	1001223 (Rear)
Back Cap	1000487
Back Cap Tool	1001251

## PALISADES IMPACT AIR 45 Type (Requires Sable Back Cap).

Turbine Steel	1000158
Conversion Kit	1001158
Rebuild Kit	1000716
Steel Bearing	1000604M
Ceramic Bearing	1000694M
O-Ring(s)	1000902 (Rear)
	1000906 (Front)
Back Cap	1000488

## SABLE ACCESS™ ACCESS II™

Turbine Steel	2009901
Rebuild Kit	1000719
Spindle/Impeller	2009501
Spindle	2009801
Steel Bearing	1000611M
Ceramic Bearing	1000697M
O-Ring(s)	1000920
Spring Washer(s)	1001223
Back Cap	2009601 (Access)
	2009603 (Access II)
Cap Tool	2001215

## SABLE ROTAMAX II™

Turbine Steel	2009903
Rebuild Kit	1000719
Spindle/Impeller	2009503
Spindle	2009803
Steel Bearing	1000611M
Ceramic Bearing	1000697M
O-Ring(s)	1000920
Spring Washer(s)	1001223
Back Cap	2009604
Cap Tool	2001216

12 Month warranty on Sable steel & ceramic bearing turbines.  
\*6 Month warranty on Midwest Tradition Screw Type.

12 Month warranty on Sable steel & ceramic bearing turbines.

# REPLACEMENT TURBINES BY BRAND

## SABLE ACCESS L™ ACCESS KL™



Turbine Steel	2009907
Rebuild Kit	1000725
Steel Bearing	1000611M
Ceramic Bearing	1000697M
O-Ring(s)	1000916 (x1) (Front) 1000915 (x2) (Rear)
Spring Washer(s)	1001223
Back Cap	2009605
Cap Tool	2001221

## SABLE ROTAMAX L™ ROTAMAX KL™



Turbine Steel	2009908
Steel Bearing	1000612M
O-Ring(s)	1000917
Back Cap	2009606
Cap Tool	2001222

## SIEMENS 4000SL, TS1 Type



Turbine Steel	1000108
Rebuild Kit	1000710
Spindle/Impeller	1000508
Spindle	1000208
Steel Bearing	1000601M
Ceramic Bearing	1000699M
O-Ring(s)	1000921
Back Cap	1000408

## SIEMENS 4000ML, TM1 Type



Turbine Steel	1000123
Rebuild Kit	1000710
Spindle/Impeller	1000523
Spindle	1000223
Steel Bearing	1000601M
Ceramic Bearing	1000699M
O-Ring(s)	1000921

## SABLE S600NL STANDARD Ceramic



Turbine Ceramic	2009910
Rebuild Kit	1000727
Spindle/Impeller	1000556
Steel Bearing	1000605M
Ceramic Bearing	1000698M
O-Ring(s)	1000920
Spring Washer(s)	1001223
Back Cap	2009608
Cap Tool	2001214

## SABLE S700NL TORQUE Ceramic



Turbine Ceramic	2009911
Rebuild Kit	1000727
Spindle/Impeller	1000557
Steel Bearing	1000605M
Ceramic Bearing	1000698M
O-Ring(s)	1000920
Spring Washer(s)	1001223
Back Cap	2009609
Cap Tool	2001216

## SIRONA T1 Control Highspeed Type



Turbine Steel	1000154
Spindle/Impeller	1000554
Spindle	1000254
O-Ring(s)	1000921

## SIRONA T2 CONTROL Highspeed Type



Turbine Steel	1000155
Spindle/Impeller	1000555
Spindle	1000255
O-Ring(s)	1000902
Steel Bearing	1000611M (x2)

## SABLE ACCESS™ S4500KL Ceramic



Turbine Ceramic	2009912
Rebuild Kit	1000730
Spindle/Impeller	2009504
Ceramic Bearing	1000697M
O-Ring(s)	1000919
Spring Washer(s)	1001223
Back Cap	2009610
Cap Tool	2001234

## SABLE ROTAMAX™ S6500KL Ceramic Type



Turbine Ceramic	2009913
Rebuild Kit	1000730
Spindle/Impeller	2009505
Ceramic Bearing	1000697M
O-Ring(s)	1000919
Spring Washer(s)	1001223
Back Cap	2009611
Cap Tool	2001234

## SIRONA T2 & T3 RACER Type



Turbine Steel	1000163
O-Ring(s)	1000907
Steel Bearing	1000611M
Ceramic Bearing	1000697M

## SIRONA T2 Mini / ProDrive PD-LMS / Patterson Mini Head Type



Turbine Steel	1000190
O-Ring(s)	1000907
Ceramic Bearing	1000697M

## SABLE ACCESS™ PRO-NL & KL Mini



Turbine Ceramic	2009914
Back Cap	2009612
Cap Tool	2001253

## SABLE ROTAMAX™ PRO-NL & KL Torque



Turbine Ceramic	2009915
Back Cap	2009613
Cap Tool	2001253

## SIRONA T3 STANDARD PUSH BUTTON LATCH HEAD



Turbine Ceramic	1002105
Back Cap	1000492

## SIRONA T3 MINI PUSH BUTTON LATCH HEAD



Turbine Ceramic	1002106
Back Cap	1000493

## SABLE ROTAMAX™ 45 & 45KL



Turbine Ceramic	2009916
Back Cap	2009616

## STAR 430/431 Type



Turbine Steel	1000106
Conversion Kit	1001106
Rebuild Kit	1000711 (w clips) 1000712 (w/o clips)
Spindle/Impeller	1000506
Spindle	1000206
Steel Bearing	1000609M
O-Ring(s)	1000906
Back Cap	1000406
Cap Tool	1001243

## Star 430 Torque Ceramic PB Type



Turbine Ceramic	1002112
PB Back Cap	1002403



# REPLACEMENT TURBINES BY BRAND

## STAR **LUBE FREE**



Ceramic Type	
Turbine Ceramic	1000162
Spindle/Impeller	1000506
Spindle	1000206
Ceramic Bearing	1000692M (x2)
O-Ring(s)	1000906
Back Cap	1000406
Cap Tool	1001243

## VECTOR Small Head Type



Turbine Steel	1000129
Rebuild Kit	1000723
Spindle/Impeller	1000529
Spindle	1000229
Steel Bearing	1000601M
Ceramic Bearing	1000699M
O-Ring(s)	1000921
Back Cap	1000429

## W&H 300 SERIES, SABRA Type



Turbine Steel	1000113
Rebuild Kit	1000714
Spindle/Impeller	1000513
Spindle	1000204
Steel Bearing	1000605M
O-Ring(s)	1000921
Spring Washer(s)	1001223
Back Cap	1000404

## W&H SYNEA TA-96 Type



Turbine Steel	1000118
Rebuild Kit	1000713
Spindle/Impeller	1000518
Spindle	1000220
Steel Bearing	1000611M
O-Ring(s)	1000907
Spring Washer(s)	1001223
Cap Tool 2 Notch	1001215
2 Notch Back Cap	1000482

## VECTOR Large Head Triple Spray Type



Turbine Steel	1000130
Rebuild Kit	1000723
Spindle/Impeller	1000530
Spindle	1000230
Steel Bearing	1000601M
Ceramic Bearing	1000699M
O-Ring(s)	1000921
Back Cap	1000430

## VECTOR Large Head Single Spray Type



Turbine Steel	1000152
Rebuild Kit	1000723
Spindle/Impeller	1000530
Spindle	1000230
Steel Bearing	1000601M (x2)
Ceramic Bearing	1000699M (x2)
O-Ring(s)	1000921
Back Cap	1000430

## W&H TA-98L Type



Turbine Steel	1000119
Rebuild Kit	1000713
Spindle/Impeller	1000519
Spindle	1000221
Steel Bearing	1000611M
O-Ring(s)	1000907
Spring Washer(s)	1001223
Back Cap 2 Notch	1000471
2 Notch Cap Tool	1001225

## W&H TREND TC-95RM Type



Turbine Steel	1000124
Rebuild Kit	1000713
Spindle/Impeller	1000524
Spindle	1000225
Ceramic Bearing	1000697M
Steel Bearing	1000611M (x2)
O-Ring(s)	1000907

## W&H TOPAIR 795, 798, 895, 898 Type



Turbine Steel	1000104
Rebuild Kit	1000713
Spindle/Impeller	1000504
Spindle	1000204
Steel Bearing	1000611M (x2)
O-Ring(s)	1000907
Spring Washer(s)	1001223
Back Cap	1000404
Cap Tool	1001216

## W&H TOPLIGHT 796, 896 Type



Turbine Steel	1000105
Rebuild Kit	1000713
Spindle/Impeller	1000505
Spindle	1000205
Steel Bearing	1000611M
O-Ring(s)	1000907
Spring Washer(s)	1001223
Back Cap	1000405
Cap Tool	1001215

## W&H Adec Synea TA-97 LW, TA-97 LED, TK-97L, TG-97 (L, LW, LM, LN) 400 Series Type



Turbine Steel	1000164
Rebuild Kit	1000713
Spindle/Impeller	1000535
Spindle	1000241
Steel Bearing	1000611M
O-Ring(s)	1000907
Spring Washer(s)	1001223
Back Cap 3 Notch	1000473
3 Notch Cap Tool	1001226

## W&H Adec Synea TA-98 LW, TG98 TA-98 LED, TK-98L, TK98 (L, LW, LM, LN), 400 Series Type



Turbine Steel	1000172
Rebuild Kit	1000713
Spindle/Impeller	1000536
Spindle	1000240
Steel Bearing	1000611M
Ceramic Bearing	1000697M
O-Ring(s)	1000907
Spring Washer(s)	1001223
Back Cap 3 Notch	1000472
3 Notch Cap Tool	1001235

## Sable W & H Alegra TE-97 PB Type


Turbine Steel	1002110
---------------	---------

## Sable W & H Alegra TE-98 PB Type

Turbine Steel	1002111
---------------	---------

The flow table shows bearing variants available. By matching the bearing removed from the handpiece with the illustrations shown, it is possible to select the bearing type. Special note should be taken of number of shields and shield location, flanged or non flanged, special OD grooving and any other unusual features.

“S” Steel Ball Bearing  
“C” Ceramic Ball Bearing  
\* Angular Contact



NSK Pana-Air Mini Canister  
**1001301**  
NSK MODEL# NPA-M03



NSK Pana-Air Standard Canister\*  
**1001302**  
NSK MODEL# NPA-S03



NSK Pana-Air Torque Canister  
**1001303**  
NSK MODEL# NPA-T03



## NSK & OEM CANISTER PB TYPE TURBINES

NSK PANA-AIR Mini Canister  
**1001401**  
NSK MODEL# NPA-MU03



NSK PANA-AIR Standard Canister\*  
**1001402**  
NSK MODEL# NPA-SU03



NSK PANA-AIR Torque Canister  
**1001403**  
NSK MODEL# NPA-TU03



NSK PANA 2000 Mini Canister  
**1001409**  
NSK MODEL# PANA 2000-NU03



NSK PANA 2000 Torque Canister  
**1001410**  
NSK MODEL# PANA 2000-SU03



\* Also fits Champion, Healthco, Henry Schein, Darby, Naminex, Medidenta, Dental Aire.

**PUSH BUTTON TYPE**

PANA MAX Standard  
BRASSLER N5000

**1001496**  
NSK MODEL# PAX-SU03



NSK PHATELUS I Mini Twist-Free Mini  
ALPHALITE Mini KINETIC VIPER, Mini

**1001497**  
NSK MODEL# NTF-MU03



NSK PHATELUS I STANDARD TWIST-FREE  
Standard ALPHALITE Standard  
KINETIC VIPER, Standard

**1001498**  
NSK MODEL# NTF-SU03



NSK PHATELUS I Torque Twist-Free  
Torque ALPHALITE Torque  
KINETIC VIPER, Torque

**1001499**  
NSK MODEL# NTF-TU03



NSK PHATELUS II Mini Canister  
NL-65 / CH-QD

**1001404**  
NSK MODEL# NCH-MU03



NSK Phatelus II NL-75  
45 Degree Standard & 45-S  
Surgical HPC AK450L

**1001405**  
NSK MODEL# NCH-SU03



NSK PHATELUS II Torque Canister  
NL-65 / CH-QD

**1001406**  
NSK MODEL# NCH-TU03



HEALTHCO Push And Swivel  
HEALTHCO PUSH, Swivel And Light  
NSK VIP-AIR / VIP II

**1001407**  
NSK MODEL# NVA-SU03



NSK MACH-LITE XT Mini Canister  
MACH-LITE / MACH / MACH-QD  
NL-75 / NL-85 AK200L

**1001501** NSK MODEL# NMC-MU03



NSK MACH-LITE XT Standard Canister  
MACH-LITE / MACH / MACH-QD  
AK300L, NL-75 / NL-85

**1001502** NSK MODEL# NMC-SU03



NSK MACH-LITE XT Torque Canister  
MACH-LITE / MACH / MACH-QD  
NL-75 / NL-85

**1001503** NSK MODEL# NMC-TU03



MICRO MOTORS Ih730  
MINI CANISTER LHS-720

**1001504**



MICRO MOTORS LHS700, 715  
Standard Canister LHS-715 LHS-700

**1001505**



SABLE NSK CANISTERS - USA grade ABEC 7 Angular  
Contact Bearings 6 Month Manufacturers Warranty.

## NSK ECONOMY CANISTER TURBINES

NSK TI-MAX A500L,  
AK500L, NL-95M  
**1001506**  
NSK MODEL# TI-MU03



NSK TI-MAX A600L,  
AK600L, NL-95S  
**1001507**  
NSK MODEL# TI-SU03



NSK TI-MAX A700L, AK700L, NL-95T  
**1001508**  
NSK MODEL# TI-TU03



NSK S-MAX M500 TYPE  
PB Mini Canister  
**1001509**  
NSK MODEL# SMAX MU03



NSK S-MAX M600 TYPE  
PB Standard Canister  
**1001510**  
NSK MODEL# SMAX-SU03



### CHUCK TYPE

NSK PANA-AIR Mini Canister  
**1001601**  
NSK MODEL# NPA-M03



NSK PANA-AIR Standard Canister+  
**1001602**  
NSK MODEL# NPA-S03



NSK PANA-AIR Torque Canister  
**1001603**  
NSK MODEL# NPA-T03



SAFCO Super Grade Canister  
**1001604**



+ Also fits Champion, Healthco, Henry Schein, Darby, Naminex, Medidenta, Dental Aire. USA grade ABEC 7 Radial Bearings.

\* No Manufacturers Warranty offered.

## NSK ECONOMY CANISTER TURBINES

### PUSH BUTTON TYPE

NSK PANA-AIR Mini Canister  
**1001701**  
NSK MODEL# NPA-MU03



NSK PANA-AIR Standard Canister  
**1001702**  
NSK MODEL# NPA-SU03



NSK PANA-AIR Torque Canister  
**1001703**  
NSK MODEL# NPA-TU03



SAFCO Super Grade Canister  
**1001704**



### NSK & OEM BACK CAPS

NSK Chuck Type PANA-Air Mini & Standard\*  
Fits w/ Sable Part #'s 1001301, 1001302, 1001601, 1001602  
**1000431**

NSK PB Type PANA-Air Mini & Standard\*  
Fits w/ Sable Part #'s 1001401, 1001402, 1001701, 1001702  
**1000432**

NSK PB Type PHATELUS II Mini & Standard  
NSK PB Type MACH-LITE Mini & Standard  
Fits w/ Sable Part #'s 1001404, 1001405, 1001501, 1001502  
**1000433**

+ Also fits Champion, Healthco, Henry Schein, Darby, Naminex, Medidenta, Dental Aire.

## HANDPIECE REPLACEMENT BULBS

KAVO LUX Coupler  
OEM code 553 3881  
SINGLE **1000801**  
3 PACK **1000802**



NSK OEM code Y900 529  
SINGLE **1000803**  
3 PACK **1000804**



MIDWEST MW Power Optic  
In-Hose Swivel OEM code 650064  
SINGLE **1000810**  
2 PACK **1000809**



MIDWEST XGT/Stylus/Stylus ATC  
OEM code 790257  
SINGLE **1000818**



MK-DENT MK Dent QC 5016  
Coupler OEM code BU7012  
SINGLE **1000812**  
3 PACK **1000811**



W&H Type Single  
**1000821**



DABI ATLANTE MSL  
OEM code 42001-021/2  
SINGLE **1000807**  
3 PACK **1000808**



5 Hole Tubing Bulb  
Fits Sable Part #2100201 / 2100202  
**1000815**



NLC 400 (Sable NSK Type Coupler)  
**1000819**



VECTOR  
Type Single **1000814**  
Type (3 Pack) **1000813**

## LED FIBER OPTIC HANDPIECE REPLACEMENT BULBS

KAVO Type Single LED **1000850**  
3 Pack LED **1000851**



KAVO Motor Type Single LED **1000852**  
3 Pack LED **1000853**

NSK Type Single LED **1000864**  
3 Pack LED **1000865**



MIDWEST Type Single LED **1000866**  
XGT/Stylus/Stylus ATC 3 Pack LED **1000867**



W&H Type Single LED **1000854**  
3 Pack LED **1000855**



BIEN AIR Motor Type Single LED **1000858**  
3 Pack LED **1000859**



BIEN AIR Type Single LED **1000860**  
3 Pack LED **1000861**

SIRONA Type Single LED **1000856**  
3 Pack LED **1000857**

SIRONA Motor Type Single LED **1000862**  
3 Pack LED **1000863**

**SABLE LIGHT CONTROL**

- Power Source: 100-240V, 9VDC; 2.2 Amp
- Adjustable Output: DC 2V-3.5V
- Adjustable Illumination time delay:  
0 - 20 seconds for each handpiece
- 1 or 2 Optic Handpieces
- Operates All Brands Of Handpieces
- Weight: 4 oz.
- W 2.25" x L 3.5" x H 1"

**SABLE Light Control 7FT Tubing Kit** Includes:  
2300101: Sable Light Control Board, 1-2 Handpieces  
2300203 SABLE Power Transformer 100-240V; 9VDC; 2.2 Amp  
2100103 6 Hole, 7 Foot, Straight GRAY Tubing  
**2300001**

**SABLE Light Control 10FT Tubing Kit** Includes:  
2300101: Sable Light Control Board, 1-2 Handpieces  
2300203 SABLE Power Transformer 100-240V; 9VDC; 2.2 Amp  
2100105 6 Hole, 10 Foot, Straight GRAY Tubing  
**2300002**

**SABLE Light Control 5-Hole Tubing Kit** Includes:  
2300101: Sable Light Control Board, 1-2 Handpieces  
2300203 SABLE Power Transformer 100-240V; 9VDC; 2.2 Amp  
2100202 5 Hole, 7 Foot Straight GRAY Tubing with Bulb (1000815)  
**2300003**

Sable Light Control Board  
1 or 2 Optic Handpieces  
Adjustable Output: DC 2V-3.5V  
Adjustable Illumination  
Time Delay: 0 - 20 seconds  
for each handpiece  
**2300101** OEM P/N 905-8400



Sable Power Transformer  
9VDC, 100-240V, 2.2 Amp  
**2300203**

9VDC 2.2 Amp wall type military spec.  
Power transformer.  
Any STAR Dental Halogen Bulb can be  
used with our light control board.



Sable Tubing Hook-Up Kit  
**2300202**



Tubing can be substituted for different colours of equal value at no charge.



GRAY STERLING SURF DARK  
(WHITE) SURF SURF

**HANDPIECE TUBING ISO FIBER OPTIC**

**6 Hole, 7 Foot, Straight;** OEM P/N KL035009

Gray	<b>2100103</b>
Sterling (White)	<b>2100108</b>
Surf (Light Brown Or Sand)	<b>2100104</b>
Dark Surf (A-dec Type)	<b>2100110</b>

**6 Hole, 10 Foot, Straight;** OEM P/N KL035010

Gray	<b>2100105</b>
Sterling (White)	<b>2100109</b>
Surf (Light Brown Or Sand)	<b>2100106</b>
Dark Surf (A-dec Type)	<b>2100111</b>

**5 Hole, 7 Foot, Straight;**

Gray W/ Bulb	<b>2100202</b>
Sterling (White) W/ Bulb	<b>2100201</b>

**Replacement Bulb;**

Bulb	<b>1000815</b>
------	----------------

**SABLE BRAND EURO TYPE TUBING**

Silicone GRAY Tubing with 10-32 to Metric Barb

<b>5-Foot, 4-Hole</b>	<b>2100301</b>
<b>5-Foot, 6-Hole</b>	<b>2100302</b>
<b>7-Foot, 4-Hole</b>	<b>2100304</b>
<b>7-Foot, 6-Hole</b>	<b>2100305</b>



**SABLE BRAND ADEC TYPE TUBING**

52 1/4", ADEC type Silicone GRAY Tubing with Quick Disconnect Fittings


<b>4-Hole</b>	<b>2100401</b>
<b>6-Hole</b>	<b>2100402</b>





**SABLE BUCK-BOOST TRANSFORMERS**

- Steps up or steps down incoming Voltage.
- Made in Canada. CSA & UL approved.
- Variety of applications such as compressors, vacuums, sterilizers x-rays and more.
- Include Simple step by step bilingual install directions and drawing.
- Front accessible separate high and low voltage lead wires.
- Valid for single and three phase.
- 10 Year Warranty.



Sable Buck-Boost Transformers represent an economical way to both raise supply voltage caused by line drop and equipment demands or lower voltage caused by increased system voltages.

**Select the correct buck-boost transformer:**

1. Line Voltage – The voltage that you want to buck (decrease) or boost (increase). This can be done by measuring the supply line voltage with a voltmeter.
2. Load Voltage – The voltage at which your equipment is designed to operate. This information can be found on the nameplate of the equipment that you want to operate.
3. Load kVA or Load Amps – You do not need to know both – one or the other is sufficient. This information can be found on the nameplate of the equipment.

3500KVA, (0.35 KVA), 120/240V-12/24,  
50/60Hz, CU NEMA-3R, 80 C, 4.5"W x 5.75"D x 7.25"H  
Wall Mount 14 LBS, 3500 Type **3001000**

5000KVA, (0.50 KVA), 120/240V-12/24,  
50/60 Hz, CU NEMA-3R, 80 C, 5.00"W x 4.75"D x 9.25"H  
Wall Mount 15 LBS, 5000 Type **3001010**

7500KVA, (0.75 KVA), 120/240V-12/24,  
50/60Hz, CU NEMA-3R, 80 C, 5.00"W x 4.75"D x 9.25"H  
Wall Mount 18 LBS, 7500 Type **3001020**

**NEW 3 WIRE BUCK - BOOST**

Three Wire, 5000KVA  
Ideal for single head compressors or suction units **3001015**

Three Wire, 7500KVA  
Ideal for dual head compressors or dual head suction units **3001025**



Please refer to Buck-Boost Documentation on Sable Website.

**Autoscaler 30K  
1900001**

- Rated EXCELLENT by a leading research institute.
- Accepts SABLE, Hu-Freidy, Whaledent, Premier, and Dentsply/Cavatron brand tips.
- Includes a detachable handpiece assembly for ease of service.
- Extended 6 year warranty on circuitry.



**AUTOSCALERS REPLACEMENT PARTS**

Replacement Handpiece 25K/30K Units w/serial's 9676-current	<b>1900102</b>
Solenoid & Regulator 25K/30K Units w/serial's 11732-current	<b>1900202</b>
Foot Control for Autoscaler Black (All Units)	<b>1900405</b>

**ULTRASONIC SCALING TIP**

30,000 Hz Hu Friedly Ultrasonic Insert #10 - Universal  
**1900603**



**AUTOCLAVE PARTS**

**STATIM COMPRESSORS**

Compressor  
Statim 2000  
OEM code 01-100560S  
Pelton & Crane Delta 8, 10  
**1700101**



Compressor Filter Statim 2000  
OEM code 01-100207S Pelton & Crane Delta 8, 10  
(Pkg. of 10) **1700201**



**AUTOCLAVE GASKETS**

SciCan Statim 2000 Cassette Type Gasket Seal	<b>1700301</b>
SciCan Statim 5000 Cassette Type Gasket Seal	<b>1700302</b>
MidMark M11 Door Type Gasket Seal with Dam	<b>1700307</b>



For eVac Maintenance or Between Patient Flush



**Guaranteed To Improve  
Suction Performance!**

**EPA Amalgam Rule Compliant pH6-8**

Cleans 24/7 - 1 Microbe becomes 10,000,000 in 24 hours

- The Right Microbes and Enzymes together for clean vacuum lines
- All natural & non-foaming
- Compatible & best choice for all Amalgam Separators

**48 oz. = 533 treatments**



## BIO-PURE EVACUATION SYSTEM CLEANER



### BIO-PURE EVAC CLEANER STARTER KITS

Includes: Maintenance Jar, and HVE-SE Dispenser.  
Restore 6 oz. Kit case of 6 (67 Treatments) **2801101**

### BIO-PURE EVAC SYSTEM CLEANER

Maintenance Cleaner 6 oz. case of 12 (67 Treatments) **2801001**  
Maintenance Cleaner 14 oz. case of 12 (156 Treatments) **2801002**  
Maintenance Cleaner 48 oz. case of 6 (533 Treatments) **2801003**  
Bio-Pure eVac System Dispenser case of 6 **2801004**

### BIO-PURE - LIQUID EVAC CLEANER (LSF)

Bio-Pure Liquid eVac Cleaner 32oz. (128 Treatments) **2801402**  
Bio-Pure Liquid eVac Cleaner Starter Kit (128 Treatments) **2801405**  
(**2801405: 2801402** 32 oz. Cleaner & **2801004** 48oz. Dispenser)

### BIO-PURE - BETWEEN PATIENT FLUSH (BPF)

Bio-Pure Between Patient Flush 32oz. (128 Treatments) **2801202**  
Bio-Pure Between Patient Flush 64oz. (256 Treatments) **2801203**  
Bio-Pure Between Patient Flush Starter Kit (128 Treatments) **2801206**  
(**2801206: 2801202** 32 oz. Cleaner & **2801004** 48oz. Dispenser)

## BIO-PURE EPA COMPLIANCE PLUS SOFTWARE (USA Only)

Bio-Pure EPA Compliance Plus Software **2801600**  
1 Year License  
Bio-Pure EPA Compliance Plus Software Renewal **2801601**  
1 Year License Renewal

## ULTRASONIC 3-IN-1 PRESOAK & CLEANER



Bio-Pure is a unique proprietary environmentally friendly ultrasonic cleaner safely removing blood, tissue, proteins and contaminants from instruments, burs, dentures, handpieces, aspirator tips, metal parts, bridges, diamonds, orthodontic appliances, gold work, porcelain, crowns, ultrasonic tips and much more.

### Effective!

- Triple Action Enzymes & Proprietary Surfactants.
- Contains proprietary deep cleaning ingredients that energize cleaning.
- Removes organic & inorganic deposits rapidly while protecting delicate instruments.

### Great Value!

- 3-in-1 Cleaner: Presoak Immersion Tub, Presoak Spray, Ultrasonic Unit.
- Concentrated 1 oz. makes 1.5 gal. 32 oz. bottle makes 48 gallons of solution (1:192 ratio).
- Works in warm or cold water.

### Naturally Safe!

- pH Neutral.
- Non-Corrosive, Non-Foaming & Non-Toxic.
- No Fragrance Contaminants.
- No Ammonia or Phosphates.

Bio-Pure Ultrasonic Presoak & Cleaner 32oz. **2801302**



**TURBINES**

- Sable Replacement Steel or Ceramic Bearing Turbines Includes, Spindle, Chuck, Impeller, or Bearings ..... 12 Months
- Sable Conversion Turbines Steel or Ceramic Includes
  - Spindle, Chuck, Impeller, Back Cap or Bearings ..... 12 Months
  - Manual Chuck Type Canister Turbines ..... 4 Months
  - NSK Type PB Type Canister Turbines ..... 6 Months
  - NSK Type Economy Canister Turbines ..... No Warranty
  - Midwest Screw Type Turbines ..... 6 Months

**1:5 ELECTRIC PARTS** ..... 12 Months

**TURBINE SERVICE ITEMS**

Caps, Spindles, and Spindle/Impeller Combos: ..... 12 Months

**BEARINGS**

(Loose Bearings) All High-speed Steel, Ceramic, Standard Slow Speed Bearings ..... No Warranty

**REBUILD KITS**

All rebuild kits ..... No Warranty

**O-RINGS**

All O-rings ..... Manufacturer Defect Only

**BULBS**

Handpiece Replacement Bulbs ..... Manufacturer Defect Only

**TUBINGS**

Handpiece Tubing ..... Manufacturer Defect Only

**SLOW SPEED**

- KaVo Type Premium Slow Speed products (ALL) ..... 12 Months
- Slow Speed Air Motors ..... 12 Months
- Slow Speed Contra Angles and Nosecones ..... 6 Months
- Slow Speed Packages (NSK or KaVo Type) ..... 12 Months
- Slow Speed Heads ..... 3 Months
- U Type Contra Angles ..... 3 Months
- New One Piece ..... 12 months

**HIGH SPEED**

- High Speed Handpiece ..... 12 Months
- New PRO Series High Speed Handpieces ..... 24 Months
- High Speed Couplers ..... 12 Months
- 45 Degree Handpiece ..... 12 Months
- Light Control Set ..... 12 months (on the Light Box and Transformer

**AUTOSCALER**

- South East Instruments Autoscaler ..... 1 Year Parts & Labour
- Circuitry ..... 6 Year Extended Warranty
- Handpiece Tips ..... Manufacturer Defect Only
- HP Tubing Assembly ..... Manufacturer Defect Only
- Buck - Boost Transformers ..... 10 Year Warranty

**PRODUCT RETURN PROCEDURE**

A Return Authorization is required for all returns. Please contact Sable Industries Inc. to obtain an authorization number at 1-800-368-8106.

*Customers are responsible for Return Shipping Costs. No Returns on Special Orders!*

Please indicate RA# on the outside of the carton. Please accompany returns with pertinent invoices and paperwork. Refunds will be credited to the account minus shipping costs. Products not required or ordered in error could be subject to a 25% restocking fee, all merchandise must be in new, unused condition to receive credit.

*All USA Repairs and Returns are through our USA Authorized Repair Centre:*

**Sable USA,**  
27300 Haggerty Road, Suite F-1  
Farmington Hills, MI 48331 USA  
Phone # 1-734-293-7503

**CONDITIONS**

Sable Industries Inc. reserves the right to review warranty claims, request supporting documentation, and/or substitute product or remaining portion of original cost.

Return Address: **Sable Industries Inc.**  
100 Campbell Avenue Suite 5  
Kitchener, On N2H 4X8 Canada  
P: 1-800-368-8106 F: 1-519-579-9324

**PRICING**

A minimum order of \$250.00 is required, if minimum is not met, a \$30.00 processing fee will apply. We make every effort to maintain prices for the duration of the catalogue, but we reserve the right to make price adjustments to any increase in transportation, labour, material, or other cost. We reserve the right to correct obvious errors in specifications or prices.

**EXPORT COMPLIANCE**

Any subsequent export, transfer, resale or other disposition of the products from the United States /Canada must be made in accordance with applicable law. Purchaser shall not export or re-export the products in violation of any export control laws, or without all necessary authorizations and/or licenses. Purchaser shall be solely responsible for obtaining and paying for such required licenses or approvals from the appropriate governmental agencies or authorities. Purchaser shall indemnify, defend and hold seller harmless from and against any and all losses, damages, costs and expenses (including attorneys' fees) incurred by seller by reason of any failure of purchaser to comply with such export control laws. Seller shall be relieved of all obligations to provide products to purchaser if purchaser violates any export control laws, or if any related approvals/license are suspended or revoked.



**SABLE**  
INDUSTRIES INC.



**SABLE HEAD OFFICE**

100 Campbell Avenue, Suite 5,  
Kitchener, Ontario, N2H 4X8 CANADA

Tel: 519-579-9323 / 1-800-368-8106

Fax: 519-579-9324

Email: [info@sableindustriesinc.com](mailto:info@sableindustriesinc.com)

**SABLE USA AUTHORIZED REPAIR CENTRE**

27300 Haggerty Road, Suite F-1  
Farmington Hills, MI 48331 USA

**WWW.SABLEINDUSTRIESINC.COM**

PRINTED IN CANADA LIT-0201