

**NEW**



# PHOTON EXE PHOTON EXE PLUS

Soft-Tissue Dental Diode Laser

Enhanced Precision & Redefined Excellence

Updated  
LBR Timer

Interactive  
User  
Inference

Enhanced  
Battery  
Performance

Customized  
Treatments

High-Resolution  
Large Display

30+ Additional  
Presets

Training &  
Customer Service  
Available



SCAN



Contact us or  
scan the QR code

+905-593-3605  
info@zolartek.com  
www.zolartek.com

# APPLICATIONS

- Exposure of Unerupted Teeth
- Fibroma Removal
- Gingival Toughing
- Crown Impressions
- Gingivoplasty
- Gingival Incision
- Implant Recovery
- Incision and Drainage of Abscess
- Operculectomy
- Reduction Gingival Hypertrophy
- Soft tissue Crown Lengthening

- Gingival Bleeding Index
- Excisional in Incisional Biopsies
- Vestibuloplasty
- Leukoplakia
- Oral Papillectomies
- Crown Lengthening
- Heptic Ulcer Treatment
- Fibroma Removal and more
- **LLLT: Low Level Laser Therapy Treatments (Photon Exe)**
- **Teeth Whitening (Photon Exe Plus)**

# TECHNICAL SPECIFICATIONS

<b>Settings</b>	Manual (Power, Time, Mode) More than 30 preset	<b>Technical Information</b>	<b>Photon:</b> Wavelength: 810nm Power: 3W Weight: <2kgs Dimensions: 187 × 148 × 38 mm <b>Photon Plus:</b> Wavelength: 980nm Power: 10W Weight: <2kgs Dimensions: 187 × 148 × 38 mm
<b>Modes</b>	Continuous or Pulse Pulse Interval: 0.01ms to 9.9s Pulse Duration: 0.01ms to 9.9s Duty Cycle: Variable		
<b>Power Requirement</b>	AC100 to 240V Lithium battery 16.8V - 10AH Wireless Pedal 1.5V Battery		
		<b>Laser Type</b>	Dental Diode Laser

# ACCESSORIES

LLLT Handpiece	Pain Therapy Tip	Cutting Fiber	Foot Pedal	Permanent Handpiece for Disposable tips
				

