

BASE series

3 SECTION ELECTRIC TREATMENT TABLE

This table is designed for simple examination at gynecology procedures. It suitable for a wide range of treatment and other therapy.




Product details



Product features

- Electric Hi-Lo operation by foot bar
- Sleek design with rounded frame legs for safety
- Head section has face hole and plug
- Equip with stirrups on both side
- Adjustable feet for uneven surfaces
- 75mm/3" retractable castors for safety and mobility
- Reliable TIMOTION motor
- Height range from 48-92cm/19"-36" in approx 30 seconds
- Padding thickness 5cm/2"
- Foam Density: 36kg/m³(79.41lbs/m³)
- Warranty: 5 Years Construction, 3 Years Electrics
- Upholstery: PVC vinyl that is fire retardant and meet the requirement of BS 5852, stain, mildew, and wear resistance 150,000 times
- Qualification: ISO13485/ Passed MDR CE certification / Registered with the FDA in the United States
- Weight Capacity: 250kgs/551lbs

Technical data



| Angle of head section | Angle of leg section | Wheel | Product weight | Volume | Packing size |
|-----------------------|----------------------|-------|----------------|--------------------|-------------------|
| -45°~ +40° | 0°~ +80° | φ75mm | 62kgs | 0.83m ³ | 189x70x62cm |
| -45°~ +40° | 0°~ +80° | φ3" | 136.71lbs | 0.83m ³ | 74.4"x27.6"x24.4" |

Optional accessories



Product model

BASE | EL03S

| | |
|---------------------|---------------------------------------|
| LxWxH | 195x66x(48-92)cm 76.8"x26"x19"-36" |
| Control mode | motor |
| Height adjustment | foot bar |
| Adjustable foot bar | 360° |
| Lifting speed | 30s |
| Down speed | 30s |

Color of vinyl

