

Wallach®

Endocell® Provides Simple, Safe Endometrial Sampling

Disposable Cell Sampler Screens for a Wide Range of Conditions

Endocell Disposable Endometrial Sampler

Endocell® is an easy-to-use, disposable cell sampler that requires only manual suction to safely extract a histological sample from the endometrium.

The simplicity of Endocell enables all practitioners to perform the endometrial screening procedure. For patients, Endocell is a simple, non-traumatic test that can be used in the early detection of cancer of the uterus.

- Excellent suction provided by elastomeric seal piston plunger
- No cervical dilation necessary in most cases
- Numbered sounding marks for safety and precision
- 3.1mm outside diameter
- Sterilized and individually packaged for convenience
- Flexibility promotes patient comfort



| Item # | Description |
|----------|---|
| 908014 A | Endocell box of 35, sterile |
| 908015 | Endocell 6 boxes of 35 (210 pieces), sterile |
| 908016 | Endocell 15 boxes of 35 (525 pieces), sterile |

CooperSurgical

800.243.2974 | 203.601.5200
www.coopersurgical.com/distribution
95 Corporate Drive, Trumbull, CT 06611

©2018 CooperSurgical, Inc.

82502 Rev. 4/18